

REF			SYSTEM
09005021190*	09005021500	300	cobas e 402 cobas e 801
09005021214*			

* Some kits shown may not be available in all countries.

English

System information

Short name	ACN (application code number)
ATG	10202

Please note

The measured anti-Tg value of a patient's sample can vary depending on the testing procedure used. The laboratory finding must therefore always contain a statement on the anti-Tg assay method used. Anti-Tg values determined on patient samples by different testing procedures cannot be directly compared with one another and could be the cause of erroneous medical interpretations. If there is a change in the anti-Tg assay procedure used while monitoring therapy, then the anti-Tg values obtained upon changing over to the new procedure must be confirmed by parallel measurements with both methods.

Intended use

Immunoassay for the in vitro quantitative determination of antibodies to thyroglobulin in human serum and plasma. The anti-Tg determination is used as an aid in the detection of autoimmune thyroid diseases.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

Summary

Thyroglobulin (Tg) is produced in the thyroid gland and is a main component in the lumen of the thyroid follicle. In synergy with the enzyme thyroid-specific peroxidase (TPO), Tg has an essential function in the iodination of L-tyrosine and in the formation of the thyroid hormones T4 and T3.¹ Both Tg and TPO are potentially autoantigenic.^{2,3}

Elevated serum concentrations of antibodies against Tg (Tg-autoantibodies) are found in subjects with autoimmunity-based thyroiditis.^{2,3} High concentrations of anti-Tg together with anti-TPO are present in most patients with chronic lymphocytic-infiltrative thyroiditis (Hashimoto's disease).³ The frequency of thyroglobulin antibodies is approximately 50-80 % in subjects with autoimmune-thyroiditis, including Hashimoto's disease, and approximately 30-50 % in individuals with Graves' disease.^{3,4,5,6} The anti-Tg assay can also provide useful information for monitoring the course of Hashimoto's thyroiditis and for differential diagnosis.^{3,7} This includes cases of suspected autoimmune thyroiditis of unknown origin with negative anti-TPO test results,^{8,9} and to distinguish Hashimoto's thyroiditis from nontoxic nodular goiter or from other forms of thyroiditis.⁴

Anti-Tg has also been reported as a useful surrogate diagnostic marker for differentiated thyroid cancer when serum Tg is negative,¹⁰ and for ruling out interference by Tg autoantibodies when measuring serum Tg using a Tg test.^{11,12}

Although the sensitivity of the procedure can be increased by simultaneously determining additional thyroid antibodies (anti-TPO, anti-TSHR), a negative result does not definitively rule out the presence of an autoimmune disease. The antibody titer does not correlate with the clinical activity of the disease. Titers that are elevated initially can become negative if the disease persists for a longer period of time or if remission occurs. If antibodies reappear after remission, relapse is likely.

The Elecsys Anti-Tg assay uses human antigen and monoclonal human anti-Tg antibodies.¹³

Test principle

Competition principle. Total duration of assay: 18 minutes.

- 1st incubation: 6 µL of sample are incubated with biotinylated Tg and the antibodies of the sample bind the antigen.

- 2nd incubation: After addition of anti-Tg antibodies labeled with ruthenium complex^{a)} and streptavidin-coated microparticles, the immunocomplex produced becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the **cobas** link.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The **cobas e** pack is labeled as ATG.

- M Streptavidin-coated microparticles, 1 bottle, 14.1 mL:
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Tg-biotin, 1 bottle, 19.7 mL:
Biotinylated Tg (human) 0.200 mg/L; TRIS buffer 100 mmol/L, pH 7.0; preservative.
- R2 Anti-Tg-Ab~Ru(bpy)₃²⁺, 1 bottle, 19.7 mL:
Monoclonal anti-Tg antibodies (human) labeled with ruthenium complex 0.620 mg/L; TRIS buffer 100 mmol/L, pH 7.0; preservative.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

- H317 May cause an allergic skin reaction.

Prevention:

- P261 Avoid breathing mist or vapours.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves.

Response:

- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

- P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

Elecsys Anti-Tg

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods use assays that have been approved or cleared by the FDA or that are in compliance with the legal rules of the European Union (IVDR 2017/746/EU, IVDD 98/79/EC, Annex II, List A). However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.^{14,15}

The initial thyroid glandular tissue extract containing the human thyroglobulin has shown to be free from HBsAg and antibodies to HCV and HIV.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is available via the **cobas** link.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the **cobas e** pack **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
on the analyzers	16 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

K₂-EDTA and K₃-EDTA plasma.

Criterion: Slope 0.9-1.1 + intercept within ± 20 IU/mL + coefficient of correlation ≥ 0.95 .

Stable for 4 days at 20-25 °C, 4 days at 2-8 °C, 2 months at -20 °C (± 5 °C). Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples and calibrators are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples and calibrators on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 09005030190, Anti-Tg CalSet, for 4 x 1.5 mL
- [REF] 05042666191, PreciControl ThyroAB, for 4 x 2.0 mL
- General laboratory equipment
- **cobas e** analyzer

Additional materials for **cobas e** 402 and **cobas e** 801 analyzers:

- [REF] 06908799190, ProCell II M, 2 x 2 L system solution
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 07485409001, Reservoir Cup, 8 cups to supply ProCell II M and CleanCell M
- [REF] 06908853190, PreClean II M, 2 x 2 L wash solution
- [REF] 05694302001, Assay Tip/Assay Cup tray, 6 magazines x 6 magazine stacks x 105 assay tips and 105 assay cups, 3 wasteliners
- [REF] 07485425001, Liquid Flow Cleaning Cup, 2 adaptor cups to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning Detection Unit
- [REF] 07485433001, PreWash Liquid Flow Cleaning Cup, 1 adaptor cup to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning PreWash Unit
- [REF] 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use.

Place the cooled (stored at 2-8 °C) **cobas e** pack on the reagent manager. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the **cobas e** pack.

Calibration

Traceability: This method has been standardized against the NIBSC (National Institute for Biological Standards and Control) 65/93 Standard.

The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the **cobas e** pack was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 28 days when using the same **cobas e** pack on the analyzer
- as required: e.g. quality control findings outside the defined limits

Quality control

Use Elecsys PreciControl ThyroAB or other suitable controls for routine quality control procedures.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per **cobas e** pack, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in IU/mL or kIU/L).

Limitations - interference

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations and no impact on results was observed.

Endogenous substances

Compound	Concentration tested
Bilirubin	≤ 1129 $\mu\text{mol/L}$ or ≤ 66 mg/dL

Elecsys Anti-Tg

Compound	Concentration tested
Hemoglobin	≤ 0.373 mmol/L or ≤ 600 mg/dL
Intralipid	≤ 2000 mg/dL
Biotin	≤ 4912 nmol/L or ≤ 1200 ng/mL
Rheumatoid factors	≤ 300 IU/mL

Criterion: For concentrations of 10-75 IU/mL the deviation is ≤ 11 IU/mL. For concentrations > 75 IU/mL the deviation is ≤ 15 %.

For samples ≤ 115 IU/mL no interference was observed for hemoglobin concentrations ≤ 600 mg/dL. In samples with a concentration of > 115 IU/mL a lower hemoglobin concentration may result in increased anti-Tg values.

Pharmaceutical substances

In vitro tests were performed on 17 commonly used pharmaceuticals. No interference with the assay was found.

In addition, the following special drugs were tested. No interference with the assay was found.

Special drugs

Drug	Concentration tested mg/L
Iodide	50
Carbimazole	30
Methimazole	16
Propylthiouracil	180
Perchlorate	2000
Propranolol	48
Amiodarone	40
Prednisolone	100
Hydrocortisone	200
Fluocortolone	100
Octreotide	0.300
Levothyroxine	0.250
Liothyronine	0.045
Nivolumab	96

Drug interferences are measured based on recommendations given in CLSI guidelines EP07 and EP37 and other published literature. Effects of concentrations exceeding these recommendations have not been characterized.

For Tg concentrations exceeding the normal range (> 100 ng/mL) an influence on anti-Tg concentrations of more than 15 % may occur.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

10-4000 IU/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 10 IU/mL. Values above the measuring range are reported as > 4000 IU/mL.

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank = 9 IU/mL

Limit of Detection = 10 IU/mL

Limit of Quantitation = 15 IU/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of ≤ 20 %.

Dilution

Sample dilution is not possible. The autoantibodies are heterogeneous and this gives rise to non-linear dilution phenomena.

Approximately 5 % of the pathological samples can have concentrations ≥ 4000 IU/mL.

Expected values

Studies conducted with the Elecsys Anti-Tg assay in 5 clinical centers covering a total of 391 healthy subjects (MCE Elecsys Anti-Tg assay) confirmed the threshold value of 115 IU/mL; this value corresponds to the 94th percentile.

For detailed information about reference intervals in children, adolescents and pregnant women, refer to the brochure "Reference Intervals for Children and Adults", English: [REF](#) 04640292.

This booklet also contains results of a detailed study about influencing factors on thyroid parameters in a well characterized reference group of adults. Different inclusion and exclusion criteria were applied (e.g. sonographic results (thyroid volume and density) as well as criteria according to the guidelines of the National Academy of Clinical Biochemistry - NACB).

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, samples and controls in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days ($n = 84$). The following results were obtained:

cobas e 402 and cobas e 801 analyzers					
Sample	Mean IU/mL	Repeatability		Intermediate precision	
		SD IU/mL	CV %	SD IU/mL	CV %
Human serum 1	14.2	0.639	4.5	1.12	7.9
Human serum 2	114	2.67	2.3	3.50	3.1
Human serum 3	1676	37.7	2.3	50.2	3.0
Human serum 4	1985	46.8	2.4	58.6	3.0
Human serum 5	3378	90.1	2.7	110	3.3
PC ^{b)} THYRO1	63.6	1.58	2.5	2.60	4.1
PC THYRO2	148	3.72	2.5	4.68	3.2

b) PC = PreciControl

Method comparison

a) A comparison of the Elecsys Anti-Tg assay, [REF](#) 09005021190 (cobas e 801 analyzer; y), with the Elecsys Anti-Tg assay, [REF](#) 07026919190 (cobas e 801 analyzer; x), gave the following correlations (IU/mL):

Number of samples measured: 150

Elecsys Anti-Tg

Passing/Bablok¹⁶ Linear regression
 $y = 1.01x - 5.04$ $y = 1.03x - 27.9$
 $\tau = 0.949$ $r = 0.997$

The sample concentrations were between 10.3 and 3785 IU/mL.

b) A comparison of the Elecsys Anti-Tg assay, [REF] 09005021190 (cobas e 402 analyzer; y), with the Elecsys Anti-Tg assay, [REF] 09005021190 (cobas e 801 analyzer; x), gave the following correlations (IU/mL):

Number of samples measured: 161

Passing/Bablok¹⁶ Linear regression
 $y = 1.03x - 0.385$ $y = 1.05x - 7.00$
 $\tau = 0.972$ $r = 0.999$

The sample concentrations were between 10.0 and 3675 IU/mL.

Analytical specificity

The following cross-reactivities were tested with anti-Tg concentrations of approximately 30 IU/mL and 115 IU/mL.

No influence with human autoantibodies to thyroid peroxidase (< 1500 IU/mL) was detectable.

References

- Mansourian AR. Metabolic pathways of tetraiodothyronine and triiodothyronine production by thyroid gland: a review of articles. Pak J Biol Sci 2011;14(1):1-12.
- Ruf J, Ferrand M, Durand-Gorde JM, et al. Significance of thyroglobulin antibodies cross-reactive with thyroperoxidase (TGPO antibodies) in individual patients and immunized mice. Clin Exp Immunol 1993;92(1):65-72.
- Thomas L. Thyroid function. Thyroglobulin antibodies. In: Thomas L (ed.). Deutsch: Labor und Diagnose. TH-Books, Frankfurt. 5th edition 1998:1043. English: Clinical Laboratory Diagnosis. 1st edition 1998:1021.
- Slatosky J, Shipton B, Wahba H. Thyroiditis: differential diagnosis and management. Am Fam Physician 2000;61(4):1047-1052.
- Garber JR, Cobin RH, Gharib H, et al. Clinical practice guidelines for hypothyroidism in adults: cosponsored by the American Association of Clinical Endocrinologists and The American Thyroid Association. Thyroid 2012;22(12):1200-1235.
- Iddah MA, Macharia BN. Autoimmune thyroid disorders. ISRN Endocrinol 2013:509764.
- Schmeltz LR, Blevins TC, Aronoff SL, et al. Anatabine supplementation decreases thyroglobulin antibodies in patients with chronic lymphocytic autoimmune (Hashimoto's) thyroiditis: A randomized controlled clinical trial. J Clin Endocrinol Metab 2014;99:E137-E142.
- Feldt-Rasmussen U. Analytical and clinical performance goals for testing autoantibodies to thyroperoxidase, thyroglobulin, and thyrotropin receptor. Clin Chem 1996;42(1):160-163.
- Lazarus J, Brown RS, Daumerie C, et al. 2014 European Thyroid Association guidelines for the management of subclinical hypothyroidism in pregnancy and in children. Eur Thyroid J 2014;3:76-94.
- Nam HY, Paeng JC, Chung JK, et al. Monitoring differentiated thyroid cancer patients with negative serum thyroglobulin. Diagnostic implication of TSH-stimulated antithyroglobulin antibody. Nuklearmedizin 2014;53(2):32-38.
- Spencer CA, Takeuchi M, Kazarosyan M, et al. Serum Thyroglobulin Antibodies: Prevalence, Influence on Serum Thyroglobulin Measurement, and Prognostic Significance in Patients with Differentiated Thyroid Carcinoma. J Clin Endocrin Metabol 1998;83(4):1121-1127.
- Spencer C. International Thyroid Testing Guidelines. National Academy of Clinical Biochemistry, August 2001;Section 3E,11-14.
- Prentice L, Kiso Y, Fukuma N, et al. Monoclonal Thyroglobulin Autoantibodies: Variable Region Analysis and Epitope Recognition. J Clin Endocrin Metabol 1995;80:977.

- Occupational Safety and Health Standards: Bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see navifyportal.roche.com for definition of symbols used):

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number
Rx only	For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2023, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com
 +800 5505 6606

