


**Alanine Aminotransferase acc. to IFCC II****Order information**

REF		CONTENT		Analyzer(s) on which <b>cobas c</b> pack(s) can be used
08104697190*	08104697500	Alanine Aminotransferase acc. to IFCC II (800 tests)	System-ID 2014 001	<b>cobas c 303, cobas c 503, cobas c 703</b>
08104697214*	08104697500	Alanine Aminotransferase acc. to IFCC II (800 tests)	System-ID 2014 001	<b>cobas c 303, cobas c 503, cobas c 703</b>

Materials required (but not provided):

10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 20401	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 20391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 20391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 20392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 20392	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

\* Some kits shown may not be available in all countries.

**English****System information****ALTP2: ACN 20140****Intended use**

In vitro test for the quantitative determination of alanine aminotransferase (ALT) with pyridoxal phosphate activation in human serum and plasma on **cobas c** systems.

**Summary**

Alanine aminotransferase (ALT) measurements, performed with this device, in human serum and plasma are used as an aid in diagnosis of hepatocellular injury and in monitoring chronic liver injury.

The enzyme alanine aminotransferase (ALT) is present in highest concentrations in the liver, in the cytosol of the hepatocytes, although it is also found in the kidney, and, in much smaller quantities, in heart and skeletal muscle cells.<sup>1</sup> ALT catalyzes the transfer of amino groups from L-alanine to  $\alpha$ -ketoglutarate, resulting in L-glutamate and pyruvate. This is a critical process of the tricarboxylic acid cycle, in which the coenzyme pyridoxal phosphate (also known as pyridoxal-5-phosphate or active vitamin B6) is required. When liver injury occurs, ALT is released from injured liver cells and causes a significant serum elevation.<sup>1</sup>

Measurement of ALT activity is therefore used for the diagnosis of hepatic diseases such as acute and chronic viral hepatitis, nonalcoholic fatty liver disease (NAFLD), alcohol-related liver disease, ischemic hepatopathy, autoimmune hepatitis, biliary injury, suspected malignant infiltration, cholestasis.<sup>1</sup> Serum elevations of ALT activity are rarely observed in conditions other than parenchymal liver disease.<sup>2</sup> In addition, ALT testing is recommended for monitoring chronic hepatitis status and progression.<sup>3</sup>

Although both serum aspartate aminotransferase (AST) and ALT become elevated whenever disease processes affect liver cell integrity, evidence suggests that ALT is a more specific marker of hepatic injury than AST. Moreover, elevations of ALT activity persist longer than elevations of AST activity.<sup>1,4</sup>

In patients with vitamin B6 deficiency (insufficient endogenous pyridoxal phosphate), serum aminotransferase activity may be decreased. The addition of pyridoxal phosphate to this assay causes an increase in aminotransferase activity (activation higher for AST than for ALT) and prevents falsely low aminotransferase test results in these samples.<sup>2</sup>

**Test principle**

This assay follows the recommendations of the IFCC, but was optimized for performance and stability.<sup>5</sup>

ALT catalyzes the transfer of an amino group between L-alanine and 2-oxoglutarate to form pyruvate and L-glutamate. The pyruvate then reacts with NADH in the presence of lactate dehydrogenase (LDH) to form L-lactate and NAD<sup>+</sup>. Pyridoxal phosphate serves as a coenzyme in the amino transfer reaction. It ensures full enzyme activation.



The rate of the NADH oxidation is directly proportional to the catalytic ALT activity. It is determined by measuring the decrease in absorbance.

**Reagents - working solutions**

**R1** TRIS buffer: 230 mmol/L, pH 7.15 (37 °C); L-alanine: 1140 mmol/L; LDH (microorganisms):  $\geq 94 \mu\text{kat/L}$ ; pyridoxamine phosphate: 0.23 mmol/L; albumin (bovine): 0.25 %; stabilizers; preservative

**R3** NADH:  $\geq 0.71 \text{ mmol/L}$ ; 2-oxoglutarate: 96 mmol/L; preservative

R1 is in position B and R3 is in position C.

**Precautions and warnings**

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

**Reagent handling**

Ready for use

**Storage and stability**

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer:

12 weeks

**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin and K<sub>2</sub>- and K<sub>3</sub>-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

**Alanine Aminotransferase acc. to IFCC II**

Stability: 4 days at 15-25 °C  
7 days at 2-8 °C

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section  
General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****Test definition**

Reporting time	10 min		
Wavelength (sub/main)	700/340 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	52 µL	48 µL	
R3	15 µL	–	
<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	4.5 µL	–	–
Decreased	4.5 µL	10 µL	90 µL
Increased	4.5 µL	–	–

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

**Calibration**

Calibrators	S1: H <sub>2</sub> O
	S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	Automatic full calibration
	- after reagent lot change
	Full calibration
	- as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the original IFCC formulation using calibrated pipettes together with a manual photometer providing absolute values and the substrate-specific absorptivity,  $\epsilon$ .<sup>5</sup>

**Quality control**

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 12 weeks. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

**cobas c** systems automatically calculate the analyte activity of each sample in the unit U/L (µkat/L).

Conversion factor: U/L x 0.0167 = µkat/L

**Limitations - interference**

Criterion: Recovery within  $\pm 4.0$  U/L of initial values of samples  $\leq 40$  U/L and  $\pm 10$  % for samples  $> 40$  U/L.

Icterus:<sup>6</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:<sup>6</sup> No significant interference up to an H index of 100 (approximate hemoglobin concentration: 62.2 µmol/L or 100 mg/dL). Contamination with erythrocytes will elevate results, because the analyte level in erythrocytes is higher than in normal sera. The level of interference may be variable depending on the content of analyte in the lysed erythrocytes.

Lipemia (Intralipid):<sup>6</sup> No significant interference up to an L index of 500. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Lipemic samples may cause > Abs flagging.

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>7,8</sup>

Drug interferences are measured based on recommendations given in CLSI guidelines EP07 and EP37 and other published literature. Effects of concentrations exceeding these recommendations have not been characterized.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>9</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

**Limits and ranges****Measuring range**

5-700 U/L (0.08-11.7 µkat/L)

Determine samples having higher activities via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 10.

**Lower limits of measurement***Limit of Blank, Limit of Detection and Limit of Quantitation*

Limit of Blank	= 5 U/L (0.08 µkat/L)
Limit of Detection	= 5 U/L (0.08 µkat/L)
Limit of Quantitation	= 5 U/L (0.08 µkat/L)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the activity below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low activity samples.

The Limit of Detection corresponds to the lowest analyte activity which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte activity that can be reproducibly measured with a total error of 20 %. It has been determined using low activity alanine aminotransferase samples.

**Expected values****U/L**

Acc. to IFCC/Standard Method 94 with pyridoxal phosphate activation measured at 37 °C:<sup>10</sup>

Males: 10-50 U/L

Females: 10-35 U/L

Consensus values with pyridoxal phosphate activation:<sup>11</sup>

Males: up to 50 U/L

Females: up to 35 U/L

**µkat/L\***

Acc. to IFCC/Standard Method 94 with pyridoxal phosphate activation measured at 37 °C:<sup>10</sup>

Males: 0.17-0.84 µkat/L

Females: 0.17-0.58 µkat/L

Consensus values with pyridoxal phosphate activation:<sup>11</sup>

Males: up to 0.84 µkat/L

Females: up to 0.58 µkat/L

\*calculated by unit conversion factor

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

**Precision**

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c 503** analyzer.

<i>Repeatability</i>	<i>Mean</i> U/L	<i>SD</i> U/L	<i>CV</i> %
PCCC <sup>1a)</sup>	49.4	0.534	1.1
PCCC <sup>2b)</sup>	127	3.04	2.4
Human serum 1	12.8	0.474	3.7
Human serum 2	30.8	0.601	2.0
Human serum 3	54.7	0.965	1.8
Human serum 4	359	2.45	0.7
Human serum 5	630	2.81	0.4
<i>Intermediate precision</i>	<i>Mean</i> U/L	<i>SD</i> U/L	<i>CV</i> %
PCCC <sup>1a)</sup>	49.2	1.80	3.7
PCCC <sup>2b)</sup>	127	4.96	3.9
Human serum 1	12.8	0.611	4.8
Human serum 2	30.8	0.818	2.7
Human serum 3	54.7	1.58	2.9
Human serum 4	359	3.28	0.9
Human serum 5	638	5.07	0.8

a) PreciControl ClinChem Multi 1

b) PreciControl ClinChem Multi 2

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s) and **cobas c 703** analyzer(s).

**Method comparison**

ALT values for human serum and plasma samples obtained on a **cobas c 503** analyzer (y) were compared with those determined using the test ALTLP on a **cobas c 501** analyzer (x).

Sample size (n) = 100

Passing/Bablok <sup>12</sup>	Linear regression
$y = 0.993x + 1.52$ U/L	$y = 0.988x + 1.71$ U/L
$\tau = 0.988$	$r = 1.000$

The sample activities were between 8.9 and 683 U/L.

ALT values for human serum and plasma samples obtained on a **cobas c 303** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 503** analyzer (x).

Sample size (n) = 50

Passing/Bablok <sup>12</sup>	Linear regression
$y = 1.036x - 0.787$ U/L	$y = 1.039x - 1.61$ U/L
$\tau = 0.997$	$r = 1.000$

The sample activities were between 27.0 and 635 U/L.

ALT values for human serum and plasma samples obtained on a **cobas c 703** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 503** analyzer (x).

Sample size (n) = 65

Passing/Bablok <sup>12</sup>	Linear regression
$y = 1.004x - 0.634$ U/L	$y = 1.003x - 0.635$ U/L
$\tau = 0.974$	$r = 1.000$

The sample concentrations were between 6.07 and 647 U/L.

**References**

- Kim WR, Flamm SL, Di Bisceglie AM, et al. Public Policy Committee of the American Association for the Study of Liver Disease. Serum activity of alanine aminotransferase (ALT) as an indicator of health and disease. *Hepatology* 2008 Apr;47(4):1363-1370. doi: 10.1002/hep.22109.
- Panteghini M, Bais R. Amino Transferases. In: Burtis CA, Ashwood ER, Bruns DE. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 5th ed. 2012;573-576.
- Guidelines for the Prevention, Care and Treatment of Persons with Chronic Hepatitis B Infection. Geneva: World Health Organization; 2015 Mar.
- Kwo PY, Cohen SM, Lim JK. ACG Clinical Guideline: Evaluation of Abnormal Liver Chemistries. *Am J Gastroenterol* 2017 Jan;112(1):18-35. doi: 10.1038/ajg.2016.517.
- Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 ° – Part 4. Reference Procedure for the Measurement of Catalytic Activity Concentration of Alanine Aminotransferase. *Clin Chem Lab Med* 2002;40(7):718-724. doi: 10.1515/CCLM.2002.124.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.

# ALTP2

## Alanine Aminotransferase acc. to IFCC II



- 10 Klauke R, Schmidt E, Lorentz K. Recommendations for carrying out standard ECCLS procedures (1988) for the catalytic concentrations of creatine kinase, aspartate aminotransferase, alanine aminotransferase and  $\gamma$ -glutamyltransferase at 37 °C. Eur J Clin Chem Clin Biochem 1993;31:901-909.
- 11 Thomas L, Müller M, Schumann G, et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29(5):301-308.
- 12 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [navifyportal.roche.com](http://navifyportal.roche.com) for definition of symbols used):

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number
Rx only	For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)  
 +800 5505 6606

