## Tina-quant Complement C3c ver.2



#### Order information

08105537500V7 (

REF	Ĩ			Analyzer(s) on which <b>cobas c</b> pack(s) can be used
08105537190*	08105537500	Tina-quant Complement C3c ver.2 (150 tests)	System-ID 2032 001	<b>cobas c</b> 303, <b>cobas c</b> 503, <b>cobas c</b> 703
08105537214*	08105537500	Tina-quant Complement C3c ver.2 (150 tests)	System-ID 2032 001	<b>cobas c</b> 303, <b>cobas c</b> 503, <b>cobas c</b> 703

#### Materials required (but not provided):

11355279216	Calibrator f.a.s. Proteins (5 × 1 mL)	Code 20656	
10557897122	Precinorm Protein (3 × 1 mL)	Code 20302	
11333127122	Precipath Protein (3 × 1 mL)	Code 20303	
05117003190	PreciControl ClinChem Multi 1 (20 × 5 mL)	Code 20391	
05947626190	PreciControl ClinChem Multi 1 (4 × 5 mL)	Code 20391	
05117216190	PreciControl ClinChem Multi 2 (20 × 5 mL)	Code 20392	
05947774190	PreciControl ClinChem Multi 2 (4 × 5 mL)	Code 20392	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

08063494190 | Diluent NaCl 9 % (123 mL) \* Some kits shown may not be available in all countries.

#### English

#### System information

C3C-2: ACN 20320

#### Intended use

In vitro test for the quantitative determination of Complement C3c in human serum and plasma on cobas c systems.

#### Summarv

Measurements of Complement C3c performed with this assay, in human serum and plasma, can be used as an aid in assessing possible complement system imbalance, associated with or observed during a number of underlying disease states or pathological conditions, including inflammatory and infectious diseases.

### Complement component C3 (C3), is one of the plasma and membrane-associated components of the human complement system.1 Complement components can be activated via 3 pathways - the Classical, Alternative, and Lectin pathways - all converging at C3 activation via proteolytic processing.<sup>1,2</sup> This cleavage event releases C3a, an anaphylatoxin, and C3b, an opsonin and necessary component of the C5 convertase, which in turn, is responsible for the assembly of the terminal pathway and formation of the membrane attack complex.<sup>1,2</sup> C3b is tightly regulated both functionally and by further proteolytic processing. C3c represents a stable, terminal cleavage product of C3b, and its generation requires Factor I, Complement Receptor 1 and other cofactors of Factor I such as Factor H and Membrane Cofactor Protein.<sup>1,2</sup> Therefore, the concentration of C3c can be evaluated as a parameter for diagnostic workup or for monitoring of complement system activation.<sup>3,4</sup> Measurement of C3c is most clinically valuable when performed alongside the measurement of other complement proteins and other diagnostic tests.<sup>3,4</sup> Increased C3c levels are indicative of complement activation, 3,4 however a normal C3c level on its own should not exclude C3 activation.

Complement activation is observed in a number of autoimmune, inflammatory conditions and infectious diseases,<sup>4,5</sup> for example systemic lupus erythematosus (SLE), rheumatoid arthritis, subacute bacterial endocarditis, viremia, parasitic infections and bacterial sepsis. 4,5,6,7,8,9,10,11,12 Considerable C3 activation can also be found in patients with acquired partial lipodystrophy and C3 glomerulopathy, and C3c levels may be elevated.  $^{13,14,15,16}$ 

#### Test principle<sup>17</sup>

#### Immunoturbidimetric assay

Human C3c forms a precipitate with a specific antiserum which is determined turbidimetrically.

#### **Reagents - working solutions**

R1 TRIS buffer: 100 mmol/L, pH 8.0; polyethylene glycol: 3.0 %; preservative

**R**3 Anti-human C3c antibody (goat): dependent on titer; TRIS buffer: 33 mmol/L; preservative

R1 is in position B and R3 is in position C.

#### Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste: Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures. Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal. Safety data sheet available for professional user on request.

#### Reagent handling

Ready for use

#### Storage and stability

Shelf life at 2-8 °C:	See expiration date on <b>cobas c</b> pack label.
On board in use and refrigerated on the	22 wooks

On-board in use and refrigerated on the 22 weeks analyzer:

#### Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable. Serum

Plasma: Li-heparin plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

Stability:18 4 days at 20-25 °C 8 days at 4-8 °C 8 days at -20 °C (± 5 °C)

Freeze only once.

## 08105537500V7 ( Tina-guant Complement C3c ver.2



The degree of fragmentation of C3 to C3c ex vivo depends on the age and storage conditions of the sample. For fresh samples the values obtained are found to be up to 25 % lower than those obtained for aged samples depending on the extent to which fragmentation has occurred.<sup>19</sup>

#### Materials provided

See "Reagents – working solutions" section for reagents.

#### Materials required (but not provided)

See "Order information" section

General laboratory equipment

#### Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

#### Application for serum and plasma

#### **Test definition**

Reporting time	10 min		
Wavelength (sub/main)	700/340 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	75 µL	-	
R3	14 µL	17 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	8.3 μL	5 µL	100 µL
Decreased	8.3 μL	2 µL	82 µL
Increased	8.3 µL	5 µL	100 µL

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

#### Calibration

Calibrators	S1: H <sub>2</sub> O
	S2: C.f.a.s. Proteins
Calibration mode	Non-linear
Calibration frequency	Automatic full calibration - after reagent lot change
	Full calibration - as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the reference preparation of the IRMM (Institute for Reference Materials and Measurements) BCR470/CRM470 (RPPHS - Reference Preparation for Proteins in Human Serum).20

#### Quality control

For guality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 22 weeks. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits

Follow the applicable government regulations and local guidelines for guality control.

#### Calculation

cobas c systems automatically calculate the analyte concentration of each sample in the unit g/L (mg/dL).

Conversion factor:  $g/L \times 100 = mg/dL$ 

#### Limitations - interference

Criterion: Recovery within ± 10 % of initial values at a C3c level of 0.9 g/L.

Icterus:<sup>21</sup> No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:<sup>21</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 621 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):<sup>21</sup> No significant interference up to an L index of 2000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 1200 IU/mL.

High-dose hook effect: No false result occurs up to a C3c concentration of 12.5 g/L.

Drugs: No interference was found at therapeutic concentrations using common drug panels.22,23

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>24</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

#### ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on cobas c systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

#### Limits and ranges

#### Measuring range

0.04-5.0 g/L (4-500 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:2 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 2.

#### Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank	= 0.04 g/L (4 mg/dL)
Limit of Detection	= 0.04 g/L (4 mg/dL)
Limit of Quantitation	= 0.05 g/L (5 mg/dL)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \ge 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration C3c samples.

#### Expected values<sup>25</sup>

0.9-1.8 g/L (90-180 mg/dL\*) \*calculated by unit conversion facto



Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

#### Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

#### Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c** 503 analyzer.

Repeatability	Mean	SD	CV
	g/L	g/L	%
PCCC1 <sup>a)</sup>	0.904	0.00711	0.8
PCCC2 <sup>b)</sup>	1.55	0.0111	0.7
Human serum 1	0.101	0.00328	3.2
Human serum 2	0.879	0.00697	0.8
Human serum 3	1.73	0.0104	0.6
Human serum 4	2.31	0.0285	1.2
Human serum 5	4.12	0.0316	0.8
Intermediate precision	Mean a/L	SD a/L	CV %
PCCC1 <sup>a)</sup>	0.880	0.00951	1.1
PCCC2 <sup>b)</sup>	1.55	0.0174	1.1
Human serum 1	0.101	0.00461	4.6
Human serum 2	0.838	0.00991	1.2
Human serum 3	1.79	0.0168	0.9
Human serum 4	2.31	0.0313	1.4
Human serum 5	4 15	0.0476	1.1

a) PreciControl ClinChem Multi 1

b) PreciControl ClinChem Multi 2

The data obtained on **cobas c** 503 analyzer(s) are representative for **cobas c** 303 analyzer(s) and **cobas c** 703 analyzer(s).

#### Method comparison

C3c values for human serum and plasma samples obtained on a **cobas c** 503 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 69

Passing/Bablok <sup>26</sup>	Linear regression
y = 1.013x + 0.0232 g/L	y = 0.983x + 0.0530 g/L
т = 0.981	r = 1.000

The sample concentrations were between 0.0560 and 4.94 g/L. C3c values for human serum and plasma samples obtained on a **cobas c** 303 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 501 analyzer (x).

Sample size (n) = 68

Passing/Bablok <sup>26</sup>	Linear regression
y = 0.983x + 0.0332 g/L	y = 0.959x + 0.0633 g/L
т = 0.983	r = 0.999

The sample concentrations were between 0.0410 and 4.93 g/L.



C3c values for human serum and plasma samples obtained on a **cobas c** 703 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 503 analyzer (x).

Sample size (n) = 68

Passing/Bablok <sup>26</sup>	Linear regression
y = 0.999x - 0.0241 g/L	y = 1.000x - 0.0237 g/L
т = 0.987	r = 0.999

The sample concentrations were between 0.0464 and 4.88 g/L.

#### References

- 1 Dunkelberger JR, Song WC. Complement and its role in innate and adaptive immune responses. Cell Res 2010 Jan;20(1):34-50.
- 2 Mathern DR, Heeger PS. Molecules Great and Small: The Complement System. Clin J Am Soc Nephrol 2015 Sep 4;10(9):1636-1650.
- 3 Nilsson B, Ekdahl KN. Complement diagnostics: concepts, indications, and practical guidelines. Clin Dev Immunol 2012;2012:962702.
- 4 Ekdahl KN, Persson B, Mohlin C, et al. Interpretation of Serological Complement Biomarkers in Disease. Front Immunol 2018 Oct 24;9:2237.
- Hebert LA, Cosio FG, Neff JC. Diagnostic significance of hypocomplementemia. Kidney Int 1991 May;39(5):811-821.
- 6 Petri M, Orbai AM, Alarcón GS, et al. Derivation and validation of the Systemic Lupus International Collaborating Clinics classification criteria for systemic lupus erythematosus. Arthritis Rheum 2012 Aug;64(8):2677-2686.
- 7 Keeling SO, Alabdurubalnabi Z, Avina-Zubieta A, et al. Canadian Rheumatology Association Recommendations for the Assessment and Monitoring of Systemic Lupus Erythematosus. J Rheumatol 2018 Oct;45(10):1426-1439.
- 8 Aringer M, Costenbader K, Daikh D, et al. 2019 European League Against Rheumatism/American College of Rheumatology classification criteria for systemic lupus erythematosus. Ann Rheum Dis 2019 Sep;78(9):1151-1159.
- 9 Li W, Li H, Song W, et al. Differential diagnosis of systemic lupus erythematosus and rheumatoid arthritis with complements C3 and C4 and C-reactive protein. Exp Ther Med 2013 Nov;6(5):1271-1276.
- 10 Messias-Reason IJ, Hayashi SY, Nisihara RM, et al. Complement activation in infective endocarditis: correlation with extracardiac manifestations and prognosis. Clin Exp Immunol 2002 Feb;127(2):310-315.
- 11 El-Fatah Fahmy Hanno A, Mohiedeen KM, Deghedy A, et al. Serum complements C3 and C4 in chronic HCV infection and their correlation with response to pegylated interferon and ribavirin treatment. Arab J Gastroenterol 2014 Jun;15(2):58-62.
- 12 Ren J, Zhao Y, Yuan Y, et al. Complement depletion deteriorates clinical outcomes of severe abdominal sepsis: a conspirator of infection and coagulopathy in crime? PLoS One 2012;7(10):e47095.
- 13 Mastrangelo A, Serafinelli J, Giani M, et al. Clinical and Pathophysiological Insights Into Immunological Mediated Glomerular Diseases in Childhood. Front Pediatr 2020 May 12;8:205.
- 14 Sethi S, Fervenza FC, Zhang Y et al. C3 glomerulonephritis: clinicopathological findings, complement abnormalities, glomerular proteomic profile, treatment, and follow-up. Kidney Int 2012 Aug;82(4):465-473.
- 15 Sethi S, Nester CM, Smith RJ. Membranoproliferative glomerulonephritis and C3 glomerulopathy: resolving the confusion. Kidney Int 2012 Mar;81(5):434-441.
- 16 Corvillo F, Akinci B. An overview of lipodystrophy and the role of the complement system. Mol Immunol 2019 Aug;112:223-232.
- 17 Müller-Eberhard HJ. Complement: Chemistry and pathways. In: Inflammation: Basic principles and clinical correlates. Gallin I, Goldstein IM, Snyderman R, eds. New York: Raven Press 1988;21-53.
- 18 Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.

08105537500V7.0

#### Tina-quant Complement C3c ver.2

- 19 Okumura N, Nomura M, Tada T, et al. Effects of sample storage on serum C3C assay by nephelometry. Clin Lab Sci 1990;3:54-57.
- 20 Baudner S, Bienvenu J, Blirup-Jensen S, et al. The certification of a matrix reference material for immunochemical measurement of 14 human serum proteins CRM470. Report EUR 15243 EN 1993;1-186.
- 21 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 22 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 23 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 24 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 25 Dati F, Schumann G, Thomas L, et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-520.
- 26 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

#### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

CONTENT	Contents of kit
$\longrightarrow$	Volume for reconstitution
GTIN	Global Trade Item Number
Rx only	For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM, PRECIPATH and TINA-QUANTare trademarks of Roche. All other product names and trademarks are the property of their respective owners. Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH	
Sandhofer Strasse 116	
68305 Mannheim, Germany	
www.roche.com	
+800 5505 6606	



# cobas®