

REF		CONTENT		Analyzer(s) on which <b>cobas c</b> pack(s) can be used
08057532190*	08057532500	Creatinine Jaffé Gen.2 (2500 tests)	System-ID 2047 001	<b>cobas c 303, cobas c 503, cobas c 703</b>
08057532214*	08057532500	Creatinine Jaffé Gen.2 (2500 tests)	System-ID 2047 001	<b>cobas c 303, cobas c 503, cobas c 703</b>

Materials required (but not provided):

10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 20401	
03121313122	Precinorm PUC (4 x 3 mL)	Code 20240	
03121291122	Precipath PUC (4 x 3 mL)	Code 20241	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 20391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 20391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 20392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 20392	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

\* Some kits shown may not be available in all countries.

**English****System information****CREJ2:** ACN 20470 (Serum/plasma)**CREJ2U:** ACN 20471 (Urine)**Intended use**In vitro test for the quantitative determination of creatinine in human serum, plasma and urine on **cobas c** systems.**Summary**

Creatinine measurements, performed with this assay, in human serum, plasma and urine are used as an aid in diagnosis and monitoring of renal disease and in monitoring of renal dialysis. Creatinine measurements are also used for the calculation of the fractional excretion of other urine analytes (e. g., albumin,  $\alpha$ -amylase).

Creatinine is a break-down product of creatine phosphate in muscle, and is usually produced at a fairly constant rate by the body (depending on muscle mass). It is freely filtered by the glomeruli and, under normal conditions, is not reabsorbed by the tubules to any appreciable extent. A small but significant amount is also actively secreted. Its concentration is thus, inversely related to glomerular filtration rate (GFR).<sup>1,2</sup>

The assay of creatinine in serum or plasma is the most commonly used test to assess renal function. Chronic kidney disease is a worldwide problem that carries a substantial risk for cardiovascular morbidity and death. Current guidelines define chronic kidney disease as kidney damage or decreased glomerular filtration rate (GFR) (less than 60 mL/min per 1.73 m<sup>2</sup>) for 3 months or more.<sup>2,3</sup>


Since a rise in blood creatinine is observed only with marked damage of the nephrons, it is not suited to detect early stage kidney disease. A considerably more sensitive test and better estimation of glomerular filtration rate (GFR) is given by the creatinine clearance test based on creatinine's concentration in urine and serum or plasma, and urine flow rate. For this test a precisely timed urine collection (usually 24 hours) and a blood sample are needed. However, since this test is prone to error due to the inconvenient collection of timed urine, mathematical attempts to estimate GFR (eGFR) based only on the creatinine concentration in serum or plasma have been made.<sup>4</sup> Among the various approaches suggested, three have found wide recognition: the Cockcroft and Gault, the Modification of Diet in Renal Disease (MDRD) Study equation and the CKD-EPI (Chronic Kidney Disease Epidemiology) equation. While the Cockcroft and Gault equation was derived from data in which serum creatinine was measured with the conventional Jaffé method, the MDRD study equation measured serum creatinine using the Jaffé method calibrated to an isotope dilution mass spectrometry (IDMS).<sup>5,6</sup> These estimates of GFR are useful during monitoring of renal dialysis.<sup>7,8</sup> In children, the Bedside Schwartz formula should be used.<sup>9,10,11</sup>

In addition to the diagnosis and treatment of renal disease and the monitoring of renal dialysis, creatinine measurements are used for the calculation of the fractional excretion of other urine analytes (e. g., albumin,

$\alpha$ -amylase). Numerous methods were described for determining creatinine. Automated assays established in the routine laboratory include the Jaffé alkaline picrate method in various modifications, as well as enzymatic tests.<sup>2</sup>

**Test principle**<sup>12,13,14</sup>

This kinetic colorimetric assay is based on the Jaffé method. In alkaline solution, creatinine forms a yellow-orange complex with picrate. The rate of dye formation is proportional to the creatinine concentration in the specimen. The assay uses "rate-blanking" to minimize interference by bilirubin. To correct for non-specific reaction caused by serum/plasma pseudo-creatinine chromogens, including proteins and ketones, the results for serum or plasma are corrected by -26  $\mu$ mol/L (-0.3 mg/dL).

**Reagents - working solutions**

**R1** Potassium hydroxide: 900 mmol/L; phosphate: 135 mmol/L; pH  $\geq$  13.5; preservative; stabilizer

**R3** Picric acid: 38 mmol/L; pH 6.5; non reactive buffer

R1 is in position B and R3 is in position C.

**Precautions and warnings**

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H314 Causes severe skin burns and eye damage.

**Prevention:**

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

3 weeks at -20 °C (± 5 °C)

**Response:**

- P301 + P330 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
+ P331
- P303 + P361 IF ON SKIN (or hair): Take off immediately all contaminated  
+ P353 clothing. Rinse skin with water.
- P304 + P340 IF INHALED: Remove person to fresh air and keep  
+ P310 comfortable for breathing.  
Immediately call a POISON CENTER/ doctor.
- P305 + P351 IF IN EYES: Rinse cautiously with water for several  
+ P338 minutes. Remove contact lenses, if present and easy to do.  
+ P310 Continue rinsing. Immediately call a POISON CENTER/  
doctor.

**Disposal:**

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

**Reagent handling**

Ready for use

**Storage and stability**

Shelf life at 15-25 °C: See expiration date on  
**cobas c** pack label.

On-board in use and refrigerated on the analyzer: 26 weeks

**Specimen collection and preparation<sup>15</sup>**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K<sub>2</sub>-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine.

Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14 to 47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used. If stabilizers are added to the sample, the sample index feature must not be used.

Stability in *serum/plasma*:<sup>16</sup> 7 days at 15-25 °C  
7 days at 2-8 °C  
3 months at -20 °C (± 5 °C)

Freeze only once.

Stability in *urine* (without preservative):<sup>16</sup> 2 days at 15-25 °C  
6 days at 2-8 °C  
6 months at -20 °C (± 5 °C)

Freeze only once.

Stability in *urine* (with preservative): 3 days at 15-25 °C  
8 days at 2-8 °C

Freeze only once.

Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****Test definition**

Reporting time	10 min		
Wavelength (sub/main)	570/505 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	10 µL	58 µL	
R3	13 µL	23 µL	

**Sample volumes**

	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	7.5 µL	–	–
Decreased	7.5 µL	20 µL	80 µL
Increased	7.5 µL	–	–

**Application for urine****Test definition**

Reporting time	10 min		
Wavelength (sub/main)	570/505 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	10 µL	58 µL	
R3	13 µL	23 µL	

**Sample volumes**

	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	7.5 µL	4 µL	96 µL
Decreased	7.5 µL	1.5 µL	135 µL
Increased	7.5 µL	4 µL	96 µL

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

**Calibration**

*Application for serum/plasma (ACN 20470)*

Calibrators S1: H<sub>2</sub>O  
S2: C.f.a.s.

Calibration mode Linear

Calibration frequency	Automatic full calibration - after reagent lot change
	Full calibration - every 8 weeks on-board - as required following quality control procedures

**Application for urine (ACN 20471)**

Transfer of calibration from serum/plasma application (ACN 20470)

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against ID/MS.

**Quality control**

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

Serum/plasma: PreciControl ClinChem Multi 1, PreciControl ClinChem Multi 2

Urine: Precinorm PUC, Precipath PUC

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 26 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

**cobas c** systems automatically calculate the analyte concentration of each sample in the unit  $\mu\text{mol/L}$  (mg/dL, mmol/L, mg/L).

Conversion factors:	$\mu\text{mol/L} \times 0.0113 = \text{mg/dL}$
	$\mu\text{mol/L} \times 0.001 = \text{mmol/L}$
	$\mu\text{mol/L} \times 0.113 = \text{mg/L}$

**Limitations – interference**

Criterion: Recovery within  $\pm 10\%$  of initial value at a creatinine concentration of  $80 \mu\text{mol/L}$  ( $0.90 \text{ mg/dL}$ ) in serum/plasma and  $2.5 \text{ mmol/L}$  ( $28.3 \text{ mg/dL}$ ) in urine.

**Serum/plasma**

Icterus (**CREJ2**):<sup>17</sup> No significant interference up to an I index of 5 for conjugated bilirubin and 10 for unconjugated bilirubin (approximate conjugated bilirubin concentration:  $86 \mu\text{mol/L}$  or  $5 \text{ mg/dL}$ ; approximate unconjugated bilirubin concentration:  $171 \mu\text{mol/L}$  or  $10 \text{ mg/dL}$ ).

Hemolysis:<sup>17</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration:  $621 \mu\text{mol/L}$  or  $1000 \text{ mg/dL}$ ).

Lipemia (Intralipid):<sup>17</sup> No significant interference up to an L index of 800. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Pyruvate: No significant interference from pyruvate up to a concentration of  $0.3 \text{ mmol/L}$  ( $2.6 \text{ mg/dL}$ ).

Glucose: No significant interference from glucose up to a concentration of  $25 \text{ mmol/L}$  ( $450 \text{ mg/dL}$ ).

Ascorbic acid: No significant interference from ascorbic acid up to a concentration of  $5 \text{ mmol/L}$  ( $88 \text{ mg/dL}$ ).

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>18,19</sup>

Exception: Antibiotics containing cephalosporin lead to significant false-positive values.<sup>20,21</sup> Cefoxitin causes artificially high creatinine results. Cyanokit (Hydroxocobalamin) may cause interference with results.

Values  $< 15 \mu\text{mol/L}$  ( $< 0.17 \text{ mg/dL}$ ) or negative results are reported in rare cases in children  $< 3$  years and in elderly patients. In such cases use the Creatinine plus test to assay the sample.

Do not use Creatinine Jaffé for the testing of creatinine in hemolyzed samples from neonates, infants or adults with HbF levels  $\geq 60 \text{ mg/dL}$  for

**CREJ2** applications.<sup>22</sup> In such cases, use the Creatinine plus test ( $\leq 600 \text{ mg/dL HbF}$ ) to assay the sample.

Estimation of the Glomerular Filtration Rate (GFR) on the basis of the Schwartz Formula can lead to an overestimation.<sup>23</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>24</sup>

The presence of ketone bodies can cause artificially high results in serum and plasma.

**Urine**

Icterus: No significant interference up to a conjugated bilirubin concentration of  $855 \mu\text{mol/L}$  or  $50 \text{ mg/dL}$ .

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration of  $621 \mu\text{mol/L}$  or  $1000 \text{ mg/dL}$ ).

Glucose: No significant interference from glucose up to a concentration of  $120 \text{ mmol/L}$  ( $2162 \text{ mg/dL}$ ).

Urea: No significant interference from urea up to a concentration of  $2100 \text{ mmol/L}$  ( $12612 \text{ mg/dL}$ ).

Urobilinogen: No significant interference from urobilinogen up to a concentration of  $676 \mu\text{mol/L}$  ( $40 \text{ mg/dL}$ ).

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>19</sup>

Exception: Cyanokit (Hydroxocobalamin) may cause interference with results.

High homotestic acid concentrations in urine samples lead to false results.

The presence of ketone bodies can cause artificially high results in urine.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

**Limits and ranges****Measuring range****Serum/plasma**

$15\text{--}2200 \mu\text{mol/L}$  ( $0.17\text{--}24.9 \text{ mg/dL}$ )

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 5.

**Urine**

$0.375\text{--}55 \text{ mmol/L}$  ( $4.2\text{--}622 \text{ mg/dL}$ )

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3.6 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 3.6.

**Lower limits of measurement****Limit of Blank, Limit of Detection and Limit of Quantitation****Serum/plasma (CREJ2)**

Limit of Blank =  $15 \mu\text{mol/L}$  ( $0.17 \text{ mg/dL}$ )

Limit of Detection =  $15 \mu\text{mol/L}$  ( $0.17 \text{ mg/dL}$ )

Limit of Quantitation =  $15 \mu\text{mol/L}$  ( $0.17 \text{ mg/dL}$ )

**Urine (CREJ2U)**

Limit of Blank =  $0.375 \text{ mmol/L}$  ( $4.24 \text{ mg/dL}$ )

Limit of Detection =  $0.375 \text{ mmol/L}$  ( $4.24 \text{ mg/dL}$ )

Limit of Quantitation =  $0.375 \text{ mmol/L}$  ( $4.24 \text{ mg/dL}$ )

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from  $n \geq 60$  measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration creatinine samples.

**Expected values****µmol/L***Serum/plasma**Adults<sup>25</sup>*

Females	44-80 µmol/L
Males	62-106 µmol/L

*Children<sup>26</sup>*

Neonates (premature)	25-91 µmol/L
Neonates (full term)	21-75 µmol/L
2-12 months	15-37 µmol/L
1- < 3 years	21-36 µmol/L
3- < 5 years	27-42 µmol/L
5- < 7 years	28-52 µmol/L
7- < 9 years	35-53 µmol/L
9- < 11 years	34-65 µmol/L
11- < 13 years	46-70 µmol/L
13- < 15 years	50-77 µmol/L

**mmol/L***Urine**1st morning urine<sup>25</sup>*

Females	2.47-19.2 mmol/L
Males	3.45-22.9 mmol/L

*24-hour urine<sup>27</sup>*

Females	7.0-14.0 mmol/24 h
Males	9.0-21.0 mmol/24 h

Creatinine clearance<sup>27,28</sup> 71-151 mL/min

Refer to reference for a prospective study on creatinine clearance in children.<sup>29</sup>

**mg/dL***Serum/plasma**Adults<sup>25</sup>*

Females	0.50-0.90 mg/dL
Males	0.70-1.20 mg/dL

*Children<sup>26</sup>*

Neonates (premature)	0.29-1.04 mg/dL
Neonates (full term)	0.24-0.85 mg/dL
2-12 months	0.17-0.42 mg/dL
1- < 3 years	0.24-0.41 mg/dL

3- < 5 years	0.31-0.47 mg/dL
5- < 7 years	0.32-0.59 mg/dL
7- < 9 years	0.40-0.60 mg/dL
9- < 11 years	0.39-0.73 mg/dL
11- < 13 years	0.53-0.79 mg/dL
13- < 15 years	0.57-0.87 mg/dL

*Urine**1st morning urine<sup>25</sup>*

Females	28-217 mg/dL
Males	39-259 mg/dL

*24-hour urine<sup>27</sup>*

Females	740-1570 mg/24 h
Males	1040-2350 mg/24 h

Creatinine clearance<sup>27,28</sup> 71-151 mL/min

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

**Precision**

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability ( $n = 84$ ) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c 503** analyzer.

*Serum/plasma (CREJ2)*

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>µmol/L</i>	<i>µmol/L</i>	<i>%</i>
PCCC1 <sup>a)</sup>	90.4	1.57	1.7
PCCC2 <sup>b)</sup>	347	3.87	1.1
Human serum 1	48.2	1.40	2.9
Human serum 2	71.8	1.51	2.1
Human serum 3	480	4.15	0.9
Human serum 4	1064	12.0	1.1
Human serum 5	1791	20.6	1.2

<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>µmol/L</i>	<i>µmol/L</i>	<i>%</i>
PCCC1 <sup>a)</sup>	89.2	2.33	2.6
PCCC2 <sup>b)</sup>	347	5.24	1.5
Human serum 1	48.2	1.64	3.4
Human serum 2	71.8	1.89	2.6
Human serum 3	480	7.59	1.6
Human serum 4	1064	16.1	1.5
Human serum 5	1791	30.0	1.7

a) PreciControl ClinChem Multi 1

b) PreciControl ClinChem Multi 2

*Urine (CREJ2U)*

Repeatability	Mean	SD	CV
	mmol/L	mmol/L	%
PN PUC <sup>c)</sup>	8.89	0.0922	1.0
PP PUC <sup>d)</sup>	4.56	0.0560	1.2
Human urine 1	1.19	0.0310	2.6
Human urine 2	2.41	0.0311	1.3
Human urine 3	22.5	0.210	0.9
Human urine 4	28.4	0.318	1.1
Human urine 5	49.7	0.480	1.0

Intermediate precision	Mean	SD	CV
	mmol/L	mmol/L	%
PN PUC <sup>c)</sup>	8.89	0.151	1.7
PP PUC <sup>d)</sup>	4.56	0.0824	1.8
Human urine 1	1.19	0.0341	2.9
Human urine 2	2.43	0.0398	1.6
Human urine 3	22.5	0.339	1.5
Human urine 4	28.4	0.417	1.5
Human urine 5	49.7	0.745	1.5

c) Precinorm PUC

d) Precipath PUC

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s) and **cobas c 703** analyzer(s).

**Method comparison**

Creatinine values for human serum, plasma and urine samples obtained on a **cobas c 503** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

*Serum/plasma (CREJ2)*

Sample size (n) = 71

Passing/Bablok <sup>30</sup>	Linear regression
$y = 1.012x - 3.68 \mu\text{mol/L}$	$y = 1.010x - 3.19 \mu\text{mol/L}$
$\tau = 0.980$	$r = 1.000$

The sample concentrations were between 23.2 and 2133  $\mu\text{mol/L}$ .

*Urine (CREJ2U)*

Sample size (n) = 72

Passing/Bablok <sup>30</sup>	Linear regression
$y = 1.065x - 0.0368 \text{ mmol/L}$	$y = 1.056x + 0.00514 \text{ mmol/L}$
$\tau = 0.984$	$r = 1.000$

The sample concentrations were between 0.388 and 50.8 mmol/L.

Creatinine values for human serum, plasma and urine samples obtained on a **cobas c 303** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

*Serum/plasma (CREJ2)*

Sample size (n) = 70

Passing/Bablok <sup>30</sup>	Linear regression
$y = 1.018x - 5.48 \mu\text{mol/L}$	$y = 1.015x - 4.42 \mu\text{mol/L}$
$\tau = 0.968$	$r = 1.000$

The sample concentrations were between 24.1 and 2114  $\mu\text{mol/L}$ .

*Urine (CREJ2U)*

Sample size (n) = 69

Passing/Bablok<sup>30</sup>

$$y = 1.088x - 0.0452 \text{ mmol/L}$$

$$\tau = 0.984$$

The sample concentrations were between 0.787 and 49.1 mmol/L.

Creatinine values for human serum, plasma and urine samples obtained on a **cobas c 703** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 503** analyzer (x).

*Serum/plasma (CREJ2)*

Sample size (n) = 86

Passing/Bablok<sup>30</sup>

$$y = 1.000x + 1.00 \mu\text{mol/L}$$

$$\tau = 0.983$$

## Linear regression

$$y = 1.093x - 0.0846 \text{ mmol/L}$$

$$r = 1.000$$

## Linear regression

$$y = 0.997x + 2.37 \mu\text{mol/L}$$

$$r = 1.000$$

The sample concentrations were between 25.5 and 2120  $\mu\text{mol/L}$ .

*Urine (CREJ2U)*

Sample size (n) = 89

Passing/Bablok<sup>30</sup>

$$y = 0.983x - 0.0174 \text{ mmol/L}$$

$$\tau = 0.993$$

## Linear regression

$$y = 0.984x - 0.0442 \text{ mmol/L}$$

$$r = 1.000$$

The sample concentrations were between 0.478 and 53.4 mmol/L.

**References**

- 1 Thomas C, Thomas L. Labordiagnostik von Erkrankungen der Nieren und ableitenden Harnwege. In: Thomas L, ed. Labor und Diagnose, 6th ed. Frankfurt/Main: TH-Books 2005;520-585.
- 2 Lamb E, Newman DJ, Price CP. Kidney function tests. In: Burtis CA, Ashwood ER, Bruns DE. Tietz textbook of clinical chemistry and molecular diagnostics. 4th ed. St.Louis, MO: Elsevier Saunders 2006;797-835.
- 3 KDIGO. KDIGO 2012 Clinical Practice Guideline for the Evaluation and Management of Chronic Kidney Disease <[https://kdigo.org/wp-content/uploads/2017/02/KDIGO\\_2012\\_CKD\\_GL.pdf](https://kdigo.org/wp-content/uploads/2017/02/KDIGO_2012_CKD_GL.pdf)> (2012).
- 4 Lamb EJ, Tomson CRV, Roderick PJ. Estimating kidney function in adults using formulae. Ann Clin Biochem 2005;42:321-345.
- 5 Miller WG. Editorial on Estimating glomerular filtration rate. Clin Chem Lab Med 2009;47(9):1017-1019.
- 6 Levey AS, Stevens LA, Schmid CH, et al. CKD-EPI (Chronic Kidney Disease Epidemiology Collaboration). A new equation to estimate glomerular filtration rate. Ann Intern Med 2009 May 5;150(9):604-12. doi: 10.7326/0003-4819-150-9-200905050-00006. Erratum in: Ann Intern Med 2011 Sep 20;155(6):408.
- 7 Debowska M, Wojcik-Zaluska A, Ksiazek A, et al. Phosphate, urea and creatinine clearances: haemodialysis adequacy assessed by weekly monitoring. Nephrol Dial Transplant 2015 Jan;30(1):129-36. doi: 10.1093/ndt/gfu266.
- 8 Tattersall J, Dekker F, Heimbürger O, et al. ERBP Advisory Board. When to start dialysis: updated guidance following publication of the Initiating Dialysis Early and Late (IDEAL) study. Nephrol Dial Transplant 2011 Jul;26(7):2082-6. doi: 10.1093/ndt/gfr168.
- 9 Schwartz GJ, Muñoz A, Schneider MF, et al. New Equations to Estimate GFR in Children with CKD. J Am Soc Nephrol 2009;20:629-637.
- 10 Schwartz GJ, Work DF. Measurement and Estimation of GFR in Children and Adolescents. Clin J Am Soc Nephrol 2009;4:1832-1843.
- 11 Staples A, LeBlond R, Watkins S, et al. Validation of the revised Schwartz estimating equation in a predominantly non-CKD population. Pediatr Nephrol 2010 Jul 22;25:2321-2326.
- 12 Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 1886;10:391-400.
- 13 Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the Centrifichem Clin Chem. 1971;17:696-700.

- 14 Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 1971;32:81-85.
- 15 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 16 Guder W, Fonseca-Wollheim W, Ehret W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
- 17 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 18 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 19 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 20 Ducharme MP, Smyth M, Strohs G. Drug-induced alterations in serum creatinine concentrations. Ann Pharmacotherapy 1993;27:622-633.
- 21 Kroll MH. Some observations on the reaction mechanism of Cefoxitin and Cephalothin with picrate. Microchem J 1990;42:241-249.
- 22 Mazzachi BC, Phillips JW, Peake MJ. Is the Jaffe creatinine assay suitable for neonates? Clin Biochem Revs 1998;19:82.
- 23 Filler G, Priem F, Lepage N, et al.  $\beta$ -Trace Protein, Cystatin C,  $\beta$ 2-Microglobulin, and Creatinine Compared for Detecting Impaired Glomerular Filtration Rates in Children. Clin Chem 2002;48:729-736.
- 24 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 25 Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinine Assays in Plasma and Serum and Early Morning Urine. Clin Lab 2000;53-55.
- 26 Schlebusch H, Liappis N, Kalina E, et al. High Sensitive CRP and Creatinine: Reference Intervals from Infancy to Childhood. J Lab Med 2002;26:341-346.
- 27 Junge W, Wilke B, Halabi A, et al. Determination of reference intervals for serum creatinine, creatinine excretion and creatinine clearance with an enzymatic and a modified Jaffé method. Clin Chim Acta 2004;344:137-148.
- 28 Zawta B, Delanghe J, Taes Y, et al. Arithmetic Compensation for Pseudo-Creatinine Interferences of the Creatinine Jaffé Method and its Effect on Creatinine Clearance Results. Clin Chem Part 2, Suppl S June 2001;46(6):487.
- 29 Wuyts B, Bernard D, van den Noortgate N, et al. Reevaluation of Formulas for Predicting Creatinine Clearance in Adults and Children Using Compensated Creatinine Methods. Clin Chem 2003;49:1011-1014.
- 30 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [navifyportal.roche.com](http://navifyportal.roche.com) for definition of symbols used):

CONTENT	Contents of kit
→	Volume for reconstitution
GTIN	Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

+800 5505 6606

