


Order information

REF		CONTENT		Analyzer(s) on which cobas c pack(s) can be used
08057605190	08057605500	Cardiac C-Reactive Protein (Latex) High Sensitive (600 tests)	System-ID 2048 001	cobas c 303 , cobas c 503 , cobas c 703

Materials required (but not provided):

11355279216	Calibrator f.a.s. Proteins (5 x 1 mL)	Code 20656	
20766321322	CRP T Control N (5 x 0.5 mL)	Code 20235	
10557897122	Precinorm Protein (3 x 1 mL)	Code 20302	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 20391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 20391	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

English

System information

CRPHS: ACN 20480

Intended use

In vitro test for the quantitative determination of C-reactive protein (CRP) in human serum and plasma on **cobas c** systems. Measurement of CRP is of use for the detection and evaluation of inflammatory disorders and associated diseases, infection and tissue injury. Highly sensitive measurement of CRP may also be used as an aid in the assessment of the risk of future coronary heart disease. When used as an adjunct to other laboratory evaluation methods of acute coronary syndromes, it may also be an additional independent indicator of recurrent event prognosis in patients with stable coronary disease or acute coronary syndrome.

Summary

C-reactive protein is the classic acute phase protein in inflammatory reactions.¹ It is synthesized by the liver and consists of five identical polypeptide chains that form a five-member ring having a molecular weight of 105000 daltons.^{1,2,3,4} CRP is the most sensitive of the acute phase reactants and its concentration increases rapidly during inflammatory processes.^{2,3} Complexed CRP activates the classical complement pathway. The CRP response frequently precedes clinical symptoms, including fever.^{1,3} After onset of an acute phase response the serum CRP concentration rises rapidly and extensively.^{2,3,4} The increase begins within 6 to 12 hours and the peak value is reached within 24 to 48 hours.^{1,3,5} CRP response may be less pronounced in patients suffering from liver disease.⁶

CRP assays are used to detect systemic inflammatory processes (apart from certain types of inflammation such as systemic lupus erythematosus (SLE) and Colitis ulcerosa);^{1,3,4,6} to assess treatment of bacterial infections with antibiotics;^{1,4,6,7} to detect intrauterine infections with concomitant premature amniorrhexis;^{4,6} to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa;^{3,4,6} to therapeutically monitor rheumatic disease and assess anti-inflammatory therapy;^{1,4,6} to determine the presence of post-operative complications at an early stage, such as infected wounds, thrombosis and pneumonia, and to distinguish between infection and bone marrow transplant rejection.^{1,4,6}

Sensitive CRP measurements have been used and discussed for early detection of infection in pediatrics and risk assessment of coronary heart disease.^{8,9,10,11} Several studies came to the conclusion that the highly sensitive measurement of CRP could be used as a marker to predict the risk of coronary heart disease in apparently healthy persons and as an indicator of recurrent event prognosis.^{10,12,13,14,15,16} Increases in CRP values are non-specific and should not be interpreted without a complete clinical history¹⁷. The American Heart Association and the Centers for Disease Control and Prevention have made several recommendations concerning the use of high sensitivity C-Reactive Protein (hsCRP) in cardiovascular risk assessment.^{17,18} Measurement of hsCRP may also be used as an aid in the assessment of the risk of future coronary heart disease and as a risk-enhancing factor in patients with borderline- or intermediate-risk for atherosclerotic cardiovascular disease.¹⁹ When used as an adjunct to other laboratory evaluation methods of acute coronary syndromes, it may also be an additional independent indicator of recurrent event prognosis in patients with stable coronary disease or acute coronary syndrome.^{17,20}

Testing for any risk assessment should not be performed while there is an indication of infection, systemic inflammation or trauma.^{11,17,21} Patients with

persistently unexplained hsCRP levels above 10 mg/L (95.2 nmol/L) should be evaluated for non-cardiovascular etiologies.^{13,17} When using hsCRP to assess the risk of coronary heart disease, measurements should be made on metabolically stable patients and compared to previous values.¹⁷ Optimally, the average of hsCRP results repeated two weeks apart should be used for risk assessment.¹⁷ Screening the entire adult population for hsCRP is not recommended, and hsCRP is not a substitute for traditional cardiovascular risk factors.¹⁷ Acute coronary syndrome management should not depend solely on hsCRP measurements.^{14,17} Similarly, application of secondary prevention measures should be based on global risk assessment and not solely on hsCRP measurements.¹⁷ Serial measurements of hsCRP should not be used to monitor treatment.¹⁷

Studies indicate an influence of gestational age on the kinetics of CRP in preterm infants, which may materialize as a blunted response to infection when comparing preterm and term newborns.^{22,23,24} This phenomenon, most likely due to immature liver function, may result in a lower sensitivity of CRP in the diagnosis of neonatal sepsis in preterm compared to term newborns.²⁵ In adult patients with advanced liver dysfunction, CRP levels are reduced in response to acute infection, however production is nevertheless maintained.²⁶ Although the liver is considered the main source of CRP, serum levels are not significantly lower in patients with cirrhosis than in other patients, and the predictive performance for infection is similar for patients with and without cirrhosis.²⁷

Various assay methods are available for CRP determination, such as nephelometry and turbidimetry.^{28,29} The Roche CRP assay is based on the principle of particle-enhanced immunological agglutination.

Test principle^{28,29}

Particle enhanced immunoturbidimetric assay.

Human CRP agglutinates with latex particles coated with monoclonal anti-CRP antibodies. The precipitate is determined turbidimetrically.

Reagents - working solutions

- R1** TRIS buffer with bovine serum albumin and immunoglobulins (mouse); preservative; stabilizers
- R3** Latex particles coated with anti-CRP (mouse) in glycine buffer; preservative; stabilizers

R1 is in position B and R3 is in position C.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 26 weeks

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.
Serum

Plasma: Li-heparin and K₂-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability:³⁰ 11 days at 15-25 °C
2 months at 2-8 °C
3 years at -20 °C (± 5 °C)

Freeze only once.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma**Test definition**

Reporting time	10 min		
Wavelength (sub/main)	- /546 nm		
Reagent pipetting	Diluent (H ₂ O)		
R1	62 µL	31 µL	
R3	21 µL	15 µL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	4.5 µL	-	-
Decreased	4.5 µL	6 µL	84 µL
Increased	4.5 µL	-	-

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

Calibration

Calibrators S1: H₂O
S2-S6: C.f.a.s. Proteins

Calibration mode Non-linear

Calibration frequency

Full calibration

- after reagent lot change
- as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against the reference preparation of the IRMM (Institute for Reference Materials and Measurements) BCR470/CRM470 (RPPHS - Reference Preparation for Proteins in Human Serum).³¹

Quality control

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 26 weeks. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

cobas c systems automatically calculate the analyte concentration of each sample in the unit mg/L (nmol/L, mg/dL).

Conversion factors: mg/L × 9.52 = nmol/L
mg/L × 0.1 = mg/dL

Limitations - interference

Criterion: Recovery within ± 10 % of initial value at CRP concentrations of approximately 1.0 mg/L.

Icterus:³² No significant interference up to an I index of 60 for conjugated bilirubin and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dL or 1026 µmol/L).

Hemolysis:³² No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 622 µmol/L or 1000 mg/dL).

Lipemia (Intralipid):³² No significant interference up to an L index of 600. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 200 IU/mL.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{33,34}

Therapeutic drugs: Significantly decreased CRP values may be obtained from samples taken from patients who have been treated with carboxypenicillins.

High dose hook-effect: No false result occurs up to a CRP concentration of 1000 mg/L.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.³⁵

Although measures were taken to minimize interference caused by human anti-mouse antibodies, erroneous findings may be obtained from samples taken from patients who have been treated with monoclonal mouse antibodies or have received them for diagnostic purposes.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

Limits and ranges**Measuring range**

0.15-20.0 mg/L (1.43-190 nmol/L)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:15 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 15.

Lower limits of measurement*Limit of Blank, Limit of Detection and Limit of Quantitation*

Limit of Blank	= 0.15 mg/L (1.43 nmol/L)
Limit of Detection	= 0.15 mg/L (1.43 nmol/L)
Limit of Quantitation	= 0.30 mg/L (2.86 nmol/L)

The Limit of Blank, the Limit of Detection and the Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration C-reactive protein samples.

Expected values**mg/L**

Consensus reference interval for adults:³⁶

IFCC/CRM 470 < 5.0 mg/L*

*calculated by unit conversion factor

The CDC/AHA recommended the following hsCRP cut-off points (tertiles) for CVD risk assessment:^{17,37}

hsCRP level (mg/L)	Relative risk
< 1.0	low
1.0-3.0	average
> 3.0	high

Patients with higher hsCRP concentrations are more likely to develop myocardial infarction and severe peripheral vascular disease.

5-95 % reference intervals of neonates and children:³⁸

Neonates (0-3 weeks):	0.1-4.1 mg/L
Children (2 months-15 years):	0.1-2.8 mg/L

nmol/L*

Consensus reference interval for adults:³⁶

IFCC/CRM 470 < 47.6 nmol/L

The CDC/AHA recommended the following hsCRP cut-off points (tertiles) for CVD risk assessment:^{17,37}

hsCRP level (nmol/L)	Relative risk
< 9.52	low
9.52-28.6	average
> 28.6	high

Patients with higher hsCRP concentrations are more likely to develop myocardial infarction and severe peripheral vascular disease.

5-95 % reference intervals of neonates and children:³⁸

Neonates (0-3 weeks):	0.95-39.0 nmol/L
Children (2 months-15 years):	0.95-26.7 nmol/L

*calculated by unit conversion factor

It is important to monitor the CRP concentration during the acute phase of the illness.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Increases in CRP values are non-specific and should not be interpreted without a complete clinical history.

When using hsCRP to assess the risk of coronary heart disease, measurements should be made on metabolically stable patients and compared to previous values. Optimally, the average of hsCRP results repeated two weeks apart should be used for risk assessment. Measurements should be compared to previous values. When the results are being used for risk assessment, patients with persistently unexplained hsCRP levels of above 10 mg/L (95.2 nmol/L) should be evaluated for non-cardiovascular origins. Testing for any risk assessment should not be performed while there is indication of infection, systemic inflammation or trauma.¹⁷

Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability ($n = 84$) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c 503** analyzer.

Repeatability	Mean mg/L	SD mg/L	CV %
PCCC1 ^{a)}	8.85	0.0329	0.4
Precinorm Protein	12.6	0.0715	0.6
Human serum 1	0.306	0.0237	7.7
Human serum 2	0.807	0.0261	3.2
Human serum 3	4.88	0.0304	0.6
Human serum 4	9.34	0.0332	0.4
Human serum 5	17.2	0.0919	0.5
Intermediate precision	Mean mg/L	SD mg/L	CV %
PCCC1 ^{a)}	9.03	0.0667	0.7
Precinorm Protein	12.6	0.0857	0.7
Human serum 1	0.306	0.0258	8.4
Human serum 2	0.807	0.0269	3.3
Human serum 3	4.88	0.0410	0.8
Human serum 4	9.73	0.0608	0.6
Human serum 5	17.2	0.140	0.8

a) PreciControl ClinChem Multi 1

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s) and **cobas c 703** analyzer(s).

Method comparison

CRP values for human serum and plasma samples obtained on a **cobas c 503** analyzer (y) were compared to those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 83

Passing/Bablok ³⁹	Linear regression
$y = 1.005x - 0.0232 \text{ mg/L}$	$y = 1.017x - 0.0836 \text{ mg/L}$

$\tau = 0.986$ $r = 0.999$

The sample concentrations were between 0.211 and 19.2 mg/L.

CRP values for human serum and plasma samples obtained on a **cobas c 303** analyzer (y) were compared to those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 67

Passing/Bablok³⁹

Linear regression

 $y = 0.980x + 0.0376 \text{ mg/L}$ $y = 0.977x + 0.0325 \text{ mg/L}$ $\tau = 0.992$ $r = 1.000$

The sample concentrations were between 0.200 and 19.6 mg/L.

CRP values for human serum and plasma samples obtained on a **cobas c 703** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 503** analyzer (x).

Sample size (n) = 70

Passing/Bablok³⁹

Linear regression

 $y = 0.997x + 0.0790 \text{ mg/L}$ $y = 0.997x + 0.0719 \text{ mg/L}$ $\tau = 0.983$ $r = 1.000$

The sample concentrations were between 0.240 and 18.6 mg/L.

References

- Thomas L. Labor und Diagnose, 7. Auflage, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main 2008;1010-1021.
- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer Verlag 1995;234-236.
- Burtis CA, Ashwood ER, eds. Tietz Fundamentals of Clinical Chemistry, 5th ed. Pa: WB Saunders Co 2001;332-333.
- Young B, Gleeson M, Cripps AW. C-reactive protein: A critical review. Pathology 1991;23:118-124.
- Sproston NR, Ashworth JJ. Role of C-Reactive Protein at Sites of Inflammation and Infection. Front Immunol 2018 Apr 13;9:754.
- Thomas L, Messenger M. Pathobiochemie und Labordiagnostik der Entzündung. Lab med 1993;17:179-194.
- Wasunna A, Whitelaw A, Gallimore R, et al. C-reactive protein and bacterial infection in preterm infants. Eur J Pediatr 1990 Mar;149(6):424-427.
- Kuller LH, Tracy RP, Shaten J, et al. Relation of c-reactive protein and coronary heart disease in the MRFIT nested case control study. Am J Epidemiol 1996;144:537-547.
- Ridker PM, Glynn RJ, Hennekens CH, et al. C-Reactive Protein Adds to the Predictive Value of Total and HDL Cholesterol in Determining Risk of First Myocardial Infarction. Circulation 1998;97:2007-2011.
- Ridker PM, Cushman M, Stampfer MJ, et al. Plasma Concentration of C-Reactive Protein and Risk of Developing Peripheral Vascular Disease. Circulation 1998;97:425-428.
- Järvisalo MJ, Harmoinen A, Hakanen M, et al. Elevated Serum C-Reactive Protein Levels and Early Arterial Changes in Healthy Children. Arterioscler Thromb Vasc Biol, (August) 2002;1323-1328.
- Ridker PM, Cushman M, Stampfer MJ, et al. Inflammation, Aspirin, and the Risk of Cardiovascular Disease in Apparently Healthy Men. N Eng J Med 1997;336(14):973-979.
- Danesh J, Wheeler JG, Hirschfield GM, et al. C-Reactive Protein and Other Circulating Markers of Inflammation in the Prediction of Coronary Heart Disease. N Eng J Med 2004;350(14):1387-1397.
- Ridker PM, Hennekens CH, Buring JE, et al. C-Reactive Protein and Other Markers of Inflammation in the Prediction of Cardiovascular Disease in Women. N Engl J Med 2000;342(12):836-843.
- Tracy RP, Lemaitre RN, Psaty BM, et al. Relationship of C-Reactive Protein to Risk of Cardiovascular Disease in the Elderly. Arterioscler Thromb Vasc Biol 1997;17:1121-1127.
- Almagor M, Keren A, Banai S. Increased C-Reactive Protein Level after Coronary Stent Implantation in Patients with Stable Coronary Artery Disease. American Heart Journal 2003;145 (2):248-253.
- Pearson TA, Mensah GA, Alexander RW, et al. Markers of inflammation and cardiovascular disease: application to clinical and public health practice: A statement for healthcare professionals from the Centers for Disease Control and Prevention and the American Heart Association. Circulation 2003 Jan 28;107(3):499-511.
- Arnett DK, Blumenthal RS, Albert MA, et al. ACC/AHA Guideline on the Primary Prevention of Cardiovascular Disease: A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines. J Am Coll Cardiol 2019 Sep 10;74(10):e177-e232.
- Grundy SM, Stone NJ, Bailey AL, et al. 2018 AHA/ACC/AACV-PR/AAPA/ABC/ACPM/ADA/AGS/APhA/ASPC/NLA/PCNA Guideline on the Management of Blood Cholesterol. J Am Coll Cardiol 2019;73:3168-3209.
- Arroyo-Espiguero R, Viana-Llamas MC, Silva-Obregón A, et al. The Role of C-reactive Protein in Patient Risk Stratification and Treatment. Eur Cardiol 2021;16:e28.
- Katritsis D, Korovesis S, Gazitzoglou E, et al. C-Reactive Protein Concentrations and Angiographic Characteristics of Coronary Lesions. Clin Chem 2001;47(5):882-886.
- Hofer N, Müller W, Resch B. Non-infectious conditions and gestational age influence C-reactive protein values in newborns during the first 3 days of life. Clin Chem Lab Med 2011 Feb;49(2):297-302.
- Turner MA, Power S, Emmerson AJB. Gestational age and the C reactive protein response. Arch Dis Child Fetal Neonatal Ed. 2004 May;89(3):F272-3.
- Doellner H, Arntzen KJ, Haereid PE, et al. Interleukin-6 concentrations in neonates evaluated for sepsis. J Pediatr 1998 Feb;132(2):295-9.
- Kawamura M, Nishida H. The usefulness of serial C-reactive protein measurement in managing neonatal infection. Acta Paediatr 1995 Jan;84(1):10-3.
- Park WB, Lee KD, Lee CS, et al. Production of C-reactive protein in Escherichia coli-infected patients with liver dysfunction due to liver cirrhosis. Diagn Microbiol Infect Dis 2005 Apr;51(4):227-30.
- Bota DP, Nuffelen MV, Zakariah AN, et al. Serum levels of C-reactive protein and procalcitonin in critically ill patients with cirrhosis of the liver. J Lab Clin Med 2005 Dec;146(6):347-51.
- Price CP, Trull AK, Berry D, et al. Development and validation of a particle-enhanced turbidimetric immunoassay for C-reactive protein. J Immunol Methods 1987;99:205-211.
- Eda S, Kaufmann J, Roos W, et al. Development of a New Microparticle-Enhanced Turbidimetric Assay for C-reactive Protein with Superior Features in Analytical Sensitivity and Dynamic Range. J Clin Lab Anal 1998;12:137-144.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- Baudner S, Bienvenu J, Blirup-Jensen S, et al. The certification of a matrix reference material for immunochemical measurement of 14 human serum proteins CRM470. Report EUR 15243 EN 1993;1-186.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.

CRPHS

Cardiac C-Reactive Protein (Latex) High Sensitive

- 36 Dati F, Schumann G, Thomas L, et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-520.
- 37 Ridker PM. Clinical Application of C-Reactive Protein for Cardiovascular Disease Detection and Prevention. Circulation 2003;107:363-369.
- 38 Schlebusch H, Liappis N, Kalina E, et al. High Sensitive CRP and Creatinine: Reference Intervals from Infancy to Childhood. J Lab Med 2002;26:341-346.
- 39 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety & Performance Report can be found here:
<https://ec.europa.eu/tools/eudamed>

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see navifyportal.roche.com for definition of symbols used):

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number
Rx only	For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, PRECICONTROL and PRECINORM are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68905 Mannheim
www.roche.com
 +800 5505 6606

