

Elecsys IGF-1

cobas®

| REF | | SYSTEM |
|--------------|-----|---------------------------------------------------------------------|
| 07475896 190 | 100 | MODULAR ANALYTICS E170 cobas e 411 cobas e 601 cobas e 602 |

English

System information

For **cobas e 411** analyzer: test number 1060
For MODULAR ANALYTICS E170, **cobas e 601** and **cobas e 602** analyzers: Application Code Number 555

Intended use

Immunoassay for the in vitro quantitative determination of insulin-like growth factor-1 (IGF-1) in human serum and plasma. The IGF-1 determination is intended to be used as an aid in the assessment of growth disorders in conjunction with other clinical and laboratory findings.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

Summary

IGF-1, a 70 amino-acid polypeptide with a molecular mass of 7.5 kDa,¹ is ubiquitously expressed in every tissue but it is mainly synthesized and secreted by the liver (~75 % of circulating IGF-1) and regulated by growth hormone (GH).² Around 80 % of IGF-1 in the circulation is bound in a ternary complex with IGFBP-3 (Insulin-like growth factor binding-protein 3) and ALS (Acid-labile subunit). The half-life of IGF-1 in this complex is around one hour. Another 20 % of IGF-1 is bound to IGFBP-3 without ALS. Only 1 % of IGF-1 is not bound at all with a half-life of only a few minutes.³

IGF-1 (also known as somatomedin)⁴ was the first established marker to screen for growth hormone deficiency (GHD).⁵ GH is secreted in pulses peaking every 60-90 minutes and has a short half-life. Additionally, GH levels are affected by external factors (e.g. exercise, fasting). In contrast, IGF-1 levels are more robust and as a consequence, the determination of IGF-1 is widely used as a first step in diagnosis of both GH deficiency and excess.⁶

Short stature in children is mainly caused by conditions that affect the growth plates. In case no reason is found, the diagnosis is idiopathic short stature (ISS). GHD is one such condition that affects the growth plates. In this context, IGF-1 is one of several laboratory parameters recommended in guidelines to identify the cause of short stature in children.⁷ In combination with other assessments an IGF-1 value around the mean-value of age or upper half of normal range of IGF-1 makes a GHD unlikely and no further testing would be required. Low IGF-1 concentrations (< 2 SD) would indicate a GHD with a high likelihood and should be confirmed with a GH-stimulation test. A GH-stimulation test would also be indicated with IGF-1 serum levels in the lower half of the normal range combined with clinical manifestations of GHD.⁸

GHD is also observed in adults. Interpretation of IGF-1 levels in the context of adults with GHD is different from short stature children. In adults a normal IGF-1 level does not exclude GHD. A very low IGF-1 level (< 2 SD) in patients with highly suspected GHD, or with long-lasting adult-onset, or multiple or total hypopituitarism may be considered as GHD without a GH-stimulation test.^{9,10}

The determination of IGF-1 is also recommended for screening of growth disorders by GH-excess like acromegaly.¹¹

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: Complexed antigen in the sample (10 µL) and diluted HCl react to cleave IGF-1 from IGFBP-3 and ALS.
- 2nd incubation: A biotinylated monoclonal IGF-1-specific antibody and a monoclonal IGF-1-specific antibody labeled with a ruthenium complex^{a)} react to form a sandwich complex. After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.

- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode or e-barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The reagent rackpack is labeled as IGF-1.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Diluted HCl (gray cap), 1 bottle, 10 mL: pH 1.5.
- R2 Anti-IGF-1-Ab-biotin, anti-IGF-1-Ab-Ru(bpy)₃²⁺ (black cap), 1 bottle, 10 mL: Biotinylated monoclonal anti-IGF-1 antibody (mouse) 0.6 µg/mL; monoclonal anti-IGF-1 antibody (mouse) labeled with ruthenium complex 1.5 µg/mL; phosphate buffer 100 mmol/L, pH 7.8; preservative.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

- H317 May cause an allergic skin reaction.

Prevention:

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves.

Response:

- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
- P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

Elecsys IGF-1



P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

| Stability: | |
|-------------------------|----------------------------------|
| unopened at 2-8 °C | up to the stated expiration date |
| after opening at 2-8 °C | 28 days |
| on the analyzers | 28 days |

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K₂-EDTA and K₃-EDTA plasma.

Plasma tubes containing separating gel can be used.

Criterion: Slope 0.9-1.1 + intercept within ± 7 ng/mL + coefficient of correlation ≥ 0.95 .

Stable for 24 hours at 15-25 °C, 48 hours at 2-8 °C, 28 days at -20 °C (± 5 °C). Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- REF 07475969190, CalSet IGF-1, for 4 x 1.0 mL
- REF 07476108190, PreciControl Growth, for 4 x 3.0 mL
- REF 05192943190, Diluent Universal 2, 2 x 36 mL sample diluent
- General laboratory equipment
- MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for **cobas e** 411 analyzer:

- REF 11662988122, ProCell, 6 x 380 mL system buffer
- REF 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- REF 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- REF 11933159001, Adapter for SysClean
- REF 11706802001, AssayCup, 60 x 60 reaction cups
- REF 11706799001, AssayTip, 30 x 120 pipette tips

- REF 11800507001, Clean-Liner
- Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- REF 04880340190, ProCell M, 2 x 2 L system buffer
- REF 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- REF 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- REF 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- REF 03004899190, PreClean M, 5 x 600 mL detection cleaning solution
- REF 12102137001, AssayTip/AssayCup, 48 magazines x 84 reaction cups or pipette tips, waste bags
- REF 03023150001, WasteLiner, waste bags
- REF 03027651001, SysClean Adapter M

Accessories for all analyzers:

- REF 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution
- REF 11298500160, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution (for USA)

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers: PreClean M solution is necessary.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized against the WHO International Standard 02/254.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 8 weeks when using the same reagent lot
- after 7 days when using the same reagent kit on the analyzer
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Growth.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Elecsys IGF-1



Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in ng/mL, µg/L or nmol/L).

Conversion factors:

$$\text{ng/mL} \times 1 = \mu\text{g/L}$$

$$\text{ng/mL} \times 0.131 = \text{nmol/L}$$

Limitations - interference

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations and no impact on results was observed.

Endogenous substances

| Compound | Concentration tested |
|--------------------|-------------------------------|
| Bilirubin | ≤ 1129 µmol/L or ≤ 66 mg/dL |
| Hemoglobin | ≤ 0.311 mmol/L or ≤ 500 mg/dL |
| Intralipid | ≤ 2000 mg/dL |
| Biotin | ≤ 205 nmol/L or ≤ 50 ng/mL |
| Rheumatoid factors | ≤ 1200 IU/mL |
| IgG | ≤ 3.3 g/dL |
| IgA | ≤ 0.5 g/dL |
| IgM | ≤ 1.0 g/dL |
| Albumin | ≤ 7.0 g/dL |

Criterion: Recovery within ± 4 ng/mL for IGF-1 concentrations ≤ 40 ng/mL or ± 10 % for concentrations > 40 ng/mL of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

There is no high-dose hook effect at IGF-1 concentrations up to 20000 ng/mL.

Pharmaceutical substances

In vitro tests were performed on 16 commonly used pharmaceuticals. No interference with the assay was found.

In addition, the following special growth disorder drugs were tested. No interference with the assay was found.

Special growth disorder drugs

| Drug | Concentration tested mg/L |
|--------------|---------------------------|
| Somatotropin | 3.0 |
| Octreotide | 1.5 |
| Pegvisomant | 80 |

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

7-1600 ng/mL (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 7 ng/mL. Values above the measuring range are reported as > 1600 ng/mL (or up to 16000 ng/mL for 10-fold diluted samples).

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank = 3.5 ng/mL

Limit of Detection = 7 ng/mL

Limit of Quantitation = 15 ng/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from n ≥ 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of ≤ 20 %.

Dilution

Samples with IGF-1 concentrations above the measuring range can be diluted with Diluent Universal 2. The recommended dilution is 1:10 (either automatically by the analyzers or manually). The concentration of the diluted sample must be > 140 ng/mL.

After manual dilution, multiply the result by the dilution factor.

After dilution by the analyzers, the software automatically takes the dilution into account when calculating the sample concentration.

Expected values

Expected values were obtained in a clinical study (CIM RD002173) that enrolled over 3000 female and over 3500 male subjects, including over 1400 subjects ≤ 17 years old.

See "Table for expected values" section for details.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, samples and controls in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

| cobas e 411 analyzer | | | | | |
|---------------------------|------------|---------------|------|------------------------|------|
| | | Repeatability | | Intermediate precision | |
| Sample | Mean ng/mL | SD ng/mL | CV % | SD ng/mL | CV % |
| Human serum 1 | 12.7 | 1.24 | 9.8 | 1.56 | 12.2 |
| Human serum 2 | 173 | 3.86 | 2.2 | 6.67 | 3.8 |
| Human serum 3 | 705 | 13.4 | 1.9 | 27.9 | 4.0 |
| Human serum 4 | 1594 | 49.6 | 3.1 | 98.4 | 6.2 |
| PC ^{b)} Growth 1 | 50.2 | 1.27 | 2.5 | 2.01 | 4.0 |
| PC Growth 2 | 358 | 8.26 | 2.3 | 14.0 | 3.9 |

b) PreciControl

| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---------------------------------------------------------------|------------|---------------|------|------------------------|------|
| | | Repeatability | | Intermediate precision | |
| Sample | Mean ng/mL | SD ng/mL | CV % | SD ng/mL | CV % |
| Human serum 1 | 14.5 | 0.298 | 2.0 | 0.518 | 3.6 |
| Human serum 2 | 171 | 2.35 | 1.4 | 4.65 | 2.7 |
| Human serum 3 | 693 | 11.4 | 1.6 | 17.5 | 2.5 |
| Human serum 4 | 1519 | 52.6 | 3.5 | 70.2 | 4.6 |
| PC Growth 1 | 50.0 | 0.508 | 1.0 | 1.42 | 2.8 |

Elecsys IGF-1



| MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers | | | | | |
|---------------------------------------------------------------|---------------|---------------|---------|------------------------|---------|
| | | Repeatability | | Intermediate precision | |
| Sample | Mean ng/mL | SD ng/mL | CV % | SD ng/mL | CV % |
| PC Growth 2 | 350 | 5.45 | 1.6 | 10.3 | 2.9 |

Method comparison

A comparison of the Elecsys IGF-1 assay (y) with a commercially available IGF-1 assay (x) using clinical samples gave the following correlations (ng/mL):

Number of samples measured: 146

Passing/Bablok¹²

Linear regression

$y = 1.13x - 14.0$

$y = 1.11x - 12.0$

$\tau = 0.969$

$r = 0.995$

The sample concentrations were between 11.2 ng/mL and 1200 ng/mL.

Analytical specificity

No significant cross-reactivity was found for the following substances:

| Substances | Concentration tested |
|------------|----------------------|
| IGF-2 | 4000 ng/mL |
| IGFBP-3 | 20000 ng/mL |
| Insulin | 1000 mIU/mL |
| Proinsulin | 1000 nmol/L |

Table for expected values

Table of age-dependent expected values; the values represent the indicated quantiles (2.5 %, 50 % and 97.5 %) for each age.

| Male subjects | | | | |
|----------------|----|------------------|-----------------|-------------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 0.25 | 41 | 12.0 | 39.4 | 94.1 |
| 0.5 | 44 | 11.8 | 40.9 | 94.6 |
| 1 | 59 | 11.8 | 44.2 | 96.4 |
| 2 | 38 | 13.9 | 51.7 | 104 |
| 3 | 28 | 18.9 | 60.5 | 116 |
| 4 | 29 | 26.8 | 70.6 | 134 |
| 5 | 34 | 36.6 | 81.9 | 156 |
| 6 | 51 | 47.1 | 94.5 | 184 |
| 7 | 34 | 57.5 | 108 | 216 |
| 8 | 58 | 67.5 | 123 | 254 |
| 9 | 61 | 76.9 | 141 | 296 |
| 10 | 51 | 85.7 | 164 | 343 |
| 11 | 49 | 93.9 | 194 | 392 |
| 12 | 47 | 101 | 231 | 434 |
| 13 | 42 | 108 | 270 | 467 |
| 14 | 35 | 115 | 304 | 489 |
| 15 | 15 | 120 | 327 | 501 |
| 16 | 13 | 125 | 339 | 503 |
| 17 | 4 | 129 | 340 | 495 |
| 18 | 1 | 132 | 331 | 476 |
| 19 | 2 | 134 | 312 | 450 |
| 20 | 4 | 136 | 291 | 421 |

| Male subjects | | | | |
|----------------|-----|------------------|-----------------|-------------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 21 | 10 | 137 | 272 | 394 |
| 22 | 10 | 137 | 254 | 370 |
| 23 | 16 | 136 | 238 | 348 |
| 24 | 19 | 135 | 225 | 328 |
| 25 | 25 | 132 | 213 | 310 |
| 26 | 15 | 130 | 203 | 295 |
| 27 | 19 | 128 | 194 | 282 |
| 28 | 16 | 125 | 188 | 271 |
| 29 | 18 | 123 | 183 | 263 |
| 30 | 18 | 120 | 180 | 257 |
| 31 | 17 | 118 | 176 | 253 |
| 32 | 16 | 116 | 173 | 250 |
| 33 | 15 | 114 | 170 | 247 |
| 34 | 21 | 111 | 166 | 244 |
| 35 | 14 | 109 | 163 | 242 |
| 36 | 16 | 107 | 160 | 239 |
| 37 | 16 | 105 | 158 | 236 |
| 38 | 19 | 103 | 155 | 234 |
| 39 | 18 | 101 | 152 | 231 |
| 40 | 39 | 98.5 | 150 | 229 |
| 41 | 92 | 96.4 | 148 | 226 |
| 42 | 93 | 94.4 | 146 | 223 |
| 43 | 101 | 92.4 | 144 | 221 |
| 44 | 99 | 90.5 | 142 | 218 |
| 45 | 75 | 88.5 | 140 | 216 |
| 46 | 100 | 86.5 | 139 | 214 |
| 47 | 98 | 84.6 | 137 | 211 |
| 48 | 79 | 82.6 | 136 | 209 |
| 49 | 88 | 80.6 | 135 | 207 |
| 50 | 97 | 78.7 | 133 | 205 |
| 51 | 61 | 76.7 | 132 | 203 |
| 52 | 78 | 74.8 | 130 | 201 |
| 53 | 76 | 72.8 | 129 | 200 |
| 54 | 54 | 70.9 | 127 | 198 |
| 55 | 62 | 68.9 | 126 | 196 |
| 56 | 44 | 67.0 | 124 | 195 |
| 57 | 63 | 65.3 | 122 | 194 |
| 58 | 70 | 63.7 | 121 | 193 |
| 59 | 70 | 62.3 | 119 | 192 |
| 60 | 61 | 61.1 | 118 | 191 |
| 61 | 58 | 60.0 | 117 | 190 |
| 62 | 85 | 59.2 | 116 | 189 |
| 63 | 62 | 58.5 | 116 | 188 |
| 64 | 64 | 57.9 | 115 | 188 |
| 65 | 46 | 57.4 | 115 | 187 |
| 66 | 57 | 56.8 | 115 | 186 |

Elecsys IGF-1



| Male subjects | | | | |
|---------------|----|---------------|--------------|----------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 67 | 53 | 56.3 | 115 | 186 |
| 68 | 58 | 55.8 | 115 | 185 |
| 69 | 68 | 55.2 | 114 | 185 |
| 70 | 68 | 54.7 | 114 | 185 |
| 71 | 68 | 54.1 | 113 | 184 |
| 72 | 64 | 53.6 | 111 | 184 |
| 73 | 72 | 53.0 | 110 | 184 |
| 74 | 40 | 52.4 | 108 | 184 |
| 75 | 39 | 51.9 | 106 | 184 |
| 76 | 32 | 51.3 | 104 | 184 |
| 77 | 27 | 50.7 | 102 | 184 |
| 78 | 19 | 50.2 | 99.0 | 184 |
| 79 | 14 | 49.6 | 96.1 | 184 |
| 80 | 0 | - | - | - |

| Female subjects | | | | |
|-----------------|----|---------------|--------------|----------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 0.25 | 28 | 13.8 | 48.8 | 86.4 |
| 0.5 | 35 | 15.4 | 50.9 | 92.0 |
| 1 | 37 | 18.7 | 55.3 | 104 |
| 2 | 34 | 26.1 | 65.0 | 128 |
| 3 | 48 | 34.2 | 76.0 | 155 |
| 4 | 42 | 43.2 | 88.2 | 185 |
| 5 | 50 | 53.0 | 102 | 216 |
| 6 | 49 | 63.6 | 116 | 250 |
| 7 | 37 | 75.0 | 133 | 286 |
| 8 | 47 | 87.3 | 154 | 324 |
| 9 | 39 | 99.9 | 180 | 363 |
| 10 | 42 | 112 | 210 | 398 |
| 11 | 50 | 123 | 244 | 427 |
| 12 | 54 | 132 | 278 | 451 |
| 13 | 38 | 140 | 306 | 468 |
| 14 | 38 | 146 | 325 | 480 |
| 15 | 21 | 151 | 331 | 485 |
| 16 | 11 | 154 | 324 | 485 |
| 17 | 14 | 156 | 305 | 479 |
| 18 | 5 | 156 | 283 | 466 |
| 19 | 3 | 155 | 261 | 449 |
| 20 | 13 | 152 | 243 | 429 |
| 21 | 7 | 148 | 227 | 410 |
| 22 | 7 | 143 | 214 | 392 |
| 23 | 15 | 138 | 203 | 375 |
| 24 | 16 | 134 | 195 | 359 |
| 25 | 15 | 130 | 189 | 343 |
| 26 | 18 | 126 | 185 | 329 |
| 27 | 13 | 122 | 182 | 315 |

| Female subjects | | | | |
|-----------------|----|---------------|--------------|----------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 28 | 13 | 118 | 179 | 303 |
| 29 | 14 | 115 | 176 | 292 |
| 30 | 10 | 112 | 173 | 281 |
| 31 | 12 | 109 | 171 | 271 |
| 32 | 10 | 107 | 169 | 263 |
| 33 | 7 | 104 | 167 | 255 |
| 34 | 10 | 102 | 165 | 248 |
| 35 | 11 | 100 | 163 | 242 |
| 36 | 9 | 98.3 | 160 | 238 |
| 37 | 14 | 96.5 | 158 | 234 |
| 38 | 15 | 94.8 | 155 | 231 |
| 39 | 6 | 93.1 | 153 | 228 |
| 40 | 51 | 91.4 | 150 | 227 |
| 41 | 74 | 89.8 | 147 | 225 |
| 42 | 88 | 88.1 | 145 | 224 |
| 43 | 79 | 86.5 | 142 | 222 |
| 44 | 71 | 84.9 | 139 | 221 |
| 45 | 72 | 83.3 | 136 | 220 |
| 46 | 53 | 81.8 | 132 | 219 |
| 47 | 70 | 80.2 | 130 | 218 |
| 48 | 69 | 78.7 | 127 | 218 |
| 49 | 94 | 77.2 | 125 | 217 |
| 50 | 59 | 75.7 | 123 | 215 |
| 51 | 47 | 74.3 | 121 | 214 |
| 52 | 52 | 72.8 | 120 | 212 |
| 53 | 48 | 71.4 | 119 | 210 |
| 54 | 44 | 70.0 | 118 | 207 |
| 55 | 68 | 68.6 | 117 | 204 |
| 56 | 46 | 67.3 | 117 | 201 |
| 57 | 55 | 65.9 | 116 | 198 |
| 58 | 51 | 64.6 | 115 | 194 |
| 59 | 36 | 63.3 | 114 | 190 |
| 60 | 59 | 62.0 | 113 | 186 |
| 61 | 60 | 60.7 | 112 | 182 |
| 62 | 55 | 59.5 | 111 | 179 |
| 63 | 57 | 58.3 | 110 | 176 |
| 64 | 47 | 57.3 | 109 | 173 |
| 65 | 40 | 56.3 | 108 | 170 |
| 66 | 50 | 55.5 | 106 | 168 |
| 67 | 41 | 54.8 | 105 | 166 |
| 68 | 71 | 54.2 | 104 | 164 |
| 69 | 45 | 53.8 | 102 | 163 |
| 70 | 48 | 53.5 | 101 | 162 |
| 71 | 59 | 53.3 | 99.8 | 161 |
| 72 | 47 | 53.2 | 98.7 | 160 |
| 73 | 44 | 53.2 | 97.6 | 160 |

Elecsys IGF-1

cobas®

| Female subjects | | | | |
|-----------------|----|---------------|--------------|----------------|
| Age (years) | N | 2.5 % (ng/mL) | 50 % (ng/mL) | 97.5 % (ng/mL) |
| 74 | 33 | 53.3 | 96.7 | 160 |
| 75 | 24 | 53.5 | 95.8 | 160 |
| 76 | 24 | 53.7 | 95.1 | 161 |
| 77 | 20 | 54.0 | 94.4 | 162 |
| 78 | 25 | 54.3 | 93.9 | 163 |
| 79 | 10 | 54.7 | 93.4 | 164 |
| 80 | 3 | 55.1 | 93.0 | 166 |

References

- Nieto-Estevez V, Defferali C, Vicario-Abejon C. IGF-I: A Key Factor that Regulates Neurogenesis and Synaptogenesis from Embryonic to Adult Stage of the Brain. *Front Neurosci* 2016;10:1-9.
- Aguirre GA, Rodriguez De Ita JR, de la Garzia RG, et al. Insulin-like growth factor-1 deficiency. *J Transl Med* 2016;14:3-26.
- Ranke MB. Insulin-like growth factor binding-protein-3 (IGFBP-3). *Best Pract Res Clin Endocrinol Metab* 2015;29:701-711.
- Binoux M, Hossenlopp P, Hardouin S, et al. Somatomedin (Insulin-Like Growth Factor)-Binding Proteins. *Horm Res* 1986;24:141-151.
- Moore DC, Rogelio HA, Smith EK, et al. Plasma Somatomedin-C as a Screening Test for Growth Hormone Deficiency in Children and Adolescents. *Horm Res* 1982;16:49-55.
- Ketha H, Singh RJ. Clinical assays for quantitation of insulin-like-growth-factor-1 (IGF1). *Methods* 2015;81:93-98.
- Wit JM, Clayton PE, Rogol AD, et al. Idiopathic short stature: Definition, epidemiology, and diagnostic evaluation. *Growth Horm IGF Res* 2008;18:89-110.
- Cohen P, Rogol AD, Deal CL, et al. Consensus Statement on the Diagnosis and Treatment of Children with Idiopathic Short Stature: A Summary of the Growth Hormone Research Society, the Lawson Wilkins Pediatric Endocrine Society, and the European Society for Paediatric Endocrinology Workshop. *J Clin Endocrinol Metab* 2008;93:4210-4217.
- Ho KKY. Consensus guidelines for the diagnosis and treatment of adults with GH deficiency II: a statement of the GH Research Society in association with the European Society for Pediatric Endocrinology, Lawson Wilkins Society, European Society of Endocrinology, Japan Endocrine Society, and Endocrine Society of Australia. *Eur J Endocrinol* 2007;157:695-700.
- Ghigo E, Aimaretti G, Corneli G. Diagnosis of adult GH deficiency. *Growth Hormone & IGF Res* 2008;18:1-16.
- Katznelson L, Laws ER, Melmed S, et al. Acromegaly: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab* 2014;99:3933-3951.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog. Roche.com for definition of symbols used):

| | |
|------------|-----------------------------------------------------|
| CONTENT | Contents of kit |
| SYSTEM | Analyzers/Instruments on which reagents can be used |
| REAGENT | Reagent |
| CALIBRATOR | Calibrator |
| → | Volume after reconstitution or mixing |
| GTIN | Global Trade Item Number |

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2020, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

+800 5505 6606



Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336