

REF			SYSTEM
07027559190*			cobas e 402
07027559214*	07027559500	100	cobas e 801

* Some kits shown may not be available in all countries.

English

System information

Short name	ACN (application code number)
INSULIN	10059

Intended use

Immunoassay for the in vitro quantitative determination of human insulin in human serum and plasma. The determination of insulin is utilized in the diagnosis and therapy of various disorders of carbohydrate metabolism, including diabetes mellitus and hypoglycemia.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on **cobas e** immunoassay analyzers.

Summary

Insulin is a 51-residue peptide hormone with a molecular weight of 5808 Da. It is secreted by the β -cells of the islets of Langerhans in the pancreas, and passes into circulation via the portal vein and the liver. Insulin is generally released in pulses.^{1,2}

The biologically active insulin molecule is monomeric and consists of two polypeptide chains, the 21 amino acid α -chain and the 30 amino acid β -chain joined by disulphide bridges. Insulin is the biosynthetic product of the single-chain precursor proinsulin, which is subsequently cleaved to give proinsulin.^{2,3,4,5} Specific proteases further cleave proinsulin to produce insulin and the connecting (C)-peptide which pass into the bloodstream simultaneously in equimolar concentrations. Circulating insulin has a half-life of 3-5 minutes and is preferentially retained and degraded in the liver. Therefore only about half of the insulin reaches the systemic circulation. Inactivation or excretion of proinsulin and C-peptide mainly takes place in the kidney and virtually none of the C-peptide is retained in the liver. As a result, C-peptide has a higher plasma concentration than insulin.⁶

The amino acid sequence of insulin is extremely well conserved, with the result that prior to the development of genetically engineered human insulin it was possible to successfully use porcine or bovine insulin in the therapy of diabetes mellitus.⁷

The action of insulin is mediated by specific receptors and primarily consists of facilitation of glucose uptake by the cells of the liver, fatty tissue and musculature; this is the basis of its hypoglycemic action.^{2,8}

Serum insulin determinations are mainly performed on patients with symptoms of hypoglycemia and may be useful in classifying the different types of diabetes.^{9,10} They are used to ascertain the glucose/insulin quotients and for clarification of questions concerning insulin secretion and β -cell function, e.g. in the evaluation of oral glucose tolerance tests or hunger provocation tests.¹¹

A disorder in insulin metabolism can have a significant impact on a number of metabolic processes. Low concentrations of free, biologically active insulin can lead to the development of diabetes mellitus. Possible causes of this include destruction of the β -cells (type I diabetes), reduced activity of insulin or reduced pancreatic synthesis (type II), circulating antibodies to insulin, delayed release of insulin or the absence (or inadequacy) of insulin receptors.¹²

Conversely, autonomous, non-regulated insulin secretion is generally the cause of hypoglycemia. This condition is brought about by inhibition of gluconeogenesis, e.g. as a result of severe hepatic or renal failure, islet cell adenoma, or carcinoma. Hypoglycemia can, however, also be facilitated intentionally or unintentionally (factitious hypoglycemia).^{10,13}

In certain individuals with reduced glucose tolerance, the metabolic state deteriorates towards diabetes mellitus over a period of time. Reduced glucose tolerance during pregnancy always requires treatment. The clearly elevated risk of mortality for the fetus necessitates intensive monitoring.¹²

The Elecsys Insulin assay employs two monoclonal antibodies which are specific for human insulin.

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: Insulin from 12 μ L sample, a biotinylated monoclonal insulin-specific antibody, and a monoclonal insulin-specific antibody labeled with a ruthenium complex^{a)} form a sandwich complex.
- 2nd incubation: After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell II M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the **cobas** link.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The **cobas e** pack is labeled as INSULIN.

- M Streptavidin-coated microparticles, 1 bottle, 5.8 mL:
Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-insulin-Ab~biotin, 1 bottle, 10.3 mL:
Biotinylated monoclonal anti-insulin antibody (mouse) 1 mg/L; MES^{b)} buffer 50 mmol/L, pH 6.0; preservative.
- R2 Anti-insulin-Ab~Ru(bpy)₃²⁺, 1 bottle, 9.5 mL:
Monoclonal anti-insulin antibody (mouse) labeled with ruthenium complex 1.75 mg/L; MES buffer 50 mmol/L, pH 6.0; preservative.

b) MES = 2-morpholino-ethane sulfonic acid

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Warning

H317 May cause an allergic skin reaction.

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

Response:

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is available via the **cobas** link.

Storage and stability

Store at 2-8 °C.

Do not freeze.

Store the **cobas e** pack **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
on the analyzers	16 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable.

Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K₂-EDTA and K₃-EDTA plasma.

Criterion: Slope 0.9-1.1 + intercept within $\pm 0.8 \mu\text{U/mL}$ + coefficient of correlation ≥ 0.95 .

Stable for 4 hours at 20-25 °C, 2 days at 2-8 °C, 6 months at -20 °C (± 5 °C). Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples and calibrators are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples and calibrators on the analyzers should be analyzed/measured within 2 hours.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 12017504122, Insulin CalSet, for 4 x 1.0 mL
- [REF] 05341787190, PreciControl Multimarker, for 6 x 2.0 mL or [REF] 11731416190, PreciControl Universal, for 4 x 3.0 mL

▪ General laboratory equipment

▪ **cobas e** analyzer

Additional materials for **cobas e** 402 and **cobas e** 801 analyzers:

- [REF] 06908799190, ProCell II M, 2 x 2 L system solution
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 07485409001, Reservoir Cup, 8 cups to supply ProCell II M and CleanCell M
- [REF] 06908853190, PreClean II M, 2 x 2 L wash solution
- [REF] 05694302001, Assay Tip/Assay Cup tray, 6 magazines x 6 magazine stacks x 105 assay tips and 105 assay cups, 3 wasteliners

- [REF] 07485425001, Liquid Flow Cleaning Cup, 2 adaptor cups to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning Detection Unit
- [REF] 07485433001, PreWash Liquid Flow Cleaning Cup, 1 adaptor cup to supply ISE Cleaning Solution/Elecsys SysClean for Liquid Flow Cleaning PreWash Unit
- [REF] 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use.

Place the cooled (stored at 2-8 °C) **cobas e** pack on the reagent manager. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the **cobas e** pack.

Calibration

Traceability: This method has been standardized using the 1st IRP WHO Reference Standard 66/304 (NIBSC).

The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration frequency: Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the **cobas e** pack was registered on the analyzer).

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 12 weeks when using the same reagent lot
- after 28 days when using the same **cobas e** pack on the analyzer
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Multimarker or PreciControl Universal. In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per **cobas e** pack, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Please note: Commercial controls may contain insulin of animal origin. When assessing results, the corresponding cross-reactivity of this test must be taken into account; see under "Analytical specificity".

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in $\mu\text{U/mL}$ or pmol/L).

$$\begin{aligned} \text{Conversion factors:} \quad & \mu\text{U/mL} \times 6.945 = \text{pmol/L} \\ & \text{pmol/L} \times 0.144 = \mu\text{U/mL} \end{aligned}$$

Limitations - interference

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations and no impact on results was observed.

Endogenous substances

Compound	Concentration tested
Bilirubin	$\leq 1539 \mu\text{mol/L}$ or $\leq 90 \text{ mg/dL}$
Intralipid	$\leq 1800 \text{ mg/dL}$
Biotin	$\leq 246 \text{ nmol/L}$ or $\leq 60 \text{ ng/mL}$

Compound	Concentration tested
Rheumatoid factors	≤ 1200 IU/mL

Criterion: For concentrations of 0.4-2 µU/mL the deviation is ≤ 0.5 µU/mL. For concentrations > 2 µU/mL the deviation is ≤ 10 %.

Hemolysis interferes, as insulin-degrading peptidases are released from erythrocytes.¹⁴

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

There is no high-dose hook effect at insulin concentrations up to 20000 µU/mL or 138900 pmol/L.

Pharmaceutical substances

In vitro tests were performed on 16 commonly used pharmaceuticals. Of these, only acetylcysteine at therapeutic dosage levels showed interference with the assay (insulin values depressed).

In addition, the following special drugs were tested. No interference with the assay was found.

Special drugs

Drug	Concentration tested mg/L
Euglucon	10.5
Tolbutamide	3

Samples from patients treated with bovine, porcine or human insulin sometimes contain anti-insulin antibodies.^{15,9} Insulin bound to these antibodies is at least partially recognized by the antibodies used in the Elecsys Insulin assay.¹⁶

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

0.4-1000 µU/mL or 2.78-6945 pmol/L (defined by the Limit of Detection and the maximum of the master curve). Values below the Limit of Detection are reported as < 0.4 µU/mL (< 2.78 pmol/L). Values above the measuring range are reported as > 1000 µU/mL (> 6945 pmol/L).

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank = 0.2 µU/mL (1.39 pmol/L)

Limit of Detection = 0.4 µU/mL (2.78 pmol/L)

Limit of Quantitation = 1 µU/mL (6.95 pmol/L)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from n ≥ 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with an intermediate precision CV of ≤ 20 %.

Dilution

Not necessary due to the broad measuring range.

Expected values

Studies with the Elecsys Insulin assay conducted in a clinical center in Germany with samples from 57 healthy, fasting individuals gave the following results (5th-95th percentile range):

2.6-24.9 µU/mL (17.8-173 pmol/L)

Status: Elecsys Insulin MCE, study No.: B99P027 of 29 March 2001.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, pooled human sera and controls in a protocol (EP05-A3) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

cobas e 402 and cobas e 801 analyzers					
Sample	Mean µU/mL	Repeatability		Intermediate precision	
		SD µU/mL	CV %	SD µU/mL	CV %
Human serum 1	22.1	0.310	1.4	0.466	2.1
Human serum 2	3.25	0.141	4.3	0.172	5.3
Human serum 3	49.7	0.403	0.8	0.719	1.4
Human serum 4	505	5.49	1.1	8.19	1.6
Human serum 5	973	10.4	1.1	14.2	1.5
PC ^{c)} Multimarker 1	20.2	0.280	1.4	0.400	2.0
PC Multimarker 2	66.8	0.692	1.0	1.06	1.6

c) PC = PreciControl

cobas e 402 and cobas e 801 analyzers					
Sample	Mean pmol/L	Repeatability		Intermediate precision	
		SD pmol/L	CV %	SD pmol/L	CV %
Human serum 1	153	2.15	1.4	3.24	2.1
Human serum 2	22.6	0.979	4.3	1.19	5.3
Human serum 3	345	2.80	0.8	4.99	1.4
Human serum 4	3507	38.1	1.1	56.9	1.6
Human serum 5	6757	72.2	1.1	98.6	1.5
PC Multimarker 1	140	1.94	1.4	2.78	2.0
PC Multimarker 2	464	4.81	1.0	7.36	1.6

Method comparison

a) A comparison of the Elecsys Insulin assay, [REF] 07027559190 (cobas e 801 analyzer; y) with the Elecsys Insulin assay, [REF] 12017547122 (cobas e 601 analyzer; x) gave the following correlations (µU/mL):

Number of serum samples measured: 164

Passing/Bablok¹⁷ Linear regression
 $y = 0.988x - 0.048$ $y = 0.973x + 1.09$
 $r = 0.993$ $r = 1.00$

The sample concentrations were between 0.924 and 989 µU/mL.

b) A comparison of the Elecsys Insulin assay, [REF] 07027559190 (cobas e 402 analyzer; y) with the Elecsys Insulin assay, [REF] 07027559190 (cobas e 801 analyzer; x) gave the following correlations (µU/mL):

Number of serum samples measured: 198

Passing/Bablok¹⁷ Linear regression
 $y = 0.996x + 0.133$ $y = 0.993x + 0.173$

Elecsys Insulin

 $\tau = 0.984$
 $r = 1.00$

The sample concentrations were between 0.694 and 957 $\mu\text{U/mL}$.

Analytical specificity

For the monoclonal antibodies used, the following cross-reactivities were found:

Cross-reactant	Concentration tested	Cross-reactivity %
Bovine insulin	20000 pmol/L	9.2
Porcine insulin	10000 pmol/L	22.2
Human proinsulin	111083 pmol/L	0.36
C-peptide	33109 pmol/L	n. d. ^{d)}
Glucagon	288 pmol/L	n. d.
Somatostatin	60 pmol/L	n. d.
Insulin-like growth factor I	10000 pmol/L	n. d.

d) n. d. = not detectable

Results for cross-reactivity with recombinant insulin analogs in a number of insulin methods have been published for example by two groups in France and the USA.^{16,18,19} The following results were published by Owen et al.¹⁸ for the Elecsys Insulin assay:

Insulin lispro, insulin aspart, and insulin glargine were each tested in concentrations of 30, 100, 300, and 1000 mIU/L in the absence of insulin. The results obtained were below the detection limit of the Elecsys Insulin assay ($< 0.4 \mu\text{U/mL}$ or $< 2.78 \text{ pmol/L}$) at all the concentrations tested.

Moreover, these results also correlate with those published earlier by Sapin et al. for insulin lispro.¹⁶

References

- Lang DA, Matthews DR, Peto J, et al. Cyclic oscillations of basal plasma glucose and insulin concentrations in human beings. *N Engl J Med* 1979;301:1023-1027.
- Steiner DF. Adventures with insulin in the islets of Langerhans. *J Biol Chem* 2011;286(20):17399-17421.
- Weiss MA. Diabetes mellitus due to the toxic misfolding of proinsulin variants. *FEBS Lett* 2013;587(13):1942-1950.
- Menting JG, Whittaker J, Margetts MB, et al. How insulin engages its primary binding site on the insulin receptor. *Nature* 2013;493(7431):241-245.
- Yang Y, Hua QX, Liu J, et al. Solution structure of proinsulin connecting domain flexibility and prohormone processing. *J Biol Chem* 2010;285(11):7847-7851.
- Ghorbani A, Shafiee-Nick R. Pathological consequences of C-peptide deficiency in insulin-dependent diabetes mellitus. *World J Diabetes* 2015;6(1):145-150.
- Pickup J. Human insulin. *Br Med J* 1986;292(6514):155-157.
- Sonksen P, Sonksen J. Insulin: understanding its action in health and disease. *Br J Anaesth* 2000;85:69-79.
- Clark PM. Assays for insulin, proinsulin(s) and C-peptide. *Ann Clin Biochem* 1999;36(5):541-564.
- Nalbantoğlu Elmas Ö, Demir K, Soyulu N, et al. Importance of insulin immunoassays in the diagnosis of factitious hypoglycemia. *J Clin Res Pediatr Endocrinol* 2014;6(4):258-261.
- Cersosimo E, Solis-Herrera C, Trautmann ME, et al. Assessment of pancreatic β -cell function: review of methods and clinical applications. *Curr Diabetes Rev* 2014;10:2-42.
- Kharroubi AT, Darwish HM. Diabetes mellitus: The epidemic of the century. *World J Diabetes*. 2015;6(6):850-867.
- Iglesias P, Díez JJ. Management of endocrine disease: a clinical update on tumor-induced hypoglycaemia. *Eur J Endocrinol* 2014;170(4):R147-157.
- Chevenne D, Letailleur A, Trivin F, et al. Effect of Hemolysis on the Concentration of Insulin in Serum Determined by RIA and IRMA. *Clin Chem* 1998;44(2):354-356.
- Tietz NW. *Clinical Guide To Laboratory Tests*. 3rd ed. Philadelphia, Pa: WB Saunders Co, 1995:366-367.
- Sapin R, Le Galudec V, Gasser F, et al. Elecsys Insulin Assay: Free Insulin Determination and the Absence of Cross-Reactivity with Insulin Lispro. *Clin Chem* 2001;47:602-605.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
- Owen WE, Roberts WL. Letter to the Editor: Cross-Reactivity of Three Recombinant Insulin Analogs with Five Commercial Insulin Immunoassays. *Clin Chem* 2004;50(1):257-259.
- Sapin R. Review: Insulin Assays: Previously Known and New Analytical Features. *Clin Lab* 2003;49(3-4):113-121.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog.roche.com for definition of symbols used):

	Contents of kit
	Analyzers/Instruments on which reagents can be used
	Reagent
	Calibrator
	Volume for reconstitution
	Global Trade Item Number

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2022, Roche Diagnostics

0123

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

+800 5505 6606

