


IRON2

Iron Gen.2

Order information**cobas**[®]

REF		CONTENT		Analyzer(s) on which cobas c pack(s) can be used
08057931190*	08057931500	Iron Gen.2 (700 tests)	System-ID 2077 001	cobas c 303, cobas c 503, cobas c 703
08057931214*	08057931500	Iron Gen.2 (700 tests)	System-ID 2077 001	cobas c 303, cobas c 503, cobas c 703

Materials required (but not provided):

10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 20401	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 20391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 20391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 20392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 20392	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

* Some kits shown may not be available in all countries.

English**System information**

IRON2: ACN 20770

Intended useIn vitro test for the quantitative determination of iron in human serum and plasma on **cobas c** systems.**Summary**

Iron measurements performed with this assay in human serum and plasma are used as an aid in diagnosis and monitoring of iron deficiency and iron overload disorders.

Iron is essential for many metabolic and biochemical processes. Similar to other micronutrients in the human body, iron is supplied with food. Ingested iron is mainly absorbed in the form of Fe²⁺ in the duodenum and proximal jejunum. The trivalent form and the heme-bound Fe³⁺-component of iron in food has to be reduced by duodenal cytochrome B. About 1-2 mg of iron is absorbed and lost daily. Upon reaching the mucosal cells, Fe²⁺ ions become bound to transport proteins. In the cellular phase iron is either stored in cellular ferritin or transported to the circulation. Iron export into the circulation requires Fe²⁺ oxidation to Fe³⁺ by hephaestin (on cellular membrane) or ceruloplasmin (in the circulation), for loading onto transferrin. Circulating Fe ions are transported by transferrin-iron complexes. A maximum of 2 Fe³⁺ ions per protein molecule can be transported.¹

Serum iron fluctuates with dietary intake and normal diurnal variation. Clinically, dysregulation of serum/plasma iron levels can be divided into iron deficiency and iron overload.^{1,2} Iron deficiency disorders can be due to increased demands (e.g. growth, pregnancy), limited external supply (e.g. malnutrition, inappropriate diet, malabsorption), increased loss (e.g. hemorrhage, hemodialysis, blood donation), or other conditions, such as chronic kidney disease resulting in renal anemia, inflammatory bowel disease, heart failure, obesity, bone marrow disease.^{2,3} Iron deficiency occurs in several stages, defined by the extent of depletion, first of iron stores and then of iron available for hemoglobin synthesis. In the first stage, iron stores can be completely depleted without causing anemia. Further iron loss causes anemia (iron deficiency anemia, IDA), which is initially normocytic, with a normal absolute reticulocyte count. Deeper deficiency results in classic anemia findings with hypochromic (low mean corpuscular hemoglobin) and microcytic (low mean corpuscular volume) red blood cells.^{3,4,5} Another type of anemia is macrocytic anemia (elevated mean corpuscular volume), which is not directly due to iron deficiency but is rather related to other causes, such as vitamin B12 and folate deficiency, bone marrow disorders (myelodysplasia), use of certain medications, alcohol abuse, liver disease, marked reticulocytosis, and hypothyroidism. Iron measurements can help define different causes of anemia.^{6,7}

Iron overload disorders normally result in increased serum/plasma iron concentration, and can be due to a number of underlying conditions, most commonly hereditary haemochromatosis (excess iron derived from increased gastrointestinal absorption due to inactivating mutations in components of the hepcidin pathway) and thalassemia (increased concentrations of iron mainly caused from regular red blood cell transfusions and to a lesser extent by increased iron absorption).^{2,8}


Numerous photometric methods have been described for the determination of iron. All have the following in common:

- Liberation of Fe³⁺ ions from the transferrin complex using acids or detergents.
- Reduction of Fe³⁺ ions to Fe²⁺ ions.
- Reaction of the Fe²⁺ ions to give a colored complex.¹

The method described here is based on the FerroZine method without deproteinization.

Test principle

Colorimetric assay.



Under acidic conditions, iron is liberated from transferrin. Lipemic samples are clarified by the detergent. Ascorbate reduces the released Fe³⁺ ions to Fe²⁺ ions which then react with FerroZine to form a colored complex. The color intensity is directly proportional to the iron concentration and can be measured photometrically.

Reagents - working solutions**R1** Citric acid: 200 mmol/L; thiourea: 115 mmol/L; detergent**R3** Sodium ascorbate: 150 mmol/L; FerroZine: 6 mmol/L; preservative

R1 is in position B and R3 is in position C.

Precautions and warnings

For in vitro diagnostic use for laboratory professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

IRON2

Iron Gen.2

cobas®



Danger

H318 Causes serious eye damage.

Prevention:

P280 Wear eye protection/ face protection.

Response:

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
 + P338 Continue rinsing. Immediately call a POISON CENTER/ doctor.
 + P310

Hazardous components:

- Poly(oxy-1,2-ethanediyl), .alpha.-isotridecyl-.omega.-hydroxy-

EUH 208 Contains DIAZOLIDINYL UREA. May produce an allergic reaction.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

When removing the **cobas c** pack from the instrument during use, please immediately store at 2-8 °C.Do not shake the **cobas c** pack to avoid foaming.**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma. Do not use EDTA or oxalate plasma.

Separate serum or plasma from the clot or cells within 1 hour.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability:^{9,10} 7 days at 15-25 °C

3 weeks at 2-8 °C

several years at -20 °C (± 5 °C)

Freeze only once.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma**Test definition**

Reporting time	10 min		
Wavelength (sub/main)	700/570 nm		
Reagent pipetting		Diluent (H ₂ O)	
R1	75 µL	–	
R3	15 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	6.4 µL	–	–
Decreased	9.0 µL	25 µL	50 µL
Increased	6.4 µL	–	–

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

Calibration

Calibrators	S1: H ₂ O
	S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	Full calibration
	- after reagent lot change
	1-point recalibration using S1
	- after cobas c pack green change
	- every 7 days on-board

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference material (SRM 937).

Quality control

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 12 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation**cobas c** systems automatically calculate the analyte concentration of each sample in the unit µmol/L (µg/dL, mg/L).

Conversion factors: µmol/L x 5.59 = µg/dL

µmol/L x 0.0559 = mg/L

Limitations - interference

Criterion: Recovery within ± 2.7 µmol/L of initial values of samples ≤ 26.9 µmol/L and within ± 10 % for samples > 26.9 µmol/L.

Icterus:¹¹ No significant interference up to an I index of 60 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:¹¹ No significant interference up to an H index of 200 (approximate hemoglobin concentration: 125 µmol/L or 200 mg/dL). Higher hemoglobin concentrations lead to artificially increased values due to contamination of the sample with hemoglobin-bound iron.

Lipemia (Intralipid):¹¹ No significant interference up to an L index of 1500. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{12, 13}

In patients treated with iron supplements or metal-binding drugs, the drug-bound iron may not properly react in the test, resulting in artificially low values.

In the presence of high ferritin concentrations > 1200 µg/L the assumption that serum iron is almost completely bound to transferrin is not valid anymore. Therefore, such iron results should not be used to calculate Total Iron Binding Capacity (TIBC) or percent transferrin saturation (% SAT).¹⁴

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.¹⁵

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOH/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

Limits and ranges

Measuring range

0.90-179 µmol/L (5.00-1000 µg/dL, 0.05-10.0 mg/L)

Determine samples having higher concentrations via the rerun function. For samples with higher concentrations, the rerun function decreases the sample volume by a factor of 2.1. The results are automatically multiplied by this factor.

Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank = 0.9 µmol/L (5.03 µg/dL)

Limit of Detection = 0.9 µmol/L (5.03 µg/dL)

Limit of Quantitation = 0.9 µmol/L (5.03 µg/dL)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration iron samples.

Expected values¹⁶

µmol/L

Adults: 5.83-34.5 µmol/L

µg/dL

Adults: 33-193 µg/dL

The concentration of iron in serum/plasma is dependent on ingestion of iron and is subject to circadian variations.¹⁷

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability ($n = 84$) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c 503** analyzer.

Repeatability	Mean µmol/L	SD µmol/L	CV %
PCCC ^{1a)}	18.6	0.111	0.6
PCCC ^{2b)}	41.4	0.163	0.4
Human serum 1	2.37	0.0817	3.4
Human serum 2	6.01	0.0830	1.4
Human serum 3	35.1	0.135	0.4
Human serum 4	89.2	0.307	0.3
Human serum 5	158	0.655	0.4

Intermediate precision	Mean µmol/L	SD µmol/L	CV %
PCCC ^{1a)}	18.6	0.212	1.1
PCCC ^{2b)}	41.6	0.369	0.9
Human serum 1	2.32	0.120	5.2
Human serum 2	5.95	0.149	2.5
Human serum 3	35.1	0.187	0.5
Human serum 4	89.2	0.337	0.4
Human serum 5	158	0.673	0.4

a) PreciControl ClinChem Multi 1

b) PreciControl ClinChem Multi 2

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s) and **cobas c 703** analyzer(s).

Method comparison

Iron values for human serum and plasma samples obtained on a **cobas c 503** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 74

Passing/Bablok ¹⁸	Linear regression
$y = 1.004x + 0.0354 \text{ µmol/L}$	$y = 1.003x + 0.00110 \text{ µmol/L}$
$r = 0.985$	$r = 1.000$

The sample concentrations were between 1.20 and 169 µmol/L.

Iron values for human serum and plasma samples obtained on a **cobas c 303** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 98

Passing/Bablok ¹⁸	Linear regression
$y = 1.011x - 0.0750 \text{ µmol/L}$	$y = 1.011x - 0.0772 \text{ µmol/L}$
$r = 0.993$	$r = 1.000$

The sample concentrations were between 1.72 and 172 µmol/L.

IRON2

Iron Gen.2

Iron values for human serum and plasma samples obtained on a **cobas c** 703 analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c** 503 analyzer (x).

Sample size (n) = 75

Passing/Bablok ¹⁸	Linear regression
$y = 1.000x + 0.0000 \mu\text{mol/L}$	$y = 1.001x - 0.0494 \mu\text{mol/L}$
$\tau = 0.998$	$r = 1.000$

The sample concentrations were between 2.55 and 176 $\mu\text{mol/L}$.

References

- Swinkels DW. Iron Metabolism. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 40, p. 418-418.e40.
- Muñoz M, García-Erce JA, Remacha ÁF. Disorders of iron metabolism. Part II: iron deficiency and iron overload. J Clin Pathol 2011 Apr;64(4):287-296.
- Camaschella C. Iron deficiency: new insights into diagnosis and treatment. Hematology Am Soc Hematol Educ Program. 2015;2015:8-13. doi: 10.1182/asheducation-2015.1.8.
- Cook JD, Skikne BS. Iron deficiency: definition and diagnosis. J Intern Med 1989 Nov;226(5):349-355.
- Camaschella C. Iron-deficiency anemia. N Engl J Med 2015 May 7;372(19):1832-1843.
- Kujovich JL. Evaluation of Anemia. Obstet Gynecol Clin North Am. 2016 Jun;43(2):247-264.
- Short MW, Domagalski JE. Iron deficiency anemia: evaluation and management. Am Fam Physician 2013 Jan 15;87(2):98-104.
- Kolnagou A, Kontoghiorghe CN, Kontoghiorghe GJ. New targeted therapies and diagnostic methods for iron overload diseases. Front Biosci (Schol Ed). 2018 Jan 1;10(1):1-20.
- Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Tietz NW, Rinker AD, Morrison SR. When Is a Serum Iron Really a Serum Iron? A Follow-up Study on the Status of Iron Measurements in Serum. Clin Chem 1996;42(1):109-111.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Löhr B, El-Samalouti V, Junge W, et al. Reference Range Study for Various Parameters on Roche Clinical Chemistry Analyzers. Clin Lab 2009;55:465-471.
- Burtis CA, Ashwood ER, Bruns DE (eds.). Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed. St Louis, Missouri; Elsevier Saunders 2006;1190.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.




A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

cobas[®]

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2024, Roche Diagnostics

 0123

 Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.roche.com
 +800 5505 6606

