

ONLINE DAT Opiates II**Order information**

REF	CONTENT	Analyzer(s) on which cobas c pack(s) can be used
04490894 190	ONLINE DAT Opiates II (200 tests)	System-ID 07 6949 5 cobas c 311, cobas c 501/502
03304671 190	Preciset DAT Plus I calibrators CAL 1-6 (6 x 5 mL)	Codes 431-436
03304680 190	Preciset DAT Plus II calibrators CAL 1-6 (6 x 5 mL)	Codes 437-442
03304698 190	C.f.a.s. DAT Qualitative Plus (6 x 5 mL)	
04590856 190	C.f.a.s. DAT Qualitative Plus Clinical (3 x 5 mL)	Code 699
03312950 190	Control Set DAT I (for 2000 ng/mL assay)	
	PreciPos DAT Set I (2 x 10 mL)	
	PreciNeg DAT Set I (2 x 10 mL)	
03312968 190	Control Set DAT II (for 300 ng/mL assay)	
	PreciPos DAT Set II (2 x 10 mL)	
	PreciNeg DAT Set II (2 x 10 mL)	
04500873 190	Control Set DAT Clinical (for 300 ng/mL assay)	
	PreciPos DAT Clinical (2 x 10 mL)	
	PreciNeg DAT Clinical (2 x 10 mL)	

English**System information**

For **cobas c 311/501** analyzers:

OP3Q2: ACN 497: for qualitative assay, 300 ng/mL

OP2Q2: ACN 495: for qualitative assay, 2000 ng/mL

OP3S2: ACN 498: for semiquantitative assay, 300 ng/mL

OP2S2: ACN 496: for semiquantitative assay, 2000 ng/mL

OP3QC: ACN 794: for qualitative assay, 300 ng/mL; using C.f.a.s. DAT Qualitative Plus Clinical

For **cobas c 502** analyzer:

OP3Q2: ACN 8497: for qualitative assay, 300 ng/mL

OP2Q2: ACN 8495: for qualitative assay, 2000 ng/mL

OP3S2: ACN 8498: for semiquantitative assay, 300 ng/mL

OP2S2: ACN 8496: for semiquantitative assay, 2000 ng/mL

OP3QC: ACN 8794: for qualitative assay, 300 ng/mL; using C.f.a.s. DAT Qualitative Plus Clinical

Intended use

Opiates II (OPI2) is an in vitro diagnostic test for the qualitative and semiquantitative detection of morphine and its metabolites in human urine on Roche/Hitachi **cobas c** systems at cutoff concentrations of 300 ng/mL and 2000 ng/mL. Semiquantitative test results may be obtained that permit laboratories to assess assay performance as part of a quality control program. Semiquantitative assays are intended to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as gas chromatography/mass spectrometry (GC-MS).

Opiates II provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. GC-MS is the preferred confirmatory method.¹ Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary

Morphine, a natural product of the opium poppy, is a narcotic analgesic used for centuries as a medicine for the relief of severe pain. Extracted from opium obtained from the poppy's resin, morphine may, in turn, be further chemically refined to heroin (the more potent, diacetylated analog of the parent drug). These chemically similar "opiates" reduce sensitivity to physical and psychological stimuli, dulling pain, fear and anxiety. Users are usually lethargic and indifferent. Accompanying effects may include constriction of the pupils, itching, constipation, nausea, vomiting, and respiratory depression. Death by overdose, usually resulting from dose miscalculation or dose-strength variability, is caused by respiratory failure.^{2,3,4}

The opiates are usually administered intravenously or subcutaneously, but may also be smoked or sniffed. Upon entering the circulation, they tend to

concentrate in the lungs, spleen, kidneys, and liver; lower concentrations are found in the body's musculature and central nervous system. A variety of pathways are involved in the body's detoxification of the opiates, including the removal of chemical side groups (dealkylation), addition of hydroxyl groups, hydrolytic breakdown, and conjugation to glucuronic acid (a common body sugar).⁵ Morphine is excreted in the urine as morphine-3-glucuronide, unchanged free morphine, and other minor metabolites. Although some opiate metabolites appear in the bile and feces, urinary excretion is the primary route of elimination.^{1,6}

The opiates produce strong physical dependence; withdrawal symptoms can begin to appear within a few hours of the last dose and may continue for 5-10 days. The addict may pursue continued opiate use as much to avoid the discomfort of withdrawal as to achieve the desired insensate euphoria.^{7,8}

Test principle

The assay is based on the kinetic interaction of microparticles in a solution (KIMS)^{9,10} as measured by changes in light transmission. In the absence of sample drug, soluble drug conjugates bind to antibody-bound microparticles, causing the formation of particle aggregates. As the aggregation reaction proceeds in the absence of sample drug, the absorbance increases.

When a urine sample contains the drug in question, this drug competes with the drug derivative conjugate for microparticle-bound antibody. Antibody bound to sample drug is no longer available to promote particle aggregation, and subsequent particle lattice formation is inhibited. The presence of sample drug diminishes the increasing absorbance in proportion to the concentration of drug in the sample. Sample drug content is determined relative to the value obtained for a known cutoff concentration of drug.¹¹

Reagents - working solutions

R1 Conjugated morphine derivative; buffer; bovine serum albumin; 0.09 % sodium azide

R2 Microparticles attached to morphine antibody (mouse monoclonal); buffer; bovine serum albumin; 0.09 % sodium azide

R1 is in position B and R2 is in position C.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent handling

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer: 8 weeks

Do not freeze.**Specimen collection and preparation**

Only the specimens listed below were tested and found acceptable.

Urine: Collect urine samples in clean glass or plastic containers. Fresh urine specimens do not require any special handling or pretreatment, but an effort should be made to keep pipetted samples free of gross debris. Samples should be within the normal physiological pH range of 5-8. No additives or preservatives are required. It is recommended that urine specimens be stored at 2-8 °C and tested within 5 days of collection.¹²

For prolonged storage, freezing of the sample is recommended.

Centrifuge highly turbid specimens before testing.

See the limitations and interferences section for details about possible sample interferences.

Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

Adulteration or dilution of the sample can cause erroneous results. If adulteration is suspected, another sample should be collected. Specimen validity testing is required for specimens collected under the *Mandatory Guidelines for Federal Workplace Drug Testing Programs*.¹³

CAUTION: Specimen dilutions should only be used to interpret results of Calc.? and Samp.? alarms, or when estimating concentration in preparation for GC-MS. Dilution results are not intended for patient values. Dilution procedures, when used, should be validated.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for urine

Deselect Automatic Rerun for these applications in the Utility menu, Application screen, Range tab.

cobas c 311 test definition - 300 ng/mL cutoff assay

	Semiquantitative	Qualitative
Assay type	2-Point End	2-Point End
Reaction time / Assay points	10 / 8-22	10 / 8-22
Wavelength (sub/main)	– /570 nm	– /570 nm
Reaction direction	Increase	Increase
Unit	ng/mL	mAbs
Reagent pipetting		Diluent (H ₂ O)
R1	100 µL	–
R2	41 µL	–
Sample volumes	Sample	Sample dilution

		Sample	Diluent (H ₂ O)
Normal	6 µL	–	–
Decreased	6 µL	–	–
Increased	6 µL	–	–

cobas c 501/502 test definition - 300 ng/mL cutoff assay

	Semiquantitative	Qualitative	
Assay type	2-Point End	2-Point End	
Reaction time / Assay points	10 / 13-31	10 / 13-31	
Wavelength (sub/main)	– /570 nm	– /570 nm	
Reaction direction	Increase	Increase	
Unit	ng/mL	mAbs	
Reagent pipetting		Diluent (H ₂ O)	
R1	100 µL	–	
R2	41 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	6 µL	–	–
Decreased	6 µL	–	–
Increased	6 µL	–	–

cobas c 311 test definition - 2000 ng/mL cutoff assay

	Semiquantitative	Qualitative	
Assay type	2-Point End	2-Point End	
Reaction time / Assay points	10 / 8-22	10 / 8-22	
Wavelength (sub/main)	– /570 nm	– /570 nm	
Reaction direction	Increase	Increase	
Unit	ng/mL	mAbs	
Reagent pipetting		Diluent (H ₂ O)	
R1	100 µL	–	
R2	41 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	2 µL	–	–
Decreased	2 µL	–	–
Increased	2 µL	–	–

cobas c 501/502 test definition - 2000 ng/mL cutoff assay

	Semiquantitative	Qualitative	
Assay type	2-Point End	2-Point End	
Reaction time / Assay points	10 / 13-31	10 / 13-31	
Wavelength (sub/main)	– /570 nm	– /570 nm	
Reaction direction	Increase	Increase	
Unit	ng/mL	mAbs	
Reagent pipetting		Diluent (H ₂ O)	
R1	100 µL	–	
R2	41 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H ₂ O)
Normal	2 µL	–	–
Decreased	2 µL	–	–

Increased 2 µL – –

Calibration

Calibrators	<i>Semiquantitative applications</i> <i>300 ng/mL cutoff assay</i> S1-6: Preciset DAT Plus II calibrators, CAL 1-6 0, 150, 300, 600, 1000, 2000 ng/mL <i>2000 ng/mL cutoff assay</i> S1-6: Preciset DAT Plus I calibrators, CAL 1-6 0, 600, 1000, 2000, 4000, 8000 ng/mL <i>Qualitative applications</i> <i>300 ng/mL cutoff assay</i> S1: C.f.a.s. DAT Qualitative Plus Clinical or Preciset DAT Plus II calibrator - CAL 3 300 ng/mL <i>2000 ng/mL cutoff assay</i> S1: Preciset DAT Plus I calibrator - CAL 4 2000 ng/mL The drug concentrations of the calibrators have been verified by GC-MS.
Calibration K Factor	For the qualitative applications, enter the K Factor as -1000 into the Calibration menu, Status screen, Calibration Result window.
Calibration mode	<i>Semiquantitative applications</i> Result Calculation Mode (RCM) ^a <i>Qualitative applications</i> Linear
Calibration frequency	Full (semiquantitative) or blank (qualitative) calibration - after reagent lot change - as required following quality control procedures

a) See Results section.

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against a primary reference method (GC-MS).

Quality control

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Drug concentrations of the Control Set DAT I, II, and Clinical have been verified by GC-MS.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Results

For the qualitative assay, the cutoff calibrator is used as a reference in distinguishing between preliminary positive and negative samples. Samples producing a positive or "0" absorbance value are considered preliminary positive. Preliminary positive samples are flagged with >Test. Samples producing a negative absorbance value are considered negative. Negative samples are preceded by a minus sign.

The semiquantitation of preliminary positive results should only be used by laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as GC-MS. It also permits the

laboratory to establish quality control procedures and assess control performance.

For the semiquantitative assay, the analyzer computer constructs a calibration curve from absorbance measurements of the standards using a 4 parameter logit-log fitting function (RCM). The logit-log function fits a smooth line through the data points. The analyzer computer uses absorbance measurements of samples to calculate drug or drug metabolite concentration by interpolation of the logit-log fitting function.

NOTE: If a result of Calc.? or Samp.? alarm is obtained, review the Reaction Monitor data for the sample and compare with the Reaction Monitor data for the highest calibrator. The most likely cause is a high concentration of the analyte in the sample, in which case the absorbance value for the sample will be less than that of the highest calibrator. Make an appropriate dilution of the sample using the 0 ng/mL calibrator and rerun the sample. A normal drug-free urine may be substituted for the 0 ng/mL calibrator if the urine and procedure have been validated by the laboratory. To ensure that the sample was not over-diluted, the diluted result, prior to multiplying by the dilution factor, must be at least half the analyte cutoff value. If the diluted result falls below half the analyte cutoff value, repeat the sample with a smaller dilution. A dilution that produces a result closest to the analyte cutoff is the most accurate estimation. To estimate the preliminary positive sample's concentration, multiply the result by the appropriate dilution factor. Dilutions should only be used to interpret results of Calc.? or Samp.? alarms, or when estimating concentration in preparation for GC-MS.

Use caution when reporting results as there are various factors that influence a urine test result, such as fluid intake and other biological factors.

As with any sensitive test for drugs of abuse on automated clinical chemistry analyzers, the possibility exists for analyte carry-over from a sample with an extremely high concentration to a normal (negative) sample which immediately follows it.

Preliminary positive results should be confirmed by another method.

Limitations - interference

See the "Specific performance data" section of this document for information on substances tested with this assay. There is the possibility that other substances and/or factors may interfere with the test and cause erroneous results (e.g., technical or procedural errors).

A preliminary positive result with this assay indicates the presence of opiates and/or their metabolites in urine. It does not measure the level of intoxication.

Interfering substances were added to drug free urine at the concentration listed below. These samples were then spiked to 300 ng/mL using a morphine stock solution. Samples were tested on a Roche/Hitachi 917 analyzer and the following results were obtained:

Substance	Concentration Tested	% Morphine Recovery
Acetone	1 %	98
Ascorbic acid	1.5 %	97
Bilirubin	0.25 mg/mL	95
Creatinine	5 mg/mL	95
Ethanol	1 %	100
Glucose	2 %	97
Hemoglobin	7.5 g/L	99
Human albumin	0.5 %	96
Oxalic acid	2 mg/mL	93
Sodium chloride	0.5 M	84
Sodium chloride	1 M	78
Urea	6 %	94

Urine levels of MgSO₄ greater than 400 mg/dL (33.2 mmol/L) were found to interfere with the assay. The results were obtained on a **cobas c 501** analyzer.

Interfering substances were added to drug free urine at the concentration listed below. These samples were then spiked to 2000 ng/mL using a morphine stock solution. Samples were tested on a Roche/Hitachi 917 analyzer and the following results were obtained:

Substance	Concentration Tested	% Morphine Recovery	Cutoff (300)	Number tested	Correct results	Confidence level
Acetone	1 %	99	0.75x	100	100	> 95 % negative reading
Ascorbic acid	1.5 %	96	1.25x	100	100	> 95 % positive reading
Bilirubin	0.25 mg/mL	98				
Creatinine	5 mg/mL	100				
Ethanol	1 %	96				
Glucose	2 %	98				
Hemoglobin	7.5 g/L	101				
Human albumin	0.5 %	96				
Oxalic acid	2 mg/mL	96				
Sodium chloride	0.5 M	95				
Sodium chloride	1 M	91				
Urea	6 %	97				

Urine levels of MgSO₄ up to 600 mg/dL (49.9 mmol/L) do not interfere with the assay. The results were obtained on a **cobas c 501** analyzer.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOH-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Expected values*Qualitative assay*

Results of this assay distinguish preliminary positive (≥ 300 ng/mL or ≥ 2000 ng/mL depending on the cutoff) from negative samples only. The amount of drug detected in a preliminary positive sample cannot be estimated.

Semiquantitative assay

Results of this assay yield only approximate cumulative concentrations of the drug and its metabolites (see "Analytical specificity" section).

Specific performance data

Representative performance data on the Roche/Hitachi analyzer are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined in an internal protocol by running a series of morphine calibrator and controls (repeatability $n = 20$, intermediate precision $n = 100$). The following results were obtained on a **cobas c 501** analyzer.

Semiquantitative precision - 300 ng/mL

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>ng/mL</i>	<i>ng/mL</i>	<i>%</i>
Level 1	225	7.1	3.1
Level 2	301	10.0	3.3
Level 3	385	12.8	3.3
<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>ng/mL</i>	<i>ng/mL</i>	<i>%</i>
Level 1	227	9.4	4.2
Level 2	305	12.0	3.9
Level 3	393	14.4	3.7

Qualitative precision - 300 ng/mL

Cutoff (300)	Number tested	Correct results	Confidence level
0.75x	100	100	> 95 % negative reading
1.25x	100	100	> 95 % positive reading

Semiquantitative precision - 2000 ng/mL

<i>Repeatability</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>ng/mL</i>	<i>ng/mL</i>	<i>%</i>
Level 1	1480	35.3	2.4
Level 2	2006	43.0	2.1
Level 3	2523	55.0	2.2
<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>ng/mL</i>	<i>ng/mL</i>	<i>%</i>
Level 1	1479	44.1	3.0
Level 2	2025	57.6	2.8
Level 3	2518	57.5	2.3

Qualitative precision - 2000 ng/mL

Cutoff (2000)	Number tested	Correct results	Confidence level
0.75x	100	100	> 95 % negative reading
1.25x	100	100	> 95 % positive reading

Accuracy

100 urine samples, obtained from a clinical laboratory where they screened negative in a drug test panel, were evaluated with the Opiates II assay. 100 % of these normal urines were negative relative to the 300 ng/mL and 2000 ng/mL cutoffs. 70 samples, obtained from a clinical laboratory where they screened preliminary positive with a commercially available immunoassay and were subsequently confirmed positive by GC-MS, were evaluated with the Opiates II assay. 100 % of these samples were positive relative to the 300 ng/mL cutoff. 54 samples, obtained from a clinical laboratory where they screened preliminary positive with a commercially available immunoassay and were subsequently confirmed positive by GC-MS, were evaluated with the Opiates II assay. 100 % of these samples were positive relative to the 2000 ng/mL cutoff. In addition, positive urine samples were diluted with drug-free urine. For each cutoff (300 ng/mL and 2000 ng/mL), 10 positive samples were diluted to obtain drug concentrations less than the respective cutoffs. For each cutoff (300 ng/mL and 2000 ng/mL), the same 10 positive samples were diluted to obtain drug concentrations greater than the respective cutoffs. Data from the accuracy studies described above that fell within the near cutoff value ranges were combined with data generated from diluted positive samples. The following results were obtained with the Opiates II assay on the Roche/Hitachi 917 analyzer relative to the GC-MS values.

Opiates II Clinical Correlation (Cutoff = 300 ng/mL)

		Negative Samples	GC-MS values (ng/mL) ^b		
			Near Cutoff		825-48247
			40-253	301-794	
Roche/Hitachi 917 analyzer	+	0	5	7	68
	-	100	8	2	0

^b GC-MS values are represented by the sum of morphine and codeine and do not include all metabolites.

Opiates II Clinical Correlation (Cutoff = 2000 ng/mL)

		Negative Samples	GC-MS values (ng/mL) ^c		
			Near Cutoff		3254-48247
			153-1982	2051-3220	
Roche/Hitachi 917 analyzer	+	0	4	18	42
	-	100	10	0	0

^c GC-MS values are represented by the sum of morphine and codeine and do not include all metabolites.

Additional clinical samples were evaluated with this assay on a **cobas c 501** analyzer and a Roche/Hitachi 917 analyzer. 100 urine samples, obtained from a clinical laboratory where they screened negative in a drug test panel, were evaluated with the Opiates II assay. 100 % of these normal urines were negative for both cutoffs relative to the Roche/Hitachi 917 analyzer. 72 urine samples for the 300 ng/mL cutoff and 48 urine samples for the 2000 ng/mL cutoff, obtained from a clinical laboratory where they screened preliminary positive with a commercially available immunoassay and were subsequently confirmed by GC-MS, were evaluated with the Opiates II assay. At the 300 ng/mL cutoff, 100 % of the samples were positive on the **cobas c 501** analyzer and 97 % of the samples were positive on the Roche/Hitachi 917 analyzer. At the 2000 ng/mL cutoff, 100 % of the samples were positive on both the **cobas c 501** analyzer and the Roche/Hitachi 917 analyzer.

Opiates II Correlation (Cutoff = 300 ng/mL)

		Roche/Hitachi 917 analyzer	
		+	-
cobas c 501 analyzer	+	70	2
	-	0	100

Opiates II Correlation (Cutoff = 2000 ng/mL)

		Roche/Hitachi 917 analyzer	
		+	-
cobas c 501 analyzer	+	48	0
	-	0	100

Analytical specificity

The specificity of this assay for structurally similar compounds was determined by generating inhibition curves for each of the compounds listed and determining the approximate quantity of each compound that is equivalent in assay reactivity to a 300 ng/mL and a 2000 ng/mL assay cutoff. The following results were obtained on a Roche/Hitachi 917 analyzer.

Compound	ng/mL Equivalent to 300 ng/mL Morphine	Approximate % Cross-reactivity
Codeine	224	134
Ethyl morphine	297	101
Diacetylmorphine	366	82
6-Acetylmorphine	386	78
Dihydrocodeine	510	59
Morphine-3-glucuronide	552	54
Dihydromorphine*	937	32
Hydrocodone	1086	28
Thebaine	1210	25
Hydromorphone	1425	21
<i>n</i> -Norcodeine	18590	2
Oxycodone	> 75000	< 0.4
Meperidine	> 100000	< 0.3
Fentanyl*	> 150000	< 0.2

* Results were obtained on a **cobas c 501** analyzer

Compound	ng/mL Equivalent to 2000 ng/mL Morphine	Approximate % Cross-reactivity
Codeine	1541	130
Ethyl morphine	2474	81

6-Acetylmorphine	2598	77
Diacetylmorphine	2915	69
Dihydrocodeine	3170	63
Morphine-3-glucuronide	3785	53
Hydrocodone	7166	28
Dihydromorphine*	7393	27
Thebaine	7579	26
Hydromorphone	10768	19
<i>n</i> -Norcodeine	99264	2
Oxycodone	> 670000	< 0.3
Meperidine	> 670000	< 0.3
Fentanyl*	> 1000000	< 0.2

* Results were obtained on a **cobas c 501** analyzer

Drug interference

The following compounds were prepared in aliquots of pooled normal human urine to yield a final concentration of 100000 ng/mL. None of these compounds gave values in the assay that were greater than 0.5 % cross-reactivity.

Acetaminophen	Ibuprofen
Acetylsalicylic acid	Imipramine
Aminopyrine	Isoproterenol
Amitriptyline	Ketamine
Amobarbital	Lidocaine
<i>d</i> -Amphetamine	LSD ^d
<i>l</i> -Amphetamine	Melanin
Ampicillin	Methadone
Ascorbic acid	<i>d</i> -Methamphetamine
Aspartame	<i>l</i> -Methamphetamine
Atropine	Methaqualone
Benzocaine	Methylphenidate
Benzoylcegonine (cocaine metabolite)	Methyprylon
Benzphetamine	Naloxone
Butabarbital	Naltrexone
Caffeine	Naproxen
Calcium hypochlorite	Niacinamide
Cannabidiol	Norethindrone
Chlordiazepoxide	<i>l</i> -Norpseudoephedrine
Chloroquine	Oxazepam
Chlorpheniramine	Penicillin G
Chlorpromazine	Pentobarbital
Cocaine	Phencyclidine
Dextromethorphan	Phenobarbital
Dextropropoxyphene	Phenothiazine
Diazepam	Phenylbutazone
Diphenhydramine	<i>d</i> -Phenylpropanolamine
Diphenylhydantoin	Phenylpropanolamine
Ecgonine	Procaine
Ecgonine methyl ester	Promethazine
<i>d</i> -Ephedrine	<i>d</i> -Pseudoephedrine
<i>d,l</i> -Ephedrine	<i>l</i> -Pseudoephedrine

/Ephedrine	Quinidine
Epinephrine	Quinine
Erythromycin	Secobarbital
Estriol	Sulindac
Fenoprofen	Tetracycline
Furosemide	Δ ⁹ THC-9-carboxylic acid ^e
Gentisic acid	Tetrahydrozoline
Glutethimide	Trifluoperazine
Guaiaicol glycerol ether	Verapamil
Hydrochlorothiazide	

d) LSD was tested at 2500 ng/mL.

e) Δ⁹ THC-9-carboxylic acid was tested at 10000 ng/mL.

The cross-reactivity for rifampin was tested with the Opiates II assay. The results obtained were 11.0 % and 15.7 % for the 300 ng/mL and 2000 ng/mL cutoffs, respectively.

References

- 1 Karch SB, ed. Drug Abuse Handbook. Boca Raton, FL: CRC Press LLC 1998.
- 2 Papadopoulos CN, Keats AS. Studies of analgesic drugs. VI. Comparative respiratory depressant activity of phenazocine and morphine. Clin Pharmacol Ther 1961;2:8.
- 3 Louria DB, Hensle T, Rose J. The major medical complications of heroin addiction. Ann Intern Med 1967;67(1):1-22.
- 4 Dripps RD, Comroe JH Jr. Clinical studies on morphine I: The immediate effect of morphine administered intravenously and intramuscularly upon the respiration of normal man. Anesthesiology 1945;6:462.
- 5 Way EL, Adler TK. The pharmacological implications of the fate of morphine and its surrogates. Pharmacol Rev 1960 Dec;12:383-446.
- 6 Boerner U. The metabolism of morphine and heroin in man. Drug Metab Rev 1975;4(1):39-73.
- 7 Blum K. Handbook of Abusable Drugs. 1st ed. New York, NY: Gardner Press Inc 1984.
- 8 Jaffe JH. Drug addiction and drug abuse. In: Goodman L, Gilman A, eds. The Pharmacological Basis of Therapeutics. 4th ed. New York: Macmillan 1970.
- 9 Armbruster DA, Schwarzhoff RH, Pierce BL, et al. Method comparison of EMIT II and ONLINE with RIA for drug screening. J Forensic Sci 1993;38:1326-1341.
- 10 Armbruster DA, Schwarzhoff RH, Hubster EC, et al. Enzyme immunoassay, kinetic microparticle immunoassay, radioimmunoassay, and fluorescence polarization immunoassay compared for drugs-of-abuse screening. Clin Chem 1993;39:2137-2146.
- 11 Bates M, Brandle J, Casaretto E, et al. An Abuscreen immunoassay for opiates in urine on the COBAS MIRA automated analyzer. Amer Acad Forensic Sci. Abstract 1991;37(6):1000.
- 12 Toxicology and Drug Testing in the Clinical Laboratory; Approved Guideline. 2nd ed. (C52-A2). Clinical and Laboratory Standards Institute 2007;27:33.
- 13 Mandatory Guidelines for Federal Workplace Drug Testing Programs. Fed Regist 2008 Nov 25;73:71858-71907.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog.roche.com for definition of symbols used):



Contents of kit

Volume after reconstitution or mixing

GTIN

Global Trade Item Number

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

ABUSCREEN, COBAS, COBAS C, ONLINE DAT and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2020, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

