


PHOS2

Phosphate (Inorganic) ver.2

cobas[®]**Order information**

REF		CONTENT		Analyzer(s) on which cobas c pack(s) can be used
03183793122	03183793500	Phosphate (Inorganic) ver.2 (250 tests)	System-ID 07 6614 3	cobas c 311 , cobas c 501/502 , COBAS INTEGRA 400 plus

Materials required (but not provided):

		cobas c 311 , cobas c 501/502	COBAS INTEGRA 400 plus
10759350190	Calibrator f.a.s. (12 x 3 mL)	Code 401	System-ID 07 3718 6
12149435122	Precinorm U plus (10 x 3 mL)	Code 300	System-ID 07 7999 7
12149443122	Precipath U plus (10 x 3 mL)	Code 301	System-ID 07 8000 6
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Code 391	System-ID 07 7469 3
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Code 391	System-ID 07 7469 3
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Code 392	System-ID 07 7470 7
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Code 392	System-ID 07 7470 7
04489357190	Diluent NaCl 9 % (50 mL)	System-ID 07 6869 3	n.a.

English**Intended use**

In vitro test for the quantitative determination of phosphorus in human serum, plasma and urine on **cobas c** and COBAS INTEGRA systems.

Summary

Phosphate measurements, performed with this assay, in human serum, plasma and urine are used as an aid in diagnosis and monitoring of phosphate imbalances such as hyper- or hypophosphatemia.

The large majority (85 %) of phosphate is contained in the skeleton combined with calcium as hydroxyapatite, about 15 % is contained in soft tissue and only < 0.1 % in the extracellular fluid. Phosphate homeostasis is a complex process involving the kidneys, intestine, and skeleton. Phosphate occurs in blood in the form of inorganic phosphate and in organically bound phosphoric acid. The small amount of extracellular organic phosphate is found almost exclusively in the form of phospholipids.¹


The ratio of phosphate to calcium in the blood is approximately 6 : 10.¹ An increase in the level of phosphate causes a decrease in the calcium level. The mechanism is influenced by interactions between parathormone and vitamin D. Hyperphosphatemia originates from excessive phosphate intake or renal reabsorption, reduced phosphate excretion or transcellular shifting.² Clinical conditions such as hypoparathyroidism, vitamin D intoxication and most commonly, renal failure with decreased glomerular phosphate filtration (like in chronic kidney disease, CKD), give rise to hyperphosphatemia.^{3,4,5} Hypophosphatemia is the result of inadequate phosphorus intake, reduced intestinal absorption, excessive urinary excretion, or redistribution of phosphate to the intracellular compartments.^{2,6} Clinical conditions such as rickets, hyperparathyroidism and Fanconi's syndrome are associated with hypophosphatemia.^{7,8,9}

The method presented here for the determination of inorganic phosphate is based on the reaction of phosphate with ammonium molybdate to form ammonium phosphomolybdate without reduction. The addition of an accelerator gives rise to a more rapid rate of reaction and the application of sample blanking yields more precise results.

Test principle¹⁰

Molybdate UV.

Inorganic phosphate forms an ammonium phosphomolybdate complex having the formula $(\text{NH}_4)_3[\text{PO}_4(\text{MoO}_3)_{12}]$ with ammonium molybdate in the presence of sulfuric acid.



The concentration of phosphomolybdate formed is directly proportional to the inorganic phosphate concentration and is measured photometrically.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

**Warning**

H290 May be corrosive to metals.

Prevention:

P234 Keep only in original packaging.

Response:

P390 Absorb spillage to prevent material damage.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Reagent handling

Ready for use

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin and K₂-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine

Collect in detergent-free containers. Acidify with hydrochloric acid after collection (pH < 3).^{11,12}

PHOS2

Phosphate (Inorganic) ver.2

Stability in serum/plasma:¹³
 24 hours at 15-25 °C
 4 days at 2-8 °C
 1 year at (-15)-(-25) °C

Freeze only once.

Stability in urine:^{11,12}
 6 months at 2-8 °C (when acidified)
 24-hour urine: Store cooled during collection.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- See "Order information" section
- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

Calculation

The systems automatically calculate the analyte concentration of each sample.

Conversion factors:
 mmol/L x 3.10 = mg/dL
 mmol/L x 31 = mg/L
 mg/L x 0.0323 = mmol/L

Expected values

Serum/plasma

Adults:¹⁴

0.81-1.45 mmol/L (2.5-4.5 mg/dL)

Children:¹⁵

Age	Male	Female
	mmol/L (mg/dL)	mmol/L (mg/dL)
1-30 d	1.25-2.25 (3.9-6.9)	1.40-2.50 (4.3-7.7)
1-12 m	1.15-2.15 (3.5-6.6)	1.20-2.10 (3.7-6.5)
1-3 y	1.00-1.95 (3.1-6.0)	1.10-1.95 (3.4-6.0)
4-6 y	1.05-1.80 (3.3-5.6)	1.05-1.80 (3.2-5.5)
7-9 y	0.95-1.75 (3.0-5.4)	1.00-1.80 (3.1-5.5)
10-12 y	1.05-1.85 (3.2-5.7)	1.05-1.70 (3.3-5.3)
13-15 y	0.95-1.65 (2.9-5.1)	0.90-1.55 (2.8-4.8)
16-18 y	0.85-1.60 (2.7-4.9)	0.80-1.55 (2.5-4.8)

Urine

1st morning urine¹⁶ 13-44 mmol/L (40-136 mg/dL)

24-hour urine¹¹ 13-42 mmol/d (0.4-1.3 g/d)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

cobas c systems

System information

For **cobas c 311** analyzer:

PHOS2: ACN 714 (serum/plasma)

SPHO2: ACN 675 (STAT, reaction time: 7: serum/plasma)

PHO2U: ACN 716 (urine)

SPH2U: ACN 656 (STAT, reaction time: 7: urine)

For **cobas c 501** analyzer:

PHOS2: ACN 714 (serum/plasma/urine)

SPHO2: ACN 675 (STAT, reaction time: 7: serum/plasma/urine)

For **cobas c 502** analyzer:

PHOS2: ACN 8714 (serum/plasma)

SPHO2: ACN 8675 (STAT, reaction time: 7: serum/plasma)

PHO2U: ACN 8716 (urine)

SPH2U: ACN 8656 (STAT, reaction time: 7: urine)

Reagents - working solutions

R1 Sulfuric acid: 0.36 mol/L; detergent

R2 Ammonium molybdate: 3.5 mmol/L; sulfuric acid: 0.36 mol/L; sodium chloride: 150 mmol/L

R1 is in position B and R2 is in position C.

Storage and stability

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

Application for serum and plasma

cobas c 311 test definition

Assay type	2-Point End		
Reaction time / Assay points	10 / 6-32 (STAT 7 / 6-32)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	-	

Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2.5 µL	-	-
Decreased	12.5 µL	15 µL	135 µL
Increased	2.5 µL	-	-

cobas c 501 test definition

Assay type	2-Point End		
Reaction time / Assay points	10 / 10-47 (STAT 7 / 10-47)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	-	

Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	2.5 µL	-	-
Decreased	12.5 µL	15 µL	135 µL

Increased 2.5 µL – –

cobas c 502 test definition

Assay type	2-Point End		
Reaction time / Assay points	10 / 10-47 (STAT 7 / 10-47)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	–	
Sample volumes	Sample	Sample dilution	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2.5 µL	–	–
Decreased	12.5 µL	15 µL	135 µL
Increased	5 µL	–	–

Application for urine**cobas c 311 test definition**

Assay type	2-Point End		
Reaction time / Assay points	10 / 6-32 (STAT 7 / 6-32)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	–	
Sample volumes	Sample	Sample dilution	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2.5 µL	15 µL	150 µL
Decreased	2.5 µL	8 µL	168 µL
Increased	2.5 µL	15 µL	150 µL

cobas c 501 test definition

Assay type	2-Point End		
Reaction time / Assay points	10 / 10-47 (STAT 7 / 10-47)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	–	
Sample volumes	Sample	Sample dilution	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2.5 µL	15 µL	150 µL
Decreased	2.5 µL	8 µL	168 µL
Increased	2.5 µL	15 µL	150 µL

cobas c 502 test definition

Assay type	2-Point End		
Reaction time / Assay points	10 / 10-47 (STAT 7 / 10-47)		
Wavelength (sub/main)	700/340 nm		
Reaction direction	Increase		
Units	mmol/L (mg/dL, mg/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	90 µL	28 µL	
R2	38 µL	–	

Sample volumes	Sample	Sample dilution	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	2.5 µL	15 µL	150 µL
Decreased	2.5 µL	8 µL	168 µL
Increased	5 µL	15 µL	150 µL

Calibration

Calibrators	S1: H ₂ O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration • after reagent lot change • as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against NERL primary reference material.

Quality control*Serum/plasma*

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

Urine

Quantitative urine controls are recommended for routine quality control.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference¹¹

Criterion: Recovery within ± 10 % of initial value at a phosphate concentration of 0.87 mmol/L (2.7 mg/dL).

Serum/plasma

Icterus:¹⁷ No significant interference up to an I index of 40 for conjugated and 60 for unconjugated bilirubin (approximate conjugated bilirubin concentration: 684 µmol/L or 40 mg/dL and approximate unconjugated bilirubin concentration: 1026 µmol/L or 60 mg/dL).

Hemolysis:¹⁷ No significant interference up to an H index of 300 (approximate hemoglobin concentration: 186 µmol/L or 300 mg/dL).

Note: This interference results from inorganic phosphates produced by the action of phosphatases on organic phosphates, both of which are released from the red cells upon hemolysis.

Lipemia (Intralipid):¹⁷ No significant interference up to an L index of 1250. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{18,19}

PHOS2

Phosphate (Inorganic) ver.2

Exception: Phospholipids contained in liposomal drug formulations (eg AmBisome) may be hydrolyzed in the test due to the acidic reaction pH and thus lead to elevated phosphate results.²⁰

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.²¹

Urine

Drugs: No interference was found at therapeutic concentrations using common drug panels.¹⁹

Criterion: Recovery within $\pm 10\%$ of initial value at a phosphate concentration of 13 mmol/L (40.3 mg/dL).

Urea: No significant interference from urea up to a concentration of 1500 mmol/L (9009 mg/dL).

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c 502** analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is required in certain cases.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

Serum/plasma

0.10-6.46 mmol/L (0.31-20.0 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:2 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 2.

Urine

1.1-92 mmol/L (3.4-285 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:2 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 2.

Lower limits of measurement

Lower detection limit of the test

Serum/plasma

0.10 mmol/L (0.31 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Urine

1.1 mmol/L (3.4 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol.

Serum/plasma:

Repeatability (n = 21), intermediate precision (3 aliquots per run, 1 run per day, 21 days);

Urine:

Repeatability (n = 21), intermediate precision (3 aliquots per run, 1 run per day, 10 days).

The following results were obtained on the **cobas c 501** analyzer:

Serum/plasma

Repeatability	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Precinorm U	1.24 (3.84)	0.01 (0.03)	0.7
Precipath U	2.05 (6.36)	0.01 (0.03)	0.6
Human serum 1	2.68 (8.31)	0.02 (0.06)	0.6
Human serum 2	1.56 (4.84)	0.01 (0.03)	0.7

Intermediate precision

	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Precinorm U	1.23 (3.81)	0.02 (0.06)	1.4
Precipath U	2.04 (6.32)	0.02 (0.06)	1.2
Human serum 3	2.67 (8.28)	0.04 (0.12)	1.4
Human serum 4	1.55 (4.81)	0.02 (0.06)	1.4

Urine

Repeatability	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Control Level 1	10.2 (31.6)	0.1 (0.3)	1.4
Control Level 2	19.9 (61.7)	0.2 (0.6)	1.2
Human urine 1	40.9 (127)	0.4 (1)	1.0
Human urine 2	6.25 (19.4)	0.08 (0.2)	1.2

Intermediate precision

	Mean	SD	CV
	mmol/L (mg/dL)	mmol/L (mg/dL)	%
Control Level 1	10.0 (31.0)	0.2 (0.6)	1.6
Control Level 2	19.6 (60.8)	0.3 (0.9)	1.7
Human urine 3	40.4 (125)	0.5 (2)	1.3
Human urine 4	6.23 (19.3)	0.12 (0.4)	2.0

The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

Method comparison

Inorganic phosphate values for human serum, plasma and urine samples obtained on a **cobas c 501** analyzer (y) were compared with those determined using the corresponding reagent on a Roche/Hitachi 917 analyzer (x).

Serum/plasma

Sample size (n) = 150

Passing/Bablok ²²	Linear regression
$y = 1.022x + 0.000$ mmol/L	$y = 1.023x - 0.002$ mmol/L
$\tau = 0.978$	$r = 1.000$

The sample concentrations were between 0.62 and 5.54 mmol/L (1.92 and 17.2 mg/dL).

Urine

Sample size (n) = 145

Passing/Bablok ²²	Linear regression
$y = 0.976x - 0.053$ mmol/L	$y = 0.974x - 0.047$ mmol/L
$\tau = 0.967$	$r = 0.999$

The sample concentrations were between 1.61 and 91.5 mmol/L (4.99 and 284 mg/dL).

PHOS2

Phosphate (Inorganic) ver.2



The data obtained on **cobas c 501** analyzer(s) are representative for **cobas c 311** analyzer(s).

COBAS INTEGRA systems

System information

PHOS2: Test ID 0-614 (serum, plasma)

PHOU2: Test ID 0-514 (urine)

Reagents - working solutions

R1 Sulfuric acid 0.36 mol/L, detergent

SR Ammonium molybdate 3.5 mmol/L, Sulfuric acid 0.36 mol/L,
Sodium chloride 150 mmol/L

R1 is in position B and SR is in position C

Storage and stability

Shelf life at 2-8 °C See expiration date on
cobas c pack label

On-board in use at 10-15 °C 12 weeks

Application for serum, plasma, and urine

Test definition

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction direction	Increase
Wavelength A/B	340/659 nm
Calc. first/last	33/63
Unit	mmol/L

Serum, plasma

Reaction mode R1-S-SR

Urine

Reaction mode D-R1-S-SR

Predilution factor 11

Pipetting parameters

Serum, plasma, and urine		Diluent (H ₂ O)
R1	90 µL	
Sample	2.5 µL	27.5 µL
SR	38 µL	
Total volume	158 µL	

Calibration

Serum, plasma, and urine

Calibrator	Calibrator f.a.s. Use deionized water as zero calibrator.
Calibration mode	Linear regression
Calibration replicate	Duplicate recommended
Calibration interval	Each lot and as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against NERL primary reference material.

For USA: This method has been standardized against NIST traceable primary reference material.

Quality control

Quality control serum, plasma PreciControl ClinChem Multi 1
Precinorm U plus

PreciControl ClinChem Multi 2
Precipath U plus

Quality control urine Quantitative urine controls are recommended for routine quality control.

Control interval 24 hours recommended

Control sequence User defined

Control after calibration Recommended

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Limitations - interference

Criterion: Recovery within ± 10 % of initial value.

Serum, plasma

Icterus:¹⁷ No significant interference up to an I index of 51 (approximate conjugated bilirubin concentration: 872 µmol/L or 51 mg/dL). No significant interference with unconjugated bilirubin.

Hemolysis:¹⁷ No significant interference up to an H index of 420 (approximate hemoglobin concentration: 261 µmol/L or 420 mg/dL).

Lipemia (Intralipid):¹⁷ No significant interference up to an Intralipid level of 1000 mg/dL. There is poor correlation between turbidity and triglycerides concentration.

Drugs: No interference was found at therapeutic concentrations using common drug panels.^{18,19} **Exception:** Phospholipids contained in liposomal drug formulations (e.g. AmBisome) may be hydrolyzed in the test due to the acidic reaction pH and thus lead to elevated phosphate results.²⁰

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.²¹

Urine

Drugs: No interference was found at therapeutic concentrations using common drug panels.¹⁹

Criterion: Recovery within ± 10 % of initial value at a phosphate concentration of 13 mmol/L (40.3 mg/dL).

Urea: No significant interference from urea up to a concentration of 1500 mmol/L (9009 mg/dL).

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the CLEAN Method Sheet for further instructions and for the latest version of the Extra wash cycle list.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

Serum/plasma

0.1-6.46 mmol/L (0.31-20 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 10.

Urine

1.1-92 mmol/L (3.41-285 mg/dL)

PHOS2

Phosphate (Inorganic) ver.2



Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 5.

Lower limits of measurement

Serum/plasma

Lower detection limit of the test:

0.1 mmol/L (0.31 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Urine

Lower detection limit of the test:

1.1 mmol/L (3.41 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 21 days). The following results were obtained on the COBAS INTEGRA 700 analyzer:

Serum and plasma

	Level 1	Level 2
Mean	1.17 mmol/L (3.63 mg/dL)	2.01 mmol/L (6.23 mg/dL)
CV repeatability	1.3 %	1.4 %
Mean	1.17 mmol/L (3.63 mg/dL)	2.00 mmol/L (6.20 mg/dL)
CV intermediate precision	2.5 %	2.4 %

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (1 aliquot per run, 1 run per day, 10 days). The following results were obtained:

Urine

	Level 1	Level 2
Mean	13.9 mmol/L (43.1 mg/dL)	27.6 mmol/L (85.6 mg/dL)
CV repeatability	1.0 %	0.7 %
Mean	13.9 mmol/L (43.1 mg/dL)	27.7 mmol/L (85.9 mg/dL)
CV intermediate precision	1.7 %	1.1 %

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

Method comparison

Inorganic phosphate values obtained on a COBAS INTEGRA 700 analyzer with the COBAS INTEGRA Phosphate (Inorganic) ver.2 reagent (y) were compared with those determined using the same reagent on a Roche/Hitachi 917 analyzer (x) and with the previous reagent (PHOS) on a COBAS INTEGRA 700 analyzer (x).

Serum and plasma

Roche/Hitachi 917 analyzer	Sample size (n) = 100
Passing/Bablok ²²	Linear regression
$y = 1.043x + 0.022$ mmol/L	$y = 1.040x + 0.025$ mmol/L

$\tau = 0.955$ $r = 1.000$

SD (md 95) = 0.040 $Sy.x = 0.018$

The sample concentrations were between 0.572 to 5.69 mmol/L (1.77 to 17.7 mg/dL).

COBAS INTEGRA 700 analyzer Sample size (n) = 96

Passing/Bablok²² Linear regression

$y = 1.029x - 0.047$ mmol/L $y = 1.040x - 0.067$ mmol/L

$\tau = 0.942$ $r = 0.999$

SD (md 95) = 0.077 $Sy.x = 0.032$

The sample concentrations were between 0.619 to 4.76 mmol/L (1.92 to 14.9 mg/dL).

Urine

Roche/Hitachi 917 analyzer Sample size (n) = 86

Passing/Bablok²² Linear regression

$y = 1.052x - 0.0235$ mmol/L $y = 1.044x - 0.028$ mmol/L

$\tau = 0.983$ $r = 1.000$

SD (md 95) = 0.743 $Sy.x = 0.349$

The sample concentrations were between 6.08 to 89.4 mmol/L (18.9 to 277 mg/dL).

COBAS INTEGRA 700 analyzer Sample size (n) = 68

Passing/Bablok²² Linear regression

$y = 1.000x - 0.399$ mmol/L $y = 1.002x - 0.405$ mmol/L

$\tau = 0.989$ $r = 1.000$

SD (md 95) = 0.396 $Sy.x = 0.180$

The sample concentrations were between 6.08 to 44.8 mmol/L (18.9 to 139 mg/dL).

The data obtained on COBAS INTEGRA 700 analyzer(s) are representative for COBAS INTEGRA 400 analyzer(s).

References

- 1 Ellervik C, Halsall DJ, Nygaard B. Thyroid Disorders. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 57, p. 806-845.e13.
- 2 Manghat P, Sodi R, Swaminathan R. Phosphate homeostasis and disorders. Ann Clin Biochem. 2014 Nov;51(Pt 6):631-56. doi: 10.1177/0004563214521399.
- 3 Bilezikian JP. Hypoparathyroidism. J Clin Endocrinol Metab. 2020 Jun 1;105(6):1722-36. doi: 10.1210/clinem/dgaa113.
- 4 Razzaque MS. Can adverse effects of excessive vitamin D supplementation occur without developing hypervitaminosis D? J Steroid Biochem Mol Biol. 2018 Jun;180:81-86. doi: 10.1016/j.jsbmb.2017.07.006.
- 5 Ketteler M, Block GA, Evenepoel P, et al. Executive summary of the 2017 KDIGO Chronic Kidney Disease-Mineral and Bone Disorder (CKD-MBD) Guideline Update: what's changed and why it matters. Kidney Int. 2017 Jul;92(1):26-36. doi: 10.1016/j.kint.2017.04.006. Erratum in: Kidney Int. 2017 Dec;92(6):1558.
- 6 Tebben PJ. Hypophosphatemia: A Practical Guide to Evaluation and Management. Endocr Pract. 2022 Oct;28(10):1091-1099. doi: 10.1016/j.eprac.2022.07.005.
- 7 Takashi Y, Kawanami D, Fukumoto S. FGF23 and Hypophosphatemic Rickets/Osteomalacia. Curr Osteoporos Rep. 2021 Dec;19(6):669-675. doi: 10.1007/s11914-021-00709-4.
- 8 Arnold A, Dennison E, Kovacs CS, et al. Hormonal regulation of biominerallization. Nat Rev Endocrinol. 2021 May;17(5):261-275. doi: 10.1038/s41574-021-00477-2.

PHOS2

Phosphate (Inorganic) ver.2

cobas[®]




- 9 Klootwijk ED, Reichold M, Unwin RJ, et al. Renal Fanconi syndrome: taking a proximal look at the nephron. *Nephrol Dial Transplant*. 2015 Sep;30(9):1456-60. doi: 10.1093/ndt/gfu377.
- 10 Henry R ed. *Clinical Chemistry: Principles and Technics*, 2nd ed. New York, NY: Harper & Row 1974;723.
- 11 Tietz NW, ed. *Clinical Guide to Laboratory Tests*, 4th ed. Philadelphia.WB Saunders Co 2006;852-855.
- 12 NCCLS GP-16A2, *Urineanalysis and Collection, Transportation and Preservation of Urine specimens*, 2nd edition 2001.
- 13 Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- 14 Burtis CA, Ashwood ER, Bruns DE (eds.). *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 4th ed. St Louis, Missouri; Elsevier Saunders 2006;2290.
- 15 Soldin JS, Bruignara C, Wong EC. *Pediatric Reference Intervals*. AACC Press. 2005, 5th ed., p. 153.
- 16 Krieg M, Gunsser KJ, Steinhagen-Thiessen E, et al. Vergleichende quantitative Analytik klinisch-chemischer Kenngrößen im 24-Stunden-Urin und Morgenerin. *J Clin Chem Clin Biochem* 1986 Nov;24(11):863-869.
- 17 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- 18 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- 19 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- 20 Lane JW, Rehak NN, Hortin GL, et al. Pseudohyperphosphatemia associated with high-dose liposomal amphotericin B therapy. *Clin Chim Acta* 2008;387:145-149.
- 21 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- 22 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see navifyportal.roche.com for definition of symbols used):

	Contents of kit
	Volume for reconstitution
	Global Trade Item Number

Rx only For USA: Federal law restricts this device to sale by or on the order of a physician.

CE 0123

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

+800 5505 6606



COBAS, NAVIFY, PRECICONTROL, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2023, Roche Diagnostics