

**Total Protein Urine/CSF Gen. 3****Order information**

REF	CONTENT	Analyzer(s) on which <b>cobas c</b> pack(s) can be used
0333825 190	Total Protein Urine/CSF Gen.3 150 tests	System-ID 07 6763 8 Roche/Hitachi <b>cobas c</b> 311, <b>cobasc</b> 501/502
03121305 122	C.f.a.s. PUC (5 x 1 mL)	Code 489
03121313 122	Precinorm PUC (4 x 3 mL)	Code 240
03121291 122	Precipath PUC (4 x 3 mL)	Code 241
04489357 190	Diluent NaCl 9 % (50 mL)	System-ID 07 6869 3

**English****System information**

For **cobas c** 311/501 analyzers:

**TPU3:** ACN 708

**TPC3:** ACN 402

For **cobas c** 502 analyzer:

**TPU3:** ACN 8708

**TPC3:** ACN 8402

**Intended use**

In vitro test for the quantitative determination of protein in human urine and cerebrospinal fluid on Roche/Hitachi **cobas c** systems.

**Summary**

Protein measurements in urine are used in the diagnosis and treatment of disease conditions such as renal or heart diseases, or thyroid disorders, which are characterized by proteinuria or albuminuria. Cerebrospinal fluid (CSF) protein measurements are used in the diagnosis and treatment of conditions such as meningitis, brain tumors and infections of the central nervous system.<sup>1</sup>

Urine is formed by ultrafiltration of plasma across the glomerular capillary wall. Proteins with a relative molecular mass > 40000 are almost completely retained, while smaller substances easily enter the glomerular filtrate. Most CSF protein originates by diffusion from plasma across the blood-CSF barrier. Elevated levels occur as a result of increased permeability of the blood-CSF barrier or with increased local synthesis of immunoglobulins.

Turbidimetric methods using trichloroacetic acid (TCA) or sulfosalicylic acid (SSA) precipitate proteins in the sample depending on their size; the resulting turbidity may be unstable and flocculate. Reagents of dye-binding methods such as Coomassie blue and pyrogallol red-molybdate react with proteins depending on their amino acid composition, but may stain glass and plastic ware. Due to their reaction mechanisms all methods, turbidimetric and colorimetric, exhibit different sensitivities to various proteins, especially to protein fragments such as Bence Jones proteins<sup>2</sup> and small proteins such as  $\alpha$ 1-microglobulin.

The Roche Diagnostics Urinary/CSF Protein assay is based on the method described by Iwata and Nishikaze,<sup>3</sup> later modified by Luxton, Patel, Keir, and Thompson.<sup>4</sup> In this method, benzethonium chloride reacts with protein in a basic medium to produce a turbidity that is more stable and evenly distributed than that observed with the SSA or TCA methodologies. This assay shows an underrecovery of  $\gamma$ -globulin compared to albumin of about 30 %,<sup>5</sup> and no interference from magnesium ions due to the addition of EDTA.

**Test principle**

Turbidimetric method.

The sample is preincubated in an alkaline solution containing EDTA, which denatures the protein and eliminates interference from magnesium ions. Benzethonium chloride is then added, producing turbidity.

**Reagents - working solutions**

**R1** Sodium hydroxide: 677 mmol/L; EDTA-Na: 74 mmol/L

**R2** Benzethonium chloride: 32 mmol/L

R1 is in position B and R2 is in position C.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:



Danger

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

H412 Harmful to aquatic life with long lasting effects.

**Prevention:**

P234 Keep only in original container.

P264 Wash skin thoroughly after handling.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

**Response:**

P301 + P330 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. + P331

P303 + P361 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. + P353

P304 + P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. + P338

P310 Immediately call a POISON CENTER or doctor/physician.

P337 + P313 If eye irritation persists: Get medical advice/attention.

P363 Wash contaminated clothing before reuse.

P390 Absorb spillage to prevent material damage.

**Storage:**

P405 Store locked up.

P406 Store in corrosive resistant stainless steel container with a resistant inner liner.

**Disposal:**

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling primarily follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

**Reagent handling**

Ready for use

**Storage and stability***TPUC3*

Shelf life at 15-25 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 6 weeks

*Diluent NaCl 9 %*

Shelf life at 2-8 °C: See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer: 12 weeks

**Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

*Urine*

Use random or 24-hour urine specimens. Use no preservatives. Refrigerate specimen during collection.

*CSF*

No special additives are required. Blood in a CSF specimen invalidates the protein value.<sup>1</sup>

Samples for urinary/CSF protein should be collected before fluorescein is given or at least 24 hours later.<sup>6</sup>

**Note:** Urine, CSF and control samples with a protein concentration above 7000 mg/L must not be measured with TPUC3 as this may clog the instrument lines.

*Stability:*<sup>7</sup>

Urine:	1 day at 15-25 °C
	7 days at 2-8 °C
	1 month at (-15)-(-25) °C
CSF:	1 day at 15-25 °C
	6 days at 2-8 °C
	> 1 year at (-15)-(-25) °C

Centrifuge samples containing precipitates before performing the assay.

Non centrifuged samples may produce elevated results.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

- See "Order information" section
- General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for urine and CSF****cobas c 311 test definition**

Assay type	2-Point End
Reaction time / Assay points	10 / 6-14
Wavelength (sub/main)	700/505 nm
Reaction direction	Increase
Units	mg/L (mg/dL, g/L)

## Reagent pipetting

R 1	100 µL	Diluent (H <sub>2</sub> O)	–
-----	--------	----------------------------	---

R 2	40 µL	–
-----	-------	---

<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>

Normal	6 µL	–	–
--------	------	---	---

Decreased	2 µL	–	–
-----------	------	---	---

Increased	6 µL	–	–
-----------	------	---	---

**cobas c 501 test definition**

Assay type	2-Point End
------------	-------------

Reaction time /	10 / 10-30
-----------------	------------

Assay points	
--------------	--

Wavelength (sub/main)	700/505 nm
-----------------------	------------

Reaction direction	Increase
--------------------	----------

Units	mg/L (mg/dL, g/L)
-------	-------------------

Reagent pipetting	Diluent (H <sub>2</sub> O)
-------------------	----------------------------

R1	100 µL	–
----	--------	---

R2	40 µL	–
----	-------	---

<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	6 µL	–	–
Decreased	2 µL	–	–
Increased	6 µL	–	–

**cobas c 502 test definition**

Assay type	2-Point End
------------	-------------

Reaction time /	10 / 10-30
-----------------	------------

Assay points	
--------------	--

Wavelength (sub/main)	700/505 nm
-----------------------	------------

Reaction direction	Increase
--------------------	----------

Units	mg/L (mg/dL, g/L)
-------	-------------------

Reagent pipetting	Diluent (H <sub>2</sub> O)
-------------------	----------------------------

R1	100 µL	–
----	--------	---

R2	40 µL	–
----	-------	---

<i>Sample volumes</i>	<i>Sample</i>	<i>Sample dilution</i>	
		<i>Sample</i>	<i>Diluent (NaCl)</i>
Normal	6 µL	–	–
Decreased	2 µL	–	–
Increased	12 µL	–	–

**Calibration**

Calibrators	S1: H <sub>2</sub> O
-------------	----------------------

	S2-S6: C.f.a.s. PUC
--	---------------------

Multiply the lot-specific C.f.a.s. PUC calibrator values by the factors given below to determine the standard concentrations for the 6-point calibration curve.

	S2: 0.025	S5: 0.250
	S3: 0.050	S6: 1.0
	S4: 0.125	
Calibration mode	RCM	
Calibration frequency	Full calibration	
	- after reagent lot change	
	- as required following quality control procedures	

Traceability:<sup>8</sup> This method has been standardized against a primary standard traceable to NIST.

**Quality control**

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factors:	mg/L x 0.1 = mg/dL
	mg/L x 0.001 = g/L

To calculate 24-hour urine protein excretion:  
mg/L x total volume (liters per 24 hours) = mg/day.

**Limitations - interference**

Criterion: Recovery within  $\pm 10\%$  of initial value at a total protein concentration of 120 mg/L (12 mg/dL; 0.12 g/L).

High dose hook-effect: Sample results with high total protein concentrations above the measuring range up to 100000 mg/L will be flagged by the instrument with > TEST or > ABS.

**Urine**

Icterus: No significant interference up to a concentration of 342  $\mu$ mol/L (20 mg/dL) for conjugated bilirubin.

Hemolysis: Hemoglobin interferes.<sup>9</sup>

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>10</sup>

Exception: Levodopa, methyl dopa and Na<sub>2</sub>-cefoxitin cause artificially high total protein results and calcium dobesilate causes artificially low protein results.

Other: Patient samples containing > 8 g/L of organically bound iodine from Radiopaque media (e.g. Hexabrix) may have falsely elevated results.

High levels of homogentisic acid can be found in urine of patients with the rare genetic disorder Alkaptonuria.<sup>11</sup> Homogentisic acid in urine samples at concentrations > 0.6 mmol/L can cause false results.

The administration of gelatin-based plasma replacements can lead to increased urine protein values.

**CSF**

Hemolysis: Hemoglobin interferes.<sup>9</sup>

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. The latest version of the carry-over evasion list can be found with the NaOHD-SMS-SmpCln1+2-SCCS Method Sheets. For further instructions refer to the operator's manual. **cobas c** 502 analyzer: All special wash programming necessary for avoiding carry-over is available via the **cobas** link, manual input is not required.

Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

**Limits and ranges****Measuring range**

40-2000 mg/L (4-200 mg/dL; 0.04-2 g/L)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 3.

**Lower limits of measurement****Lower detection limit of the test**

40 mg/L (4 mg/dL; 0.04 g/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, n = 21).

**Expected values**

Urine: <sup>12</sup>	24 h:	< 140 mg/24 h*
	random:	< 150 mg/L*
	* Values obtained from centrifuged samples	
CSF:	reference range acc. to Tietz:	150-450 mg/L (15-45 mg/dL) <sup>13</sup>
	reference range acc. to Thomas:	200-400 mg/L (20-40 mg/dL) <sup>14</sup>

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 21 days). The following results were obtained:

**Urine**

Repeatability	Mean	SD	CV
	mg/L (mg/dL)	mg/L (mg/dL)	%
Precinorm PUC	159 (15.9)	1 (0.1)	0.7
Precipath PUC	1576 (158)	8 (0.8)	0.5
Human urine 1	101 (10.1)	1 (0.1)	1.0
Human urine 2	191 (19.1)	4 (0.4)	2.2

**Intermediate precision**

	Mean	SD	CV
	mg/L (mg/dL)	mg/L (mg/dL)	%
Precinorm PUC	156 (15.6)	2 (0.2)	1.5
Precipath PUC	1482 (148)	8 (0.8)	0.5
Human urine 3	106 (10.6)	2 (0.2)	1.6
Human urine 4	154 (15.4)	1 (0.1)	0.9

**CSF**

Repeatability	Mean	SD	CV
	mg/L (mg/dL)	mg/L (mg/dL)	%
Control Level 1	281 (28.1)	4 (0.4)	1.5
Control Level 2	691 (69.1)	4 (0.4)	0.6
Human CSF 1	355 (35.5)	4 (0.4)	1.1

**Total Protein Urine/CSF Gen. 3**

Human CSF 2	517 (51.7)	5 (0.5)	1.0
<i>Intermediate precision</i>	<i>Mean</i>	<i>SD</i>	<i>CV</i>
	<i>mg/L (mg/dL)</i>	<i>mg/L (mg/dL)</i>	<i>%</i>
Control Level 1	272 (27.2)	4 (0.4)	1.6
Control Level 2	660 (66.0)	6 (0.6)	0.9
Human CSF 3	349 (34.9)	4 (0.4)	1.2
Human CSF 4	501 (50.1)	7 (0.7)	1.5

**Method comparison**

Total protein values for human urine and CSF samples obtained on a Roche/Hitachi **cobas c 501** analyzer (y) were compared with those determined with the corresponding reagent on a Roche/Hitachi 917 analyzer (x).

**Urine**

Sample size (n) = 70

Passing/Bablok <sup>15</sup>	Linear regression
$y = 0.985x + 6.23 \text{ mg/L}$	$y = 0.988x + 5.35 \text{ mg/L}$
$\tau = 0.970$	$r = 1.000$

The sample concentrations were between 47.0 and 1887 mg/L (4.70 and 189 mg/dL).

**CSF**

Sample size (n) = 86

Passing/Bablok <sup>15</sup>	Linear regression
$y = 1.015x - 7.51 \text{ mg/L}$	$y = 1.010x - 5.23 \text{ mg/L}$
$\tau = 0.975$	$r = 0.999$

The sample concentrations were between 53.0 and 1087 mg/L (5.30 and 109 mg/dL).

**References**

- 1 Tietz NW. Fundamentals of Clinical Chemistry, 3rd ed. Pa: WB Saunders Co 1987:336.
- 2 Boege F. Bence Jones-Proteine. J Lab Med 1999;23(9):477-482.
- 3 Iwata J, Nishikaze O. New micro-turbidimetric method for determination of protein in cerebrospinal fluid and urine. Clin Chem 1979;25(7):1317-1319.
- 4 Luxton RW, Patel P, Keir G, et al. A micro-method for measuring total protein in cerebrospinal fluid by using benzethonium chloride in microtiter plate wells. Clin Chem 1989;35(8):1731-1734.
- 5 Hohnadel DC, Koller A. Urine protein total. In: Pesce AJ, Kaplan LA, editors. Methods in clinical chemistry, St. Louis, Mosby 1987.
- 6 Koumantakis G. Fluorescein Interference with Urinary Creatinine and Protein Measurements. Clin Chem 1991;37/10:1799.
- 7 WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
- 8 Standard Reference Materials from NERL, traceable to NIST (National Institute of Standards and Technology).
- 9 Yilmaz FM, Yücel D. Effect of Addition of Hemolysate on Urine and Cerebrospinal Fluid Assays for Protein. Clin Chem 2006;52:152-153.
- 10 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 11 Phornphutkul C, Introne WJ, Perry MB, et al. Natural History of Alkaptonuria. N Engl J Med 2002;347(26):2111-2121.
- 12 Junge W, Wilke B, Halabi A, et al. Reference Intervals for Total Protein in Collected and Random Urine using the Benzethonium Chloride Method [Abstract]. Clin Chem 2006;52:157.
- 13 Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders Company 1995;518-523.

- 14 Thomas L. Labor und Diagnose, 6. Auflage, TH-Books Verlagsgesellschaft mbH, Frankfurt/Main 2005;930-934.
- 15 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

**Symbols**

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT

Contents of kit



Volume after reconstitution or mixing

GTIN

Global Trade Item Number

COBAS, COBAS C, PRECINORM and PRECIPATH are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

