

**ONLINE TDM Vancomycin Gen.3****Order information**

REF		CONTENT		Analyzer(s) on which <b>cobas c</b> pack(s) can be used
08445605190	08445605500	ONLINE TDM Vancomycin Gen.3 (100 tests)	System-ID 2121 002	<b>cobas c 303, cobas c 503</b>
08058849190	08445605500	ONLINE TDM Vancomycin Gen.3 (200 tests)	System-ID 2121 003	<b>cobas c 303, cobas c 503</b>

Materials required (but not provided):

03375790190	Preciset TDM I CAL A-F (1 x 5 mL) Diluent (1 x 10 mL)	Codes 20691-20696	
04521536190	TDM Control Set Level I (2 x 5 mL) Level II (2 x 5 mL) Level III (2 x 5 mL)	Code 20310 Code 20311 Code 20312	

**English****System information****VANC3:** ACN 21210 (100 tests; System ID 2121 002)**VANC30:** ACN 21211 (200 tests; System ID 2121 003)**Intended use**In vitro test for the quantitative determination of vancomycin in serum and plasma on **cobas c** systems.**Summary**

Vancomycin measurements performed with this assay, in human serum and plasma, are used for monitoring vancomycin treatment to ensure appropriate therapy.

Vancomycin is a tricyclic glycopeptide isolated from *Streptomyces orientalis* that inhibits the synthesis of the cell wall in sensitive bacteria (mainly gram-positive bacteria and some gram-negative cocci).<sup>1,2</sup> It displays concentration-independent activity<sup>2</sup> and is used because of its efficacy against methicillin-resistant staphylococci and corynebacteria.<sup>3,4</sup>Vancomycin is available as intravenous and oral formulations and is indicated in the treatment of infections such as complicated skin and soft tissue infections, community and hospital acquired pneumonia and *Clostridium difficile* infection, and for the perioperative antibacterial prophylaxis in patients that are at high risk of developing bacterial endocarditis.<sup>2,3,4</sup>Monitoring of vancomycin serum or plasma levels is used to ascertain clinical efficacy and to limit potentially dose-dependent serious side effects, such as nephro- and ototoxicity<sup>1,2,4</sup> and is recommended by guidelines and clinical societies.<sup>5,6,7,8,9,10,11</sup>**Test principle**

The assay is based on the kinetic interaction of microparticles in a solution (KIMS). Vancomycin antibody is covalently coupled to microparticles and the drug derivative is linked to a macromolecule. The kinetic interaction of microparticles in solutions, photometrically detected by turbidity measurements is induced by binding of drug-conjugate to the antibody on the microparticles and is inhibited by the presence of vancomycin in the sample. A competitive reaction takes place between the drug conjugate and vancomycin in the serum sample for binding to the vancomycin antibody on the microparticles. The resulting turbidity is indirectly proportional to the amount of drug present in the sample.

**Reagents - working solutions****R1** Vancomycin conjugate; piperazine-N,N'-bis(2-ethanesulfonic acid) (PIPES) buffer, pH 7.2; preservative; stabilizer**R3** Anti-vancomycin antibody (mouse monoclonal); latex microparticle; 3-(N-morpholino)propane sulfonic acid (MOPS) buffer, pH 7.2; stabilizer

R1 is in position B and R3 is in position C.

**Precautions and warnings**

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

**Reagent handling**

Ready for use

Carefully invert reagent container several times prior to use to ensure that the reagent components are mixed.

**Storage and stability**

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label

On-board in use and refrigerated on the analyzer:

12 weeks

**Do not freeze.****Specimen collection and preparation**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: K<sub>2</sub>- or K<sub>3</sub>-EDTA, lithium heparin.

Sample collection tubes containing separating gel have not been verified for use.

Stability:

48 hours capped at 15-25 °C

14 days capped at 2-8 °C

12 months capped at -20 °C (± 5 °C)

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Do not induce foaming of specimens. Specimens can be repeatedly frozen and thawed up to 5 times.

Invert thawed specimens several times prior to testing.

Usual sampling time varies dependent upon desired measurement of peak or trough values.<sup>12</sup>**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum and plasma****Test definition**

Reporting time	10 min		
Wavelength (sub/main)	800/600 nm		
Reagent pipetting		Diluent (H <sub>2</sub> O)	
R1	65 µL	–	
R3	46 µL	–	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (H <sub>2</sub> O)
Normal	1.3 µL	–	–
Decreased	1.3 µL	–	–
Increased	1.3 µL	–	–

For further information about the assay test definitions refer to the application parameters setting screen of the corresponding analyzer and assay.

**Calibration**

Calibrators	S1-6: Preciset TDM I calibrators
Calibration mode	Non-linear
Calibration frequency	Full calibration - after reagent lot change - every 2 weeks on-board - as required following quality control procedures

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against USP reference standards.<sup>13</sup> The calibrators are prepared to contain known quantities of vancomycin in normal human serum.

**Quality control**

For quality control, use control materials as listed in the "Order information" section.

In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. It is recommended to perform quality control always after lot calibration and subsequently at least every 12 weeks.

Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

**cobas c** systems automatically calculate the analyte concentration of each sample in the unit µg/mL (µmol/L, mg/L).

Conversion factors: <sup>14</sup>	µg/mL x 0.690 = µmol/L
	µg/mL x 1.0 = mg/L

**Limitations - interference**

Criterion: Recovery within ± 1 µg/mL (0.690 µmol/L) of initial values of samples ≤ 10 µg/mL (5.18 µmol/L) and within ± 10 % for samples > 10 µg/mL (5.18 µmol/L).

**Serum/Plasma**

Icterus:<sup>15</sup> No significant interference up to an I index of 60 for conjugated bilirubin and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dL or 1026 µmol/L).

Hemolysis:<sup>15</sup> No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dL or 622 µmol/L).

Lipemia (Intralipid):<sup>15</sup> No significant interference up to an L index of 1000. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Triglycerides: No significant interference from triglycerides up to a concentration of 1000 mg/dL (11.4 mmol/L).

Rheumatoid factors: No significant interference from rheumatoid factors up to a concentration of 1200 IU/mL.

Total protein: No significant interference from total protein in the concentration range of 2-12 g/dL.

As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample, which could cause falsely lowered results.

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>16</sup>

In very rare cases (less than 1 reported case per 1000000 tests) certain immunoglobulins can unspecifically interfere with the agglutination reaction leading to unreliable results.

**Note:** A test result flagged with ">Kin3" indicates unusual reaction kinetics. There is a high possibility that the sample contains an interfering substance which accelerates the reaction kinetics. For such samples it is not possible to report a reliable analyte concentration with this assay.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on **cobas c** systems. All special wash programming necessary for avoiding carry-over is available via the **cobas** link. The latest version of the carry-over evasion list can be found with the NaOHD/SMS/SCCS Method Sheet. For further instructions, refer to the operator's manual.

**Limits and ranges****Measuring range**

4.0-80.0 µg/mL (2.76-55.2 µmol/L)

Manually dilute samples above the measuring range 1 + 1 with the Preciset TDM I diluent (0 µg/mL) and reassay. Multiply the result by 2 to obtain the specimen value.

**Lower limits of measurement**

*Limit of Blank, Limit of Detection and Limit of Quantitation*

Limit of Blank	= 1.0 µg/mL (0.69 µmol/L)
Limit of Detection	= 1.5 µg/mL (1.04 µmol/L)
Limit of Quantitation	= 4.0 µg/mL (2.76 µmol/L)

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95<sup>th</sup> percentile value from n ≥ 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration that can be reproducibly measured with a total error of 20 %. It has been determined using low concentration vancomycin samples.

### Expected values

The practice of routine monitoring and adjustment of serum vancomycin concentrations has been the subject of intense debate for many years.<sup>17</sup> Historically trough concentrations between 5 to 10 µg/mL and peak concentrations between 20 to 40 µg/mL were generally accepted for therapeutic effectiveness.<sup>17,18,14</sup> The increased prevalence of resistant organisms, increasing vancomycin minimum inhibitory concentrations in target pathogens (particularly MRSA) and vancomycin failures have prompted more aggressive vancomycin dosing practices and recommendations.<sup>17,19</sup> Therefore, current guidelines recommend higher trough concentrations in the range of 10-15 µg/mL for uncomplicated MRSA bacteremia and even 15-20 µg/mL in cases of sustained MRSA bacteremia or endocarditis and other severe invasive MRSA infections (i.e. prosthetic joint infections, hospital-acquired pneumonia or central nervous system infections).<sup>7,17</sup> However, higher doses of vancomycin used have been associated with significantly higher vancomycin trough levels, acute renal failure and ototoxicity.<sup>19,20,21,22</sup> The decision to target increased vancomycin trough concentrations should be based on an assessment of the severity of the infection and must consider the risk associated with increased vancomycin levels.

Expected values reflect the data and information provided in the reference and do not necessarily represent therapeutic recommendations and/or dosage instructions. For therapeutic recommendations and dosage instructions refer to applicable guidelines and the full prescription information of the drug.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

### Specific performance data

Representative performance data on the analyzers are given below. These data represent the performance of the analytical procedure itself.

Results obtained in individual laboratories may differ due to heterogenous sample materials, aging of analyzer components and mixture of reagents running on the analyzer.

### Precision

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n = 84) and intermediate precision (2 aliquots per run, 2 runs per day, 21 days). Results for repeatability and intermediate precision were obtained on the **cobas c 503** analyzer.

### Serum/Plasma

Repeatability	Mean	SD	CV
	µg/mL	µg/mL	%
TDMC1 <sup>a)</sup>	7.49	0.129	1.7
TDMC2 <sup>b)</sup>	21.7	0.210	1.0
TDMC3 <sup>c)</sup>	34.6	0.304	0.9
Human serum 1	5.18	0.120	2.3
Human serum 2	8.29	0.135	1.6
Human serum 3	27.9	0.263	0.9
Human serum 4	43.1	0.472	1.1
Human serum 5	74.1	0.763	1.0
Intermediate precision	Mean	SD	CV
	µg/mL	µg/mL	%
TDMC1 <sup>a)</sup>	7.49	0.197	2.6
TDMC2 <sup>b)</sup>	21.7	0.318	1.5
TDMC3 <sup>c)</sup>	34.6	0.447	1.3
Human serum 1	5.18	0.191	3.7
Human serum 2	8.29	0.229	2.8
Human serum 3	27.9	0.414	1.5
Human serum 4	43.1	0.617	1.4
Human serum 5	72.7	1.41	1.9

a) TDM Control Set Level I

b) TDM Control Set Level II

c) TDM Control Set Level III

The data obtained on **cobas c 503** analyzer(s) are representative for **cobas c 303** analyzer(s).

### Method comparison

#### Serum/plasma

Vancomycin values for human serum samples obtained on a **cobas c 503** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 119

Passing/Bablok<sup>23</sup>

$y = 1.004x - 0.345 \mu\text{g/mL}$

$r = 0.966$

Linear regression

$y = 1.005x - 0.437 \mu\text{g/mL}$

$r = 0.996$

The sample concentrations were between 4.29 and 79.8 µg/mL.

#### Serum/plasma

Vancomycin values for human serum samples obtained on a **cobas c 303** analyzer (y) were compared with those determined using the corresponding reagent on a **cobas c 501** analyzer (x).

Sample size (n) = 113

Passing/Bablok<sup>23</sup>

$y = 1.007x - 0.427 \mu\text{g/mL}$

$r = 0.977$

Linear regression

$y = 1.022x - 0.829 \mu\text{g/mL}$

$r = 0.999$

The sample concentrations were between 4.90 and 76.6 µg/mL.

### Analytical specificity

The following compounds were tested for cross-reactivity.

Compound	Concentration	%
	Tested (µg/mL)	Cross-reactivity
Acyclovir	50	ND
Amikacin	100	ND
Amphotericin B	10	ND
Aztreonam	450	ND
Caffeine	60	ND
CDP-1	20	ND
Cefazoline	500	ND
Cefotaxime	300	ND
Chloramphenicol	60	ND
Ciprofloxacin	12	ND
Cisplatin	15	ND
Clindamycin	50	ND
Cyclosporine	3	ND
Digoxin	0.009	ND
Epinephrine	1	ND
Erythromycin	60	ND
Ethacrynic acid	1.5	ND
Flucytosine	300	ND
Furosemide	60	ND
Fusidic acid	600	ND
Gentamicin	30	ND
Imipenem	250	ND
Methicillin	250	ND

Methotrexate	455	ND
Metronidazole	150	ND
Netilmicin	30	ND
Nitroprusside	90	ND
Penicillin G	36	ND
Pentamidine	1.5	ND
Phenobarbital	150	ND
Rifampin	60	ND
Salicylate	750	ND
Sulphamethoxazole	400	ND
Theophylline	60	ND
Tobramycin	30	ND
Trimethoprim	40	ND

ND = Not detectable

Tests were performed on 16 drugs. No significant interference with the assay was found.

Acetaminophen	Heparin
Acetyl cysteine	Ibuprofen
Acetylsalicylic acid	Levodopa
Ampicillin-Na	Methyldopa + 1.5 H <sub>2</sub> O
Ascorbic acid	Metronidazole
Cefoxitin	Phenylbutazone
Cyclosporine	Rifampicin
Doxycycline (Tetracycline)	Theophylline

## References

- McLawnon RW. Chapter 10 - Guidelines for the Monitoring of Vancomycin, Aminoglycosides and Certain Antibiotics. In: Dasgupta A, editor. Therapeutic Drug Monitoring. Boston: Academic Press 2012; p. 197-218.
- Flynn Pharma Ltd. Vancocin Powder for Solution, vancomycin hydrochloride - Drug information [revised 2022 June; cited 2023 August 18]. Available from: <https://www.medicines.org.uk/emc/product/6407/smpc>.
- Milone MC, Shaw LM. Therapeutic Drugs and Their Management, In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, editors. Tietz Textbook of Laboratory Medicine, Saunders Elsevier, Philadelphia, 7th edition, 2023, chapter 42, p. 420-453.e9.
- Brown NM, Goodman AL, Horner C, et al. Treatment of methicillin-resistant Staphylococcus aureus (MRSA): updated guidelines from the UK. JAC Antimicrob Resist 2021 Feb 3;3(1):dlaa114.
- Rybak MJ, Le J, Lodise T, et al. Executive Summary: Therapeutic Monitoring of Vancomycin for Serious Methicillin-Resistant Staphylococcus aureus Infections: A Revised Consensus Guideline and Review of the American Society of Health-System Pharmacists, the Infectious Diseases Society of America, the Pediatric Infectious Diseases Society, and the Society of Infectious Diseases Pharmacists. J Pediatr Infect Dis Soc 2020 Jul 13;9(3):281-284.
- Reuter SE, Stocker SL, Alffenaar JC, et al. Optimal Practice for Vancomycin Therapeutic Drug Monitoring: Position Statement From the Anti-infectives Committee of the International Association of Therapeutic Drug Monitoring and Clinical Toxicology. Ther Drug Monit 2022 Feb 1;44(1):121-132.
- Liu C, Bayer A, Cosgrove SE, et al. Clinical practice guidelines by the Infectious Diseases Society of America for the treatment of methicillin-resistant Staphylococcus aureus infections in adults and children. Clin Infect Dis 2011;52(3):285-292.
- Matsumoto K, Takesue Y, Ohmagari N, et al. Practice guidelines for therapeutic drug monitoring of vancomycin: a consensus review of the Japanese Society of Chemotherapy and the Japanese Society of Therapeutic Drug monitoring. J Infect Chemother 2013 Jun;19(3):365-380.
- Abdul-Aziz MH, Alffenaar JC, Bassetti M, et al. Antimicrobial therapeutic drug monitoring in critically ill adult patients: a Position Paper. Intensive Care Med 2020 Jun;46(6):1127-1153.
- Delgado V, Ajmone Marsan N, de Waha S, et al. 2023 ESC Guidelines for the management of endocarditis. Eur Heart J 2023 Aug 25;ehad193.
- Baddour LM, Wilson WR, Bayer AS, et al. Infective Endocarditis in Adults: Diagnosis, Antimicrobial Therapy, and Management of Complications: A Scientific Statement for Healthcare Professionals From the American Heart Association. Circulation 2015 Oct 13;132(15):1435-1486.
- Rotschafer JC, Crossley K, Zaske DE, et al. Pharmacokinetics of vancomycin: observations in 28 patients and dosage recommendations. Antimicrob Agents Chemother 1982;22:391-394.
- USP 39-NF (U.S. Pharmacopeia National Formulary) 2016:6329-6330.
- Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th edition (Elsevier) 2006;2186.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Rybak M, Lomaestro B, Rotschafer JC, et al. Therapeutic monitoring of vancomycin in adult patients: a consensus review of the American Society of Health-System Pharmacists, the Infectious Diseases Society of America, and the Society of Infectious Diseases Pharmacists. Am J Health Syst Pharm 2009;66(1):82-98.
- Hammett-Stabler CA, Johns T. Laboratory guidelines for monitoring of antimicrobial drugs. Clin Chem 1998;44(5):1129-1140.
- Pritchard L, Baker C, Leggett J, et al. Increasing Vancomycin in Serum Trough Concentrations and Incidence of Nephrotoxicity. Am J Med 2010;123:1143-1149.
- Hidayat LK, Hsu DI, Quist R, et al. High-dose vancomycin therapy for methicillin-resistant Staphylococcus aureus infections: efficacy and toxicity. Arch Intern Med 2006;166(19):2138-2144.
- Jeffres MN, Isakow W, Doherty JA, et al. A retrospective analysis of possible renal toxicity associated with vancomycin in patients with health care-associated methicillin-resistant Staphylococcus aureus pneumonia. Clin Ther 2007;29(6):1107-1115.
- Lodise TP, Lomaestro B, Graves J, et al. Larger vancomycin doses (at least four grams per day) are associated with an increased incidence of nephrotoxicity. Antimicrob Agents Chemother 2008;52(4):1330-1336.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety & Performance Report can be found here: <https://ec.europa.eu/tools/eudamed>

## Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see [navifyportal.roche.com](http://navifyportal.roche.com) for definition of symbols used):

**CONTENT**

Contents of kit

08445605500V7.0

# VANC3

ONLINE TDM Vancomycin Gen.3

cobas®



Volume for reconstitution

GTIN

Global Trade Item Number

Rx only

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

COBAS, NAVIFY, ONLINE TDM and PRECISET are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2023, Roche Diagnostics

CE 0123



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

+800 5505 6606

